
Gebrauchsanweisung / Instructions for use / Istruzioni per l'uso / Mode d'emploi / Gebruiksaanwijzing /
 Instrucciones de uso / Käyttöohjeet / Bruksanvisning / Bruksanvisning / Οδηγίες χρήσης / Инструкции за употреба /
 Uputstvo za upotrebu / Návod k použití / Instrukcja obsługi / Instrucțiuni de utilizare / Návod na použitie /
 Navodila za uporabo / Használati utasítás / Upute za korištenje / Naudojimo instrukcija / Lietošanas instrukcija /
 Kasutusjuhend / Инструкция по применению / Kullanım kılavuzu / Instruções de uso / 使用说明 / 取扱説明書 / 사용 설명서

- DE Easy-Click-Kassette für ulricheasyINJECT Max 2M / Max 3
- EN Easy-Click-Cassette for ulricheasyINJECT Max 2M / Max 3
- IT Easy-Click-Cassetta per ulricheasyINJECT Max 2M / Max 3
- FR Easy-Click-Cassette pour ulricheasyINJECT Max 2M/Max 3
- NL Easy-Click-Cassette voor ulricheasyINJECT Max 2M / Max 3
- ES Easy-Click-Casete para ulricheasyINJECT Max 2M/Max 3
- FI Easy-Click-kasetti ulricheasyINJECT Max 2M- / Max 3 -tuotteisiin
- NO Easy-Click-kassett for ulricheasyINJECT Max 2M / Max 3
- SV Easy-Click-kassett för ulricheasyINJECT Max 2M / Max 3
- EL Easy Click-Κασέτα για ulricheasyINJECT Max 2M / Max 3
- BG Easy-Click-касета за ulricheasyINJECT Max 2M / Max 3
- SR Easy-Click-kaseta za ulricheasyINJECT Max 2M / Max 3
- CS Easy-Click-Kazeta pro ulricheasyINJECT Max 2M / Max 3
- PL Easy-Click-Kaseta do urządzenia ulricheasyINJECT Max 2M / Max 3
- RO Casetă Easy-Click pentru ulricheasyINJECT Max 2M / Max 3
- SK Kazeta Easy-Click pre injektor ulricheasyINJECT Max 2M/Max 3
- SL Easy-Click-Kaseta za ulricheasyINJECT Max 2M/Max 3
- HU Easy-Click-Kazetta az ulricheasyINJECT Max 2M / Max 3 készülékhez
- HR Easy-Click-Kaseta za ulricheasyINJECT Max 2M / Max 3
- LT Kasetė „Easy-Click“, skirta „ulricheasyINJECT Max 2M / Max 3“
- LV Easy-Click kasete injektoram ulricheasyINJECT Max 2M / Max 3
- ET Easy-Click-kassett ulricheasyINJECT Max 2M-i / Max 3 jaoks
- RU Easy-Click-Кассета для инжекторов ulricheasyINJECT Max 2M/Max 3
- TR ulricheasyINJECT Max 2M / Max 3 için Easy-Click-Aparatı
- BR Easy-Click-Cassete para ulricheasyINJECT Max 2M / Max 3
- ZH 适用于 ulricheasyINJECT Max 2M / Max 3 的卡带式管路匣
- JA MRI 用注入装置 ulricheasyINJECT Max 2M / Max 3 専用「イーजूクリックカセット」
- KO ulricheasyINJECT Max 2M / Max 3용 Easy-Click 카세트

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EN Instructions for use


About this document

These instructions for use are an integral part of the product and describe how to safely use it as intended.
→ Read the instructions for use before using the product.
→ Keep the instructions for use with the product so that they are accessible at all times.

Labeling and symbols

Labeling	Meaning
 WARNING	Warning! Hazards to persons. Non-compliance can cause death or serious injuries to persons, or damage to the product.
→	Handling instructions.

Symbols on the product and packaging

Symbol	Description	Symbol	Description
	Consult instructions for use		Do not use if package is damaged
	CE mark with identification number of the notified body		Do not use if protective caps are not properly attached
	Catalogue number		Non-pyrogenic
	Batch code		Keep dry
	Quantity		Temperature limit
	Use-by date		Humidity limitation
	Manufacturer		Keep away from sunlight
	Date of manufacture		Do not re-use! Intended for single use only for 24 hours.
	Sterilized using ethylene oxide		Latex free

Abbreviations

Abbreviation	Meaning
CEDM	Contrast Enhanced Digital Mammography
CESM	Contrast Enhanced Spectral Mammography
PET	Positron Emission Tomography

Intended use

The Easy-Click-Cassette for ulricheasyINJECT Max 2M / Max 3 is used for the venous administration of contrast media and physiological saline solutions (NaCl) during diagnostic magnetic resonance imaging examinations (MRI and MRI/PET; XD 10710 / -11 / -12 / -13) and contrast medium enhanced mammography examinations (CEDM/CESM; XD 10711 / -13). The Easy-Click-Cassette connects the injector (on the device side) to the patient tubing. The Easy-Click-Cassette can be used for any number of injections within a 24 hour period.

Contraindications

The Easy-Click-Cassettes for ulricheasyINJECT Max 2M / Max 3 are not intended for the administration of contrast media during high-pressure angiography or other applications for which the device is not indicated.

Product description

- A** Device side with Easy-Click-Cassette
- B** Patient side with patient tubing
- 1** Spike for media supply for NaCl
- 2** Spike with holder for media supply for contrast medium (Max 3 only)
- 3** Spike with holder for media supply for contrast medium
- 4** Protective cap for spike
- 5-7** Clips for sealing the media supply tubings during strip down
- 5** Clip, blue
- 6** Clip, green (Max 3 only)
- 7** Clip, green (Max 2M) / yellow (Max 3)
- 8** Locking catches
- 9** Valve diaphragm
- 10** Media flow direction
- 11** Tubing segment for roller pump
- 12** Pressure measurement cell (in the Easy-Click-Cassette)
- 13** Particle filter
- 14** SafeConnect with protective cap
- 15** RFID chip

The Easy-Click-Cassette consists of a plastic housing with an integrated tubing system (cassette) and spikes for connecting the media containers, a tubing segment as a pump segment, and a valve diaphragm. The patient tubing is connected to the Easy-Click-Cassette via SafeConnect (connector with protection against accidental contact). The Easy-Click-Cassette contains a particle filter, which retains foreign bodies, and a pressure measurement cell for measuring the pressure in the tubing system during operation. Clips, for sealing off the media supply tubings, prevent liquid from running out during replacement of the Easy-Click-Cassette. The Easy-Click-Cassette has two spikes (for the ulricheasyINJECT Max 2M injector) or three spikes (for the ulricheasyINJECT Max 3 injector). The spikes with a holder, for fastening the media containers, are available in 2 different sizes for supply of contrast media. All spikes contain an integrated bacteria filter in the supply air duct.

RFID

The Easy-Click-Cassette contains a microchip with a signature and usage data, which are checked by the device during setup and operation. The ulricheasyINJECT Max 2M / Max 3 injectors must only be operated with original disposables from ulrich medical. When it is filled for the first time, the Easy-Click-Cassette is marked as "used", and from that time forward cannot be used in any other ulricheasyINJECT Max 2M / Max 3. After the maximum service life of 24 hours, the Easy-Click-Cassette is "invalidated", and can no longer be used. During the maximum service life, any number of injections can be performed. When the maximum service life is reached, the current injection can still be completed.

Scope of delivery and area of application

Art. No.	Designation
XD 10710	Easy-Click-Cassette Max 3 holder s (for ulricheasyINJECT MRI injectors)
XD 10711	Easy-Click-Cassette Max 3 holder l (for ulricheasyINJECT MRI injectors)
XD 10712	Easy-Click-Cassette Max 2M holder s (for ulricheasyINJECT MRI injectors)
XD 10713	Easy-Click-Cassette Max 2M holder l (for ulricheasyINJECT MRI injectors)

Product intended for single use

The product is sterilized with ethylene oxide and packed in sterile packaging.

- Do not reuse the product.
- Use the product for a maximum of 24 hours.
- Do not re-sterilize the product.
- Do not use any product from open or damaged sterile packaging.
- Do not use the product after the expiry date.
- Do not use any product that has not been stored under recommended storage conditions.

Single-use product

WARNING! Infection of patient/user due to reuse and/or re-processing.

Reuse and/or re-processing cause the risk of contamination and material fatigue. Fouling, contamination, and/or impaired function can lead to injury, illness, or death.

- Do not re-process product.

Safe handling

- Before using the product, check that it is operable and in proper condition.
- In order to avoid damage due to improper use, use the product in accordance with these instructions for use and the injector's instructions for use.
- Product and accessories should only be used by persons who have the required knowledge, training and experience.

Risk of injury!

The use of tubing combinations or tubing systems not approved by ulrich medical jeopardizes the safety of the patient and/or user and compromises the proper functioning of the injector.

- Use only original products approved by ulrich medical.

Risk of injury! Injection of air, overdose, and flow reduction due to pressure regulation when using tubing extensions.

Tubing extensions jeopardize patient safety and can contaminate the tubing system.

- Do not use tubing extensions (e.g. "Heidelberg extension").

Risk of injury due to infections and contamination!

Improper handling jeopardizes patient safety.

- Handle all sterile components carefully in order to ensure sterility.
- Do not use the product if the protective caps are missing or are not firmly seated.
- After disconnecting the patient tubing, the system must immediately be equipped with new patient tubing, and the protective cap must not be removed from the Luer lock connector until it is connected to the patient.
- Replace and discard the Easy-Click-Cassette after each working day or after 24 hours at the latest.

Risk of injury due to mixing of contrast media when using software option CM-Select!

- When a different type of contrast medium is attached to the injector that is incompatible with the previously used medium: replace the Easy-Click-Cassette.

Special instructions regarding safe use of the injector system in contrast medium enhanced mammography

- Do not use Easy-Click-Cassette Max 3 / Max 2M with holder (s) (XD 10710 / XD 10712).
- Do not plug in contrast media bottles with a volume > 250 mL onto Easy-Click-Cassette Max 3 / Max 2M with holder (l) (XD 10711 / XD 10713).
- Observe flow limitation maximum of 4 mL/s (@ iodine concentration = 370 mg/mL, NaCl and contrast medium temperature = 25°C, 18 G cannula size) (higher iodine concentrations, lower temperatures and cannula sizes > 18 G can reduce max. flow rate and delivery volume).

Use

Preparation

- Remove the product from the packaging.
- Check the product for damaged parts before use.
- Do not use the product if any parts are damaged.

Connecting the product to the injector

- Connect the product only as specified in the injector's instructions for use.

Functional test

- Perform a functional test as specified in the injector's instructions for use.

Operation

- Use the product only as specified in the injector's instructions for use.
- Replace the Easy-Click-Cassette after each working day or after 24 hours at the latest.

Technical service

ulrich GmbH & Co. KG
Injector Systems Customer Service
Buchbrunnweg 12
89081 Ulm / Germany

Tel.: +49 (0)731 9654-0
Fax: +49 (0)731 9654-2808
E-mail:
service@ulrichmedical.com

Technical data

Feature	Value
Pressure tightness	22 bar max.
Material	PUR, PVC, polycarbonate (free from DEHP and latex)
Max. service life	24 hours for multiple patients

Storage

- Store the product in the sterile packaging in a dry, dark and temperature-controlled room.
- Protect the sterile packaged product from environmental factors such as dirt, dust, humidity, heat, and sunlight.

Ambient conditions

	Use	Storage
Temperature	+5°C to +40°C	+15°C to +25°C
Relative humidity	10% to 90% non-condensing	30% to 60% non-condensing

Disposal

- After use, dispose of the product in accordance with the locally applicable regulations.

Reporting of incidents

The user must report all product complaints related to the safety, efficacy or performance of the product to ulrich medical (complaint@ulrichmedical.com) or the local distributor within 72 hours.

If one or more components of the system show(s) signs of malfunction (i. e. one or more of the performance specifications is/are not met) or do/does not provide the intended performance or if there is a suspicion that this may be the case, then ulrich medical (complaint@ulrichmedical.com) or the local distributor should be notified immediately (within 24 hours).

If one or more components of the system has/have ever failed and possibly led to or contributed to the death or severe injury of a patient, then ulrich medical or the local distributor should be notified immediately.

Please use the Medical Device Vigilance Report Form from ulrich medical to report product complaints, malfunctions or product failure. It is available at: www.ulrichmedical.com/vigilancereport.

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ulrich GmbH & Co. KG | Buchbrunnenweg 12 | 89081 Ulm | Germany