
Gebrauchsanweisung / Instructions for use / Istruzioni per l'uso / Mode d'emploi / Gebruiksaanwijzing /  
 Instrucciones de uso / Käyttöohjeet / Bruksanvisning / Bruksanvisning / Οδηγίες χρήσης / Инструкции за употреба /  
 Uputstvo za upotrebu / Návod k použití / Instrukcja obsługi / Instruções de utilizar / Návod na použitie /  
 Navodila za uporabo / Használati utasítás / Upute za korištenje / Naudojimo instrukcijos / Lietošanas instrukcija /  
 Kasutusjuhend / Инструкция по применению / Kullanım talimatları / Instruções de uso / 使用说明 / 取扱説明書 /  
 사용 설명서

- DE Patientenschlauch für ulricheasyINJECT CT/MRT-Injektoren
- EN Patient tubing for ulricheasyINJECT CT/MRI injectors
- IT Tubo del paziente per iniettori ulricheasyINJECT TAC/RMN
- FR Tubulure patient pour injecteurs ulricheasyINJECT pour scanner/IRM
- NL Patiëntenslang voor ulricheasyINJECT CT/MRI injectoren
- ES Tubos del paciente para inyectores ulricheasyINJECT TC/TRM
- FI ulricheasyINJECT TT/MRI -injektorien potilasletku
- NO Pasientslange for ulricheasyINJECT CT/MR-injektorer
- SV Patientslang för ulricheasyINJECT DT/MRT-injektorer
- EL Σωλήνας ασθενούς για εγχυτές αξονικής τομογραφίας (CT)/μαγνητικής τομογραφίας (MRI) ulricheasyINJECT
- BG Пациентни маркучи за инжектори за КТ/ЯМР ulricheasyINJECT
- SR Creva za pacijenta za ubrizgivače za CT/MRI iz serije ulricheasyINJECT
- CS Hadičky pacienta pro injektory ulricheasyINJECT CT/MR
- PL Wężyk pacjenta do wstrzykiwaczy ulricheasyINJECT TK/MRI
- RO Tubulatură pentru pacient pentru injectoare ulricheasyINJECT TC/IRM
- SK Hadička pacienta pre injektory ulricheasyINJECT CT/MRI
- SL Cevke za bolnike za injektorje ulricheasyINJECT CT/MRI
- HU Páciensszerezélék az ulricheasyINJECT CT/MRI injektorokhoz
- HR Cijevi za bolesnika za injektore ulricheasyINJECT CT/MR
- LT Paciento vamzdelis, skirtas „ulricheasyINJECT“ KT / MRT injektoriams
- LV Pacienta caurulīte ulricheasyINJECT CT/MRI injektoriem
- ET Patsiendivoolik ulricheasyINJECTI KT-/MRT-injektoritele
- RU Трубка пациента для КТ/МРТ-инжекторов ulricheasyINJECT
- TR ulricheasyINJECT BT/MRG enjektörler için hasta hortumu
- BR Equipo de transferência de paciente para os injetores de TC/RM ulricheasyINJECT
- ZH 适用于 ulricheasyINJECT CT/MRI 注射器的患者导管
- JA ulricheasyINJECT CT/MRI注入装置用患者チュービングチューブ
- KO ulricheasyINJECT CT/MRI 주사기용 환자 튜브

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### Zu diesem Dokument

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## EN Instructions for use

### About this document

These instructions for use are an integral part of the product and describe how to safely use it as intended.

- Read the instructions for use before using the product.
- Keep the instructions for use with the product so that they are accessible at all times.

### Labeling and symbols

Labeling	Meaning
 <b>WARNING</b>	<b>Warning!</b> Hazards to persons. Non-compliance can cause death or serious injuries to persons, or damage to the product.
→	Handling instructions.

### Abbreviations

Abbreviation	Meaning
CT	Computed tomography
CEDM	Contrast Enhanced Digital Mammography
CESM	Contrast Enhanced Spectral Mammography
MRI	Magnetic resonance imaging
CM	Contrast media
NaCl	Sodium chloride (here: 0.9% saline solution)

### Symbols on the product and packaging

Symbol	Description	Symbol	Description
	Consult instructions for use		Sterilized using ethylene oxide
	CE mark with identification number of the notified body		Do not use if package is damaged

Symbol	Description	Symbol	Description
	Catalogue number		Do not use if protective caps are not properly attached
	Batch code		Non-pyrogenic
	Quantity		Keep dry
	Use-by date		Temperature limit
	Manufacturer		Humidity limitation
	Date of manufacture		Keep away from sunlight
	Do not re-use		Latex free

### Intended use

The patient tubing for ulricheasyINJECT CT/MRI injectors is intended for the venous administration of contrast media and physiological saline solutions (NaCl) during a CT/MRI examination and also during a contrast-enhanced mammography examination (CEDM/CESM). The patient tubing creates the connection between the Easy-Click-Cassette (on the device side) and the patient.

### Contraindications

The patient tubing for ulricheasyINJECT CT/MRI injectors is not intended for administration of contrast medium for high-pressure angiography examinations or for any other use for which it is not indicated.

### Product description

The patient tubing consists of a tubing system with two protective caps and has a male Luer lock connector with an integrated check valve on the patient side. The Easy-Click-Cassette is connected to the SafeConnect. A second check valve is integrated into the tubing. The check valves are self-closing. They open only when a medium flows in the defined direction. In case of counter-pressure or no media flow, the valves remain closed.

- A Device side with Easy-Click-Cassette
- B Patient side
- 1 SafeConnect
- 2 Protective caps for tubing ends
- 3 Male Luer lock, with the first check valve
- 4 Second check valve
- 5 Media flow direction

### RFID

The patient tubing contains a microchip with a signature and usage data, which are checked by the device during setup and operation. The ulricheasyINJECT CT/MRI injectors must only be operated with original disposables from ulrich medical.

Upon initial filling, the patient tubing is marked as "used" – and from that time cannot be used in any other ulricheasyINJECT CT/MRI injector. In the following cases, the patient tubing is "invalidated" and can no longer be used:

- after an injection volume of 400 ml
- after the maximum use period of 12 hours
- upon activation of the **New Patient** function at the injector

### Scope of delivery and area of application

Art. No.	Designation
XD 10701	Patient tubing, 2.50 m, for ulricheasyINJECT CT/MRI injectors
XD 10702	Patient tubing, 3.20 m, for ulricheasyINJECT CT/MRI injectors

### Product intended for single use

The product is sterilized with ethylene oxide and packed in sterile packaging.

- Do not reuse the product.
- Do not re-sterilize the product.
- Do not use any product from open or damaged sterile packaging.
- Do not use the product after the expiry date.
- Do not use any product that has not been stored under recommended storage conditions.

### Single-use product

-  **WARNING! Infection of patient/operator due to reuse and/or re-processing.** Reuse and/or re-processing cause the risk of contamination and material fatigue. Fouling, contamination, and/or impaired function can lead to injury, illness, or death.  
→ Do not re-process product.

### Safe handling

- Before using the product, check that it is operable and in proper condition.
- In order to avoid damage due to improper use, use the product in accordance with these instructions for use and the injector's instructions for use.
- Product and accessories should only be used by persons who have the required knowledge, training and experience.

### Risk of injury!

The use of tubing combinations or tubing systems not approved by ulrich medical jeopardizes the safety of the patient and/or operator and compromises the proper functioning of the injector.

- Use only original products approved by ulrich medical.

### Risk of injury! Injection of air, overdosing, and flow reduction due to pressure regulation when using tubing extensions.

Tubing extensions jeopardize patient safety and can contaminate the tubing system.

- Do not use tubing extensions (e.g. "Heidelberg extension").

**Risk of injury! Injection of air and spraying of contrast medium. Risk of physical damage.**

Overtightening can break the threads of the Luer lock connector and lead to ruptures.

→ Close the Luer lock connector correctly.

Being closed too tightly and being in prolonged contact with contrast media can make the Luer lock connection hard to unscrew.

→ Carefully unscrew the Luer lock connector.

**Risk of injury due to infections and contamination.**

Improper handling jeopardizes patient safety.

→ Handle all sterile components carefully in order to ensure sterility.

→ Do not use the product if the protective caps are missing or are not firmly seated.

→ Replace the patient tubing for each new patient.

→ After disconnecting the patient tubing, the system must immediately be equipped with new patient tubing, and the protective cap must not be removed from the Luer lock connector until it is connected to the patient.

→ Dispose of the patient tubing after disconnecting it from the patient, in order to prevent contamination of the Easy-Click-Cassette.

## Usage

→ Remove the product from the packaging.

→ Check the product for damaged parts before use.

→ Do not use the product if any parts are damaged.

## Connect the product to the Easy-Click-Cassette

→ Connect the product only as specified in the injector's instructions for use.

## Functional test

→ Perform a functional test as specified in the injector's instructions for use.

## Operation

→ Use the product only as specified in the injector's instructions for use.

## Technical service

ulrich GmbH & Co. KG Injector Systems Customer Service Buchbrunnenweg 12 89081 Ulm / Germany	Tel.: +49 (0)731 9654-0 Fax: +49 (0)731 9654-2809 E-mail: injectors@ulrichmedical.com
---	--

## Technical data

Feature	Value
Pressure tightness	Max. 22 bar
Material	PVC (free from DEHP and latex)
Max. service life	12 hours for one patient (single-use item)

## Storage

→ Store the product in the sterile packaging in a dry, dark and temperature-controlled room.

→ Protect the sterile packaged product from environmental factors such as dirt, dust, humidity, heat, and sunlight.

## Ambient conditions

	Use	Storage
Temperature	+5°C to +40°C	+15°C to +30°C
Relative humidity	10% to 90% non-condensing	30% to 60% non-condensing

## Disposal

→ After use, dispose of the product in accordance with the applicable local regulations.

## Reporting of incidents

→ The user must report all product complaints related to the safety, efficacy or performance of the product to ulrich medical ([complaint@ulrichmedical.com](mailto:complaint@ulrichmedical.com)) or the local distributor within 72 hours.

→ If one or more components of the system show(s) signs of malfunction (i. e. one or more of the performance specifications is/are not met) or do/does not provide the intended performance or if there is a suspicion that this may be the case, then ulrich medical ([complaint@ulrichmedical.com](mailto:complaint@ulrichmedical.com)) or the local distributor should be notified immediately (within 24 hours).

→ If one or more components of the system has/have ever failed and possibly led to or contributed to the death or severe injury of a patient, then ulrich medical or the local distributor should be notified immediately.

→ Please use the Medical Device Vigilance Report Form from ulrich medical to report product complaints, malfunctions or product failure. It is available at: [www.ulrichmedical.com/vigilancereport](http://www.ulrichmedical.com/vigilancereport).

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**ulrich GmbH & Co. KG** | Buchbrunnenweg 12 | 89081 Ulm | Germany