7-AAD Viability Dye

REF A07704 150 tests; 3 mL, 20 µL / test

TABLE OF CONTENTS

English	2
Français (FR)	6
Deutsche (DE)	11
Italiano (IT)	16
Español (ES)	21
Português Portugal (PT-PT)	26
Dansk (DK)	31
Svenska (SE)	35
Norsk (NO)	39
Suomi (FI)	43
Ελληνικά (GR)	47
日本語 (JP)	52
中文 (中国) (ZH-CN)	56
Lietuvių (LT)	60
Magyar (HU)	64
Polski (PL)	69
Čeština (CZ)	74
Slovák (SK)	78
한국어 (KR)	82
Türkçe (TR)	86
Русский (RU)	90
Eestlane (EE)	95
Hrvatski (HR)	99
Български (BG)	103
中文 (台灣) (ZH-TW)	108
Română (RO)	112
Slovenščina (SI)	116
Srpski (Latinski) (SP)	120
Latviešu (LV)	124
Українська (UA)	128
Português Brasil (PT-BR)	133
Nederlands (NL)	138
Tiếng Việt (VN)	143
Қазақша (KZ)	148
APPENDIX	153
REFERENCES	154

	Specifications 7-AAD Viability Dye
Formulation	Liquid
Volume	3 mL
Calibration	Ready for use
λ excitation	488 nM
Emission peak	655 nm (complex with double strand DNA)

7-AAD Viability Dye

REF A07704 150 tests; 3 mL, 20 μ L / test

For *In Vitro* Diagnostic Use

INTENDED USE

7-AAD Viability Dye is a chemical coloring agent in solution, which is designed for mono- or multiparametric analysis using flow cytometry. It permits non-viable cells to be identified and to be excluded from (viable) cells of interest in the analysis of human leucocytes.

PRINCIPLE

7-Amino-Actinomycin D (7-AAD) is an analogue of Actinomycin containing an amine group substituted in position 7 of the chromophore. 7-AAD inserts itself between the tops of the Cytosine/Guanine bases of DNA (1).

Staining of the DNA double strand is undertaken by incubating the sample with the 7-AAD Viability Dye. The red cells are then removed by lysis and the leucocytes are analyzed by flow cytometry.

Apoptotic necrotic and/or damaged cells are a source of interference in the analysis of viable cells using flow cytometry. Non-viable cells can be characterized and identified as they are stained by 7-AAD, whilst living cells, retaining their membranous integrity, are impermeable to 7-AAD and are unstained (i.e., they are 7-AAD negative) (2).

The flow cytometer analyzes light diffusion and the fluorescence of cells. It makes possible the localization of cells within the electronic window defined on a histogram, which correlates the orthogonal diffusion of light (Side Scatter or SS) and the fluorescence corresponding to 7-AAD staining. Other histograms combining two of the different parameters available on the cytometer are also used in the gating stage.

The fluorescence of cells thus gated is analyzed in order to exclude the positively stained events (non-viable) from the unstained ones (viable). The results can be expressed as a percentage of non-fluorescent events in relation to all the events taken into account by the gating.

INTENDED USER

This product is intended for laboratory professional use.

SAMPLES

Venous blood must be taken using sterile tubes containing an EDTA salt as the anticoagulant.

The samples should be kept at room temperature (18–25°C) and not shaken. The sample should be homogenized by gentle agitation prior to taking the test sample.

The samples must be analyzed within 24 hours of venipuncture.

WARNING AND PRECAUTIONS

1. Do not use the reagent beyond the expiry date.
2. Do not freeze.
3. Let it come to room temperature (18–25°C) before use.
4. Minimize exposure to light.
5. Avoid microbial contamination of the reagents, or false results may occur.
6. This ready for use reagent contains 0.005% (weight/volume) of 7-AAD and 1% (volume/volume) of Dimethyl Sulfoxide (DMSO).

In pure form, 7-AAD is potentially carcinogenic and DMSO is an irritant.

Though extremely diluted in the present formulation, these constituents may retain all or part of their harmful effects. Never pipet by mouth, and avoid all contact with the skin, mucosa and eyes. Use these reagents following the precautions for use (notably: the wearing of gloves, gown and protective glasses).

7. 7-AAD is a chromophore subject to degradation by exposure to light. It is stored in an opaque vial to protect it from light prior to use.

Avoid continuous exposure to light during the incubation stages and reduce light exposure of the samples once staining has been undertaken.

8. All blood samples must be considered as potentially infectious and must be handled with care (in particular: the wearing of protective gloves, gowns and goggles).
9. Blood tubes and disposable material used for handling should be disposed of in ad hoc containers intended for incineration.

GHS HAZARD CLASSIFICATION

7-AAD Viability Dye

WARNING

H316

P332+P313

Causes mild skin irritation.

If skin irritation occurs: Get medical advice/attention.

Dimethyl Sulfoxide 0.5 - 1.5%

Safety Data Sheet is available at beckman.com/techdocs

STORAGE AND STABILITY

The 7-AAD Viability Dye must be kept at between 2 and 8°C and protected from light, before and after the vial has been opened.

Closed vial shelf life per stability study: 730 days.

Stability of opened vial: the reagent is stable for 60 days.

See lot specific Certificate of Analysis at www.beckman.com.

EVIDENCE OF DETERIORATION

Any change in the physical appearance of the reagents may indicate deterioration and the reagent should not be used.

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-742-2345 (USA or Canada) or contact your local Beckman Coulter Representative.

MATERIALS NEEDED BUT NOT SUPPLIED WITH KIT:

- Sampling tubes and material necessary for sampling.
- Automatic pipettes with disposable tips for 20, 100 and 500 µL.
- Plastic haemolysis tubes.
- Red cell lysis reagent. For example: IOTest3 Lysing Solution (Ref. A07799).
- Buffer (PBS: 0.01M sodium phosphate; 0.145 M sodium chloride; pH 7.2).
- Centrifuge.
- Automatic agitator (Vortex type).
- Flow cytometer.

PREPARATION OF SAMPLES

The concentration of leucocytes in the sample must be less than 10^4 cells/µL (10^{10} /L). If necessary, dilute in PBS to bring the leucocyte concentration to 5×10^3 cells/µL (5×10^9 /L).

The procedure uses 100 µL of a sample, tube prediluted or otherwise.

NB:

- Adjustment of the intensity of fluorescence corresponding to 7-AAD can be undertaken using the staining of a sample of whole fresh blood and of non-viable stabilized cells in the same tube (see image in appendix). Under these conditions, the voltage corresponding to the detection of 7-AAD must be adjusted such that viable events (i.e., whole fresh blood) appear on the first decade of an SS histogram, depending on the 7-AAD (see example in the appendix).
- Do not use reagents containing permeating or fixating agents during the procedure in order to prevent artifactually positive staining.

PROCEDURE

1. Add 20µL of the 7-AAD Viability Dye solution to each test tube.
2. Add 100µL of the test sample (i.e., the equivalent of approximately 5×10^6 cells).
Vortex the tubes gently.
3. Incubate for 15 to 20 minutes at room temperature (18–25°C), protected from light.

- Then perform lysis of the red cells, if necessary, by following the recommendations of the lysis reagent used.

By way of example, if one wishes to use the IOTest3 Lysing Solution (Ref. A07799) add 2 mL of the working solution (1X), vortex immediately and incubate for 10 minutes at room temperature, protected from light.

If the sample does not contain red cells, do not undertake the lysis stage, but add 0.5 to 1 mL of PBS.

- Preparations should be analyzed within 1 hour.

Note: In all cases, keep the preparations between 2 and 8°C and protected from light.

PERFORMANCE

SPECIFICITY

Actinomycins are the active biological constituents of a chromophore (2-Amino-4,6Dimethylphenoxazone-3) and of cyclic pentapeptides (3). They are antibiotics of bacterial origin characterized historically in soil Ascomycetes. The actinomycins form stable complexes with double strand deoxyribonucleic acid (DNA), but do not form this type of complex either with double strand ribonucleic acid (RNA) or with RNA-DNA hybrids, or with single strand DNA or RNA.

The spectral properties of 7-AAD make it a compound which is particularly well suited to flow cytometry analysis (3). The maximum absorption of the 7-AAD/DNA complex is compatible with the blue excitation wavelength of 488 nm for cytometers fitted with an argon laser (4). The peak fluorescence emission in the deep red band (635 to 675 nm) of the 7-AAD/DNA complex permits optimal use of this probe when it is combined with Fluorescein Isothiocyanate conjugated antibodies (FITC) and with R Phycoerythrin (PE) (3). Indeed, in contrast to Propidium Iodide (PI), which is another fluorescent probe used as a DNA marker, the 7-AAD/DNA complex has a reduced overlapping emission spectrum with FITC and PE.

PRECISION

The percent positive values were determined using Whole Blood spiked with positive cells. Each sample was run 4 times, twice a day for 1 day on 2 instruments using 2 lots of 7-AAD Viability Dye reagent. Measurements (% positive) were made on Navios flow cytometer. Analysis was conducted based on the CLSI method EP5-A2: Evaluation of Precision Performance of Quantitative Measurement Methods.

Our acceptance criteria is dependent on the number of positive events measured for each population:

- If positive event <1500, CV <15%
- If positive event >1500, CV <10%

IOTest3 Lysing Solution:

Whole Blood spiked with HPBALL Cell line							
Number of positive events (Mean) = 5539							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	1.48	1.55	4.44	1.41	2.87	3.05	4.91
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

Versalyse Lysing System:

Whole Blood spiked with HPBALL Cell line							
Number of positive events (Mean) = 4034							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	3.73	1.85	8.74	1.79	3.63	7.02	9.11
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

ACCURACY

The accuracy of 7-AAD Viability Dye was assessed by comparing the results with a reference reagent as the predicate on a set of whole blood samples spiked with positive cells run on a NAVIOS Cytometer. The bias between test and reference reagent was determined based on the difference between test results. If the bias is within the allowable error range or the p-value indicates no significant difference (> 0.05), then test results of the patient samples by the two reagents are considered as equivalent.

The results obtained are summarized in the following table below:

IOTest3 Lysing Solution:

Number of donors = 25				
Positive Target	Mean Δ	Δ % Cell criteria	p-value	RESULTS
Whole Blood spiked with HPBALL Cell line	0.31	<5	0.247	PASS

Versalyse Lysing System:

Number of donors = 25				
Positive Target	Mean Δ	Δ % Cell criteria	p-value	RESULTS
Whole Blood spiked with HPBALL Cell line	-0.34	<5	0.289	PASS

LIMITATIONS

1. Flow cytometry may produce false results if the cytometer has not been aligned perfectly, if fluorescence leaks have not been correctly compensated for and if the regions have not been carefully positioned.
2. Accurate and reproducible results will be obtained as long as the procedures used are in accordance with the technical insert leaflet and compatible with good laboratory practices.
3. 7-AAD, the active substance in this reagent, has been optimized so as to offer the best specific/non-specific signal ratio. Therefore, it is important to adhere to the reagent volume/sample volume ratio in every test.
4. In the case of a hyperleucocytosis, dilute the blood in PBS so as to obtain a value of approximately 5×10^9 leucocytes/L.
5. In certain disease states, such as severe renal failure or haemoglobinopathies, lysis of red cells may be slow, incomplete or even impossible. In this case, it is advisable to isolate mononucleated cells using a density gradient (Ficoll, for example) prior to staining.

See the Appendix for examples and references.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

ADDITIONAL INFORMATION

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on In vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

REVISION HISTORY

REVISION AF:	Release date : October 2021
REVISION AW:	Release date : April 2022
Updates to comply with the Beckman Coulter Global labelling Policy and according to IVD-R (EU)2017/746 requirements:	
Added sections	Intended User, Concentration, Precision, Accuracy, Additional Information, Revision History
Added information	See sections Principle
Updated sections	Principle, Warning and Precautions, GHS Hazard Classification, Storage and Stability, Evidence of Deterioration, Limitations, Appendix
Removed sections	Methodology, Examples of Clinical Application, Expected Results, Intra-Laboratory Reproducibility, Linearity
REVISION	Release Date
AX	October 2022
Updated sections	Added Translations
REVISION	Release Date
AY	
Updated sections	Storage and stability

Symbols Key

Glossary of Symbols is available at beckman.com/techdocs (document number B60062)

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
APPENDIX

- EXAMPLES

The graph below is a biparametric representation (Fluorescence Intensity vs. Side Scatter) of a lyzed fresh whole blood sample from an healthy donor spiked with Stem-Trol Control Cells (Ref. IM3632). Since Stem-Trol Control Cells are stabilized cells, they are not viable and are stained with the the 7-AAD Viability Dye (visible in the right part of the histogram). In contrast, a majority of the fresh whole blood cells are not stained (visible in the left part of the histogram).

Acquisition is performed with a Beckman Coulter Navios flow cytometer equipped with the Navios analysis software.

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Тел. +7 (495) 228 67 00, e-mail: beckman.ru@beckman.com

IMMUNOTECH SAS A Beckman Coulter Company, 130, Avenue de Lattre de Tassigny,

BP 177, 13276 Marseille Cedex 9, France, +(33) 4 91 17 27 27

www.beckman.com