IOTest

Conjugated Antibody

CD79a-APC

REF B36287 100 tests; 1 mL, 10 µL / test

TABLE OF CONTENTS

English	2
Français (FR)	7
Deutsche (DE)	12
Italiano (IT)	17
Español (ES)	22
Português Portugal (PT-PT)	27
Dansk (DK)	32
Svenska (SE)	37
Norsk (NO)	42
Suomi (FI)	47
Ελληνικά (GR)	52
日本語 (JP)	57
中文 (中国) (ZH-CN)	62
Lietuvių (LT)	67
Magyar (HU)	72
Polski (PL)	77
Čeština (CZ)	82
Slovák (SK)	87
한국어 (KR)	92
Türkçe (TR)	97
Русский (RU)	102
Eestlane (EE)	107
Hrvatski (HR)	112
Български (BG)	117
中文 (台灣) (ZH-TW)	122
Română (RO)	127
Slovenščina (SI)	132
Srpski (Latinški) (SP)	137
Latviešu (LV)	142
Українська (UA)	147
Português Brasil (PT-BR)	152
Nederlands (NL)	157
Tiếng Việt (VN)	162
Қазақша (KZ)	167
APPENDIX	172
REFERENCES	173

	Specifications
Specificity	CD79a
Clone	HM47
Hybridoma	NS1 x balb/c
Immunogen	Synthetic peptide (amino acids 202 - 216) of human mb-1 protein containing the cyto-plasmic portion
Immunoglobulin	IgG1
Species	Mouse
Purification	Affinity Chromatography
Fluorochrome	Allophycocyanin (APC)
Molar ratio	APC / Ig: 0.5 - 1.5
λ excitation	633/638 nm
Emission Peak	660 nm
Buffer	PBS pH 7.2 plus 2 mg / mL BSA and 0.1% NaN ₃

IOTest Conjugated Antibody CD79a-APC

REF B36287 100 tests; 1 mL, 10 µL / test

For *In Vitro* Diagnostic Use

INTENDED USE

This fluorochrome-conjugated antibody allows the qualitative and non automated identification of cell populations expressing the CD79a antigen present in human biological samples using flow cytometry (see section "Samples" below).

PRINCIPLE

This test is based on the ability of specific monoclonal antibodies to bind to the antigenic determinants expressed by leucocytes.

Permeability is induced in the cytoplasmic membranes of leucocytes for the demonstration of intracellular antigenic determinants by means of monoclonal fluorescent antibodies. The leucocytes are then analyzed by flow cytometry.

The flow cytometer measures light diffusion and the fluorescence of cells. It makes possible the delimitation of the population of interest within the electronic window defined on a histogram, which correlates the orthogonal diffusion of light (Side Scatter or SS) and the diffusion of narrow-angle light (Forward Scatter or FS). Other histograms combining two of the different parameters available on the cytometer can be used as supports in the gating stage depending on the application chosen by the user.

The fluorescence of the delimited cells is analyzed in order to distinguish the positively-stained events from the unstained ones. The results are expressed as a percentage of positive events in relation to all the events acquired by the gating.

INTENDED USER

This product is intended for laboratory professional use.

CLINICAL RELEVANCE

The CD79a-APC is a CD79a antibody used to identify and characterize cells expressing the CD79a antigen by flow cytometry. This product alone cannot and is not intended to generate any diagnostic conclusion.

When used in combination with other markers, this product can be used in one or more of the following functions:

- To aid in the differential diagnosis of hematologically abnormal patients suspected of having hematopoietic neoplasm and to monitor patients with known hematopoietic neoplasm.

See following references :

1; 2; 3; 4; 5; 6; 7; 8; 9; 10

SAMPLES

Venous blood must be taken using sterile tubes containing an EDTA salt as the anticoagulant.

The samples should be kept at room temperature (18 – 25°C) and not shaken. The samples should be homogenized by gentle agitation prior to taking the test sample.

The samples must be analyzed within 24 hours of venipuncture.

CONCENTRATION

See lot specific Certificate of Analysis at www.beckman.com.

WARNING AND PRECAUTIONS

1. Do not use the reagent beyond the expiry date.
2. Do not freeze.
3. Let it come to room temperature (18 – 25°C) before use.
4. Minimize exposure to light.
5. Avoid microbial contamination of the reagents, or false results may occur.

6. Antibody solutions containing sodium azide (NaN_3) should be handled with care. Do not take internally and avoid all contact with the skin, mucosa and eyes.
Moreover, in an acid medium, sodium azide can form the potentially dangerous hydrazoic acid. If it needs to be disposed of, it is recommended that the reagent be diluted in a large volume of water before pouring it into the drainage system so as to avoid the accumulation of sodium azide in metal pipes and to prevent the risk of explosion.
7. All blood samples must be considered as potentially infectious and must be handled with care (in particular: the wearing of protective gloves, gowns and goggles).
8. Never pipet by mouth, and avoid all contact of the samples with the skin, mucosa and eyes.
9. Blood tubes and disposable material used for handling should be disposed of in ad hoc containers intended for incineration.
10. Reagents and waste should be eliminated according to local requirements.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

Safety Data Sheet is available at beckman.com/techdocs

STORAGE AND STABILITY

This reagent must be kept at between 2 and 8°C and protected from light, before and after the vial has been opened.

Closed vial shelf life per stability study: 730 days.

Stability of opened vial: the reagent is stable for 180 days.

EVIDENCE OF DETERIORATION

Any change in the physical appearance of the reagents may indicate deterioration and the reagent should not be used.

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-742-2345 (USA or Canada) or contact your local Beckman Coulter Representative.

CONTENTS

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

MATERIALS NEEDED BUT NOT SUPPLIED WITH KIT:

- Sampling tubes and material necessary for sampling.
- Automatic pipettes with disposable tips for 10, 100 and 500 μL .
- Plastic haemolysis tubes.
- To obtain optimal results, either following reagents are recommended (follow related specific procedure):
 - IntraPrep Fixation/Permeabilization reagent (Ref A07802 or A07803)
 - PerFix-nc Kit (no centrifuge assay kit), for Intra-& Extra-Cellular Staining Preparation (Ref B31167 or B31168)
- Leucocyte fixation reagent. For example : IOTest 3 Fixative Solution (Ref. A07800).
- Isotype control APC : IOTest reagent (Ref. IM2475).
- Buffer (PBS: 0.01 M sodium phosphate; 0.145 M sodium chloride; pH 7.2).
- Centrifuge.
- Automatic agitator (Vortex type).
- Flow cytometer.

PROCEDURE WITH PERFIX-NC REAGENT

Below is the recommended procedure to use with the PerFix-nc Kit (no centrifuge assay kit) for Intra-& Extra-Cellular Staining Preparation (Ref B31167 or B31168).

For each sample analyzed, in addition to the test tube, one control tube can be added in which the cells are mixed in the presence of the isotype control (Ref. IM2475).

Kit reagent preparation

PerFix-nc Buffer 1 and 2 Reagents

No reconstitution is necessary. Both reagents may be used directly from the vial.

PerFix-nc Buffer 3 Reagent preparation

Prepare extemporaneously the Final 1X Reagent.

Dilute the 10X Concentrated PerFix-nc Buffer 3 (Final 10X Solution) into deionized water: 1 volume of Buffer 3 with 9 volumes of water. Mix well before use. We recommend preparing only the volume of Final 1X Reagent necessary for the experiments of the day.

1. Pipet 50 μL of blood sample into the bottom of each appropriately labeled tube. Avoid putting some blood on the side of the tube; otherwise it will not be appropriately treated.
2. Pipet 5 or 25 μL of the Fixative Reagent to each tube depending if a low or high fixation is required (refer to PerFix-nc instruction for use for details on low/high fixation protocols).

3. Vortex immediately and incubate for 15 min at room temperature (18 – 25°C).
4. Vortex again the fixed blood and add 300 µL of the Permeabilizing Reagent to each tube; vortex immediately.
5. Add immediately to each tube 10 µL of the fluorochrome-conjugated antibodies against intracellular epitopes and surface molecules (alternatively, the antibodies can be pre-mixed into the Permeabilizing Reagent and added altogether at the end of the fixation step).
6. Vortex immediately and incubate for 15 - 30 min at room temperature.
7. Add 3 mL of the Final 1X Reagent (prepared from the 10X concentrated Final Solution) to each tube; vortex immediately; the sample is now ready for analysis on a flow cytometer.

NOTE: The preparations can be stored for 24 hours prior to cytometric analysis, it is advisable to store them at 2-8°C and protected from light.

PROCEDURE WITH INTRAPREP REAGENT

Below the recommended procedure to use with IntraPrep Fixation/Permeabilization reagent (Ref. A07802 or A07803).

For each sample analyzed, in addition to the test tube, one control tube can be added in which the cells are mixed in the presence of the isotype control (Ref. IM2475).

For a blood sample, optimal staining is obtained using a number of leucocytes between 3 and 10 x 10³ cells / µL. If the leucocyte concentration is greater than 10 x 10³ cells / µL, dilute (11).

1. Add 50 µL of blood sampled into EDTA to each tube.
2. Add to each tube 100 µL of IntraPrep reagent 1 (Fixation)
3. Vortex the tubes vigorously for 3 to 5 seconds.
4. Incubate for 15 minutes at room temperature (18 – 25°C) protected from light.
5. Add 4 mL of PBS to each tube.
6. Centrifuge for 5 minutes at 300 x g at room temperature. Remove the supernatant by aspiration.
7. Add to each tube 100 µL of IntraPrep reagent 2 (Permeabilization). Let mix by diffusion. DO NOT VORTEX.
8. Incubate for 5 minutes at room temperature without shaking.
9. Shake the tubes carefully and manually for 2 to 3 seconds.
10. Add 10 µL of specific IOTest conjugated antibody to each test tube and if necessary, 10 µL of the isotype control to each control tube.
11. Incubate for 15 to 20 minutes at room temperature (18 – 25°C), protected from light.
12. Add 4 mL of PBS and centrifuge at 300 x g for 5 minutes at room temperature.
13. Remove the supernatant by aspiration and resuspend the cell pellet in 0.5 to 1 mL of IOTest 3 Fixative Solution (Ref. A07800) at its working concentration (1X).
14. The preparations are ready for cytometric analysis.

NOTE: If the preparations are to be stored for more than 2 hours prior to cytometric analysis, it is advisable to store them at 2-8°C and protected from light. The preparations thus stored do not keep, however, for more than 24 hours.

EXPECTED VALUES

In our laboratories, the whole blood samples of 25 apparently healthy donors were treated using the reagent described above. The results obtained for the count of the positive events of interest with this reagent are given in the table below :

PerFix-nc Lysing System:

Positive Target	Number	Mean (%)	SD	CV (%)
Lymphocytes	25	10.98	4.13	37.60

IntraPrep Lysing System:

Positive Target	Number	Mean (%)	SD	CV (%)
Lymphocytes	25	10.02	4.03	40.19

These values are intended to be representative only. Each laboratory should establish its own expected values from the local population of normal donors.

PERFORMANCE

Performance data are obtained using the procedure described above on less than 24 hour-old blood samples previously collected in sterile tubes with EDTA salt as anticoagulant. Analysis is performed within 2 hours following immunostaining.

SPECIFICITY

The CD79a molecule is part of the CD79a / CD79b disulphide-linked heterodimer, non-covalently bound to surface immunoglobulins to form B cell receptors (BCR) (12). The expression of CD79a appears early in the ontogeny of B cells and its localization at the pro-B stage is therefore cytoplasmic. Later on, the CD79a forms part of the BCR. Its membrane expression persists up to the plasmocytic stage, the stage at which its localization once again becomes cytoplasmic (13).

MAb HM47 reacts with an intracytoplasmic epitope of the CD79a molecule (13). It was assigned to CD79a at the 5th HLDA Workshop on Human Leucocyte Differentiation Antigens held in Boston, USA, in 1993 (WS Code: cb017, Section B) (13).

PRECISION

The percent positive values were determined using Whole Blood. Each sample was run 4 times, twice a day for 1 day on 2 instruments using 2 lots of CD79a-APC Monoclonal Antibody Reagents. Measurements (% positive) were made on Navios flow cytometer. Analysis was conducted based on the CLSI method EP5-A2: Evaluation of Precision Performance of Quantitative Measurement Methods.

Our acceptance criteria is dependent on the number of positive events measured for each population:

- If positive event <1500, CV <15%
- If positive event >1500, CV <10%

PerFix-nc Lysing System:

Lymphocytes							
Number of positive events (Mean) = 1026							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	1.79	0.95	3.88	1.03	1.73	2.97	4.13
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

IntraPrep Lysing System:

Lymphocytes							
Number of positive events (Mean) = 880							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	1.24	2.36	4.79	1.22	2.36	3.98	5.48
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

ACCURACY

The accuracy of CD79a-APC was assessed by comparing the results with a reference reagent as the predicate on a set of whole blood samples run on a Navios flow cytometer. The bias between test and reference reagent was determined based on the difference between test results. If the bias is within the allowable error range or the p-value indicates no significant difference (> 0.05), then test results for the two reagents are considered as equivalent.

The results obtained are summarized in the following table below:

PerFix-nc Lysing System:

Number of donors = 25				
Positive Target	Mean Δ	Δ % Cell criteria	p-value	RESULTS
Lymphocytes	0.32	<3	0.005	PASS

IntraPrep Lysing System:

Number of donors = 25				
Positive Target	Mean Δ	Δ % Cell criteria	p-value	RESULTS
Lymphocytes	0.09	<3	0.659	PASS

LIMIT OF BLANK AND LIMIT OF DETECTION

A study was conducted in accordance with CLSI EP17-A2, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures. The Limit of Detection (LOD) is the lowest concentration of analyte that can be consistently detected. Results obtained are summarized in the following table:

PerFix-nc Lysing System:

Positive Target	Limit of Blank (cell/ μ L)	Limit of Detection (cell/ μ L)
Lymphocytes	0	4

IntraPrep Lysing System:

Positive Target	Limit of Blank (cell/ μ L)	Limit of Detection (cell/ μ L)
Lymphocytes	1	2

LIMITATIONS

1. Flow cytometry may produce false results if the cytometer has not been aligned perfectly, if fluorescence leaks have not been correctly compensated for and if the regions have not been carefully positioned.
2. Accurate and reproducible results will be obtained as long as the procedures used are in accordance with the technical insert leaflet and compatible with good laboratory practices.
3. The conjugated antibody of this reagent is calibrated so as to offer the best specific signal/non-specific signal ratio. Therefore, it is important to adhere to the reagent volume/sample volume ratio in every test.
4. In the case of a hyperleucocytosis, dilute the blood in PBS so as to obtain a value of approximately 5×10^9 leucocytes/L (11).
5. In certain disease states, such as severe renal failure or haemoglobinopathies, lysis of red cells may be slow, incomplete or even impossible. In this case, it is recommended to isolate mononucleated cells using a density gradient (Ficoll, for example) prior to staining (14).
6. In patients treated with anti-human monoclonal antibody therapies, detection of the specific targeted antigens may be diminished or absent due to partial or complete blocking by the therapeutic antibody.
7. The CD79a-APC results should be interpreted in light of the total clinical presentation of the patient, including: symptoms, clinical history, data from additional tests, and other appropriate information.

See the Appendix for examples and references.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

ADDITIONAL INFORMATION

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on In vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

The Summary of Safety and Performance is available from the EUDAMED database: ec.europa.eu/tools/eudamed.

REVISION HISTORY

REVISION AC:	Release date : October 2019
REVISION	Release Date
AW	February 2022
Updates to comply with the Beckman Coulter Global labelling Policy and according to IVD-R (EU)2017/746 requirements:	
Added sections	BSI 2797 Number, Intended User, Clinical Relevance, Concentration, Precision, Accuracy, Limit of blank and limit of detection, Additional Information, Revision History.
Added information	See sections Limitations
Phrasing or typographic updates	See sections Procedure, Performance, Limitations, Warning and Precautions, Storage and Stability.
Removed sections	Example of clinical applications, Reagents, Intra-laboratory reproducibility, Linearity
Updated sections	Intended Use, GHS Hazard Classification, Evidence of deterioration, Procedure, Appendix.
REVISION	Release Date
AX	
Updated sections	Add Kazakh
Updated sections	Storage and stability

Symbols Key

Glossary of Symbols is available at beckman.com/techdocs (document number B60062)

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
APPENDIX

- EXAMPLES


The graph below is a biparametric representation (Side Scatter versus Fluorescence Intensity) of a permeabilized / lysed normal whole blood sample. Whole blood samples stained with IOTest CD79a-APC Conjugated Antibody.

Acquisition is performed with a Beckman Coulter Navios flow cytometer equipped with the Navios analysis software.

PerFix-nc Lysing System

IntraPrep Lysing System

REFERENCES

1. N. Singh, H. P. P. Pati, S. Tyagi, R. Deka, "Proposed minimal panel of antibodies for cost-effectiveness and accuracy in acute leukemias immunophenotyping: Prospective study at a tertiary care center.", *Hematology*. 21, 338–342 (2016).
2. R. Paredes-Aguilera, L. Romero-Guzman, N. Lopez-Santiago, L. Burbano-Ceron, O. Camacho-Del Monte, S. Nieto-Martinez, "Flow cytometric analysis of cell-surface and intracellular antigens in the diagnosis of acute leukemia.", *Am J Hematol.* 68, 69–74 (2001).
3. G. Vasconcelos de Andrade Alves, A. L. Araújo da Cunha Fernandes, J. M. Freire, A. de Souza Paiva, R. C. de Vasconcelos, V. Soraya de Farias Sales, T. Maria de Araújo Moura Lemos, E. A. Gil, F. H. Miranda de Araújo Freire, D. G. K. C. e Silva, J. F. R. Maciel, I. G. Farkatt, G. B. C. Júnior, "Flow Cytometry Immunophenotyping Evaluation in Acute Lymphoblastic Leukemia: Correlation to Factors Affecting Clinic Outcome." *J Clin Lab Anal.* 26, 431–440 (2012).
4. S. Naeem, M. H. Bukhari, "Antigen Expression on Blast Cells and Hematological Parameters at Presentation in Acute Lymphoblastic Leukemia Patients." *J Coll Physicians Surg Pak.* 25, 407–411 (2015).
5. H.-F. Li, W.-T. Meng, Y.-Q. Jia, N.-G. Jiang, T.-T. Zeng, Y.-M. Jin, Q.-R. Huang, X. Li, H. Xu, X.-M. Mo, "Development-associated immunophenotypes reveal the heterogeneous and individualized early responses of adult B-acute lymphoblastic leukemia.", *Medicine (Baltimore)*. 2016 Aug; 95(34): e4128.
6. J.-H. Yu, J.-T. Dong, Y.-Q. Jia, N.-G. Jiang, "Individualized leukemia cell-population profiles in common B-cell acute lymphoblastic leukemia patients.", *Ai zheng.* 32, 213–223 (2013).
7. E. S. Blix, J. M. Irish, A. Husebekk, J. Delabie, L. Forfang, A. M. Tierens, J. H. Myklebust, A. Kolstad, "Phospho-specific flow cytometry identifies aberrant signaling in indolent B-cell lymphoma.", *BMC Cancer.* 12, 478 (2012).
8. N. Jiang, C. Chen, Q. Gong, K. Shields, Y. Li, Y. Chen, J. Song, T. W. McKeithan, W. C. Chan, "Flow cytometric sorting coupled with exon capture sequencing identifies somatic mutations in archival lymphoma tissues.", *Lab Invest.* 97, 1364–1374 (2017).
9. V. Deneys, A. M. Mazzon, J. L. Marques, H. Benoit, M. De Bruyère, "Reference values for peripheral blood B-lymphocyte subpopulations: a basis for multiparametric immunophenotyping of abnormal lymphocytes.", *Journal of Immunological Methods.* 253, 23–36 (2001).
10. T. M. Owaidah, A. A. Beihany, M. A. Iqbal, N. Elkum, G. T. Roberts, Cytogenetics, "Molecular and ultrastructural characteristics of biphenotypic acute leukemia identified by the EGIL scoring system.", *Leukemia.* 20, 620–626 (2006).
11. H42-A2 Vol.27 No. 16. P30. Enumeration of Immunologically Defined Cell Populations by Flow Cytometry; Approved Guideline - Second Edition.
12. Reth, M., Hombach, J., Wienands, J., Campbell, K.S., Chien, N., Justement, L.B., Cambier, J.C., "The B-cell antigen receptor complex", 1991, *Immunol. Today*, 6, 12, 196-200.
13. Engel, P., Wagner, N., Tedder, T.F., "CD79 workshop report", 1995, *Leucocyte Typing V, White Cell Differentiation Antigens*. Schlossman, S.F., et al., Eds., Oxford University Press, 667-670.
14. Ibrahim FF, Ghannam MM, Ali FM., "Effect of dialysis on erythrocyte membrane of chronically hemodialyzed patients"., 2002. Ren Fail., Nov;24(6):779-90.

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IMMUNOTECH SAS A Beckman Coulter Company, 130, Avenue de Lattre de Tassigny,

BP 177, 13276 Marseille Cedex 9, France, +(33) 4 91 17 27 27
www.beckman.com