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α: Ortho-Clinical Diagnostics; Vitros Systems
 β: Abbott; AxSYM, Architect,
 γ: Roche Diagnostics; Roche/Hitachi, Elecsys®, Modular
 γγ: Roche Diagnostics; Cobas INTEGRA®
 δ: Beckmann-Coulter; Synchron LX/CX, Immage/Array, Access
 ε: Siemens Healthcare Diagnostics; Dimension®, BN-Systems, Stratus CS
 Ω: Beckmann Coulter, Olympus-Analysers
 κ: Siemens Healthcare Diagnostics; Immulite
 λ: Bio-Rad
 μ: Siemens Healthcare Diagnostics; ADVIA Centaur /ACS 180
 σ: Siemens Healthcare Diagnostics, Enzygnost

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7 Bibliografia

1. Adcock D, Kressin D, Martar RA. Effects of 3.2 % vs. 3.8 % citrate concentration on routine coagulation testing. *Am J Clin Pathol* 1997; 107: 105-10.
2. Adcock D, Kressin D, Marlar RA. The effect of time and temperature variables on routine coagulation tests. *Blood Coagul Fibrinolysis* 1998; 9: 463-70.
3. Agnese ST, Spierito FW, Hannon WH. Evaluation of four reagents for delipidation of serum. *Clin Biochem* 1983; 16: 98-100.
4. Alber B, Mühlbayer D. Wichtige Einflussgrößen und Störfaktoren für die Beurteilung von Tumormarkerkonzentrationen. *Der Bay Int* 2004; 24: 136-40.
5. Altura BT, Shirey TL, Young CC, Dell'Orfano K, Hiti J, Welsh R, Yeh Q, Barbour RL, Altura B. Characterization of a new ion selective electrode for ionized magnesium in whole blood, plasma, serum, and aqueous samples. *Scand J Clin Lab Invest* 1994; Suppl 217: 21-36.
6. Andersson A, Lindgren A, Hultberg B. Effect of thiol oxidation and thiol export from erythrocytes on determination of redox status of homocysteine and other thiols in plasma from healthy subjects and patients with cerebral infarction. *Clin Chem* 1995; 41: 361-6.
7. Annesley TM, Hunter BC, Fidler DR, Giachero DA. Stability of tacrolimus (FK 506) and cyclosporin G in whole blood. *Therap Drug Monit* 1995; 17: 361-5.
8. Artiss JD, McEnroe RJ, Zak B. Bilirubin interference in a peroxidase-coupled procedure for creatinine eliminated by bilirubin oxidase. *Clin Chem* 1984; 30: 1389-92.
9. Artiss JD, Zak B. Problems with measurements caused by high concentrations of serum lipids. *CRC Crit Rev Clin Lab Sci* 1987; 25: 19-41.
10. Astles R, Williams CP, Sedor F. Stability of plasma lactate in vitro in the presence of antiglycolytic agents. *Clin Chem* 1994; 40: 1327-30.
11. Aufenanger J, Zawia B. Preanalytical aspects of lipoprotein measurement. *Clin Lab* 1999; 45: 535-46.
12. Balistreri WF, Shaw LM. Liver function. In: Tietz NW, ed. *Fundamentals of Clinical Chemistry*. Philadelphia: Saunders, 1987 3rd ed., 729-67.
13. Banfi G, Salvagno GL, Lippi G. The role of ethylene-diamine tetraacetic acid (EDTA) as in vitro anticoagulant for diagnostic purposes. *Clin Chem Lab Med* 2007; 45: 565-76.
14. Banks RE. Measurement of cytokines in clinical samples using immunoassays: Problems and pitfalls. *Crit Rev Clin Lab Sci* 2000; 37: 131-82.
15. Bargnoux AS, Dupuy AM, Biron-Andréani C, Schved JF, Cristol JP. Immuno-nephelometric determination of fibrinogen on citrated or heparinized plasma: comparison with functional Claus method. *Clin Lab* 2005; 51: 285-8.
16. Bauer K. Determination of free haemoglobin in serum by automated assay using 4-aminophenazone and the Cobas Bio system. *Clin Chem Clin Biochem* 1981; 19: 971-6.
17. Baum H. Kardiale Troponine - Bedeutung der Präanalytik und Analytik für die Diagnostik von Herzmuskelschädigungen. *Der Bay Int* 2004; 24: 141-6.
18. Baum H, Bookslegers P, Steinbeck G, Neumeier D. A rapid assay for the quantification of myoglobin: evaluation and diagnostic relevance in the diagnosis of acute myocardial infarction. *Eur J Clin Chem Clin Biochem* 1994; 32: 853-8.
19. Beaulieu M, Lapointe Y, Vinet B. Stability of pO₂, pCO₂ and pH in fresh blood samples stored in a plastic syringe with low heparin in relation to various blood-gas and hematological parameters. *Clin Biochem* 1999; 32: 101-7.
20. Becker B, Denzler B, Kolde HJ, Ramirez I, Rombach B, Strauß J. Fehlermöglichkeiten bei Hämostasestests, Teil 1: Präanalytik, Reagenzien. *MTA* 1994; 9: 794-800.
21. Ben Rayana MC, Burnett RW, Covington, AK, Orazio PD, Fog Anderson N, Jacobs E, et al. IFCC guidelines for sampling, measuring and reporting ionized magnesium in plasma. *Clin Chem Lab Med* 2008; 46: 21-6.
22. Bidart JM, Thuillier F, Augereau C, Chalas J, Daver A, Jacob N, et al. Kinetics of serum tumor marker concentrations and usefulness in clinical monitoring. *Clin Chem* 1999; 45: 1695-707.
23. Blank DW, Kroll MH, Ruddel ME, Elin RJ. Hemoglobin interference from in vivo hemolysis. *Clin Chem* 1985; 31: 1566-9.
24. Boink ABTJ, Buckley BM, Christiansen TF, Covington AK, Maas AHJ, Müller-Plathe O, et al. IFCC-Recommendations on sampling, transport and storage for the determination of concentration of ionized calcium in whole blood, plasma and serum. *Eur J Clin Chem Clin Biochem* 1991; 29: 767-72.
25. Bonini PA, Plebani M, Ceriotti F, Rubboli F. Errors in laboratory medicine. *Clin Chem* 2002; 48: 691-8.
26. Boomsma F, Alberts G, van Eijk L, Man in't Feld, Schalekamp MADH. Optimal collection and storage conditions for catecholamine measurement in human plasma and urine. *Clin Chem* 1993; 39: 2503-8.
27. Boyanton BL, Blick K. Stability studies of twenty-four analytes in human plasma and serum. *Clin Chem* 2002; 48: 2242-7.
28. Brady J, O'Leary N. Interference due to lipaemia in routine photometric analysis - survey of an underrated problem. *Ann Clin Biochem* 1994; 31: 281-8.
29. Burnett RW, Covington AK, Fogh-Andersen N, Külpmann WR, Maas AHJ, Müller-Plathe O, et al. Approved IFCC recommendations on whole blood sampling, transport, and storage for simultaneous determination of pH, blood gases, and electrolytes. *Eur J Clin Chem Clin Biochem* 1995; 33: 247-53.
30. Bush V, Blennerhasset J, Wells A, Dasgupta A. Stability of therapeutic drugs in serum collected in vacutainer serum separator tubes containing a new gel (SST II). *Ther Drug Monit* 2001; 23: 259-62.
31. Caliezi C, Reber G, Lämmle B, de Moerloose P, Wuillemin WA. Agreement of D-dimer results measured by a rapid ELISA (VIDAS) before and after storage during 24 h or transportation of the original whole blood samples. *Thromb Haemost* 2000; 83: 177-8.
32. Callas DD, Clark TL, Moreira PL, Lansden C, Gawryl MS, Kahn S, Bermes EW. In vitro effects of a novel haemoglobin-based oxygen carrier on routine chemistry, therapeutic drug, coagulation, hematology, and blood bank assays. *Clin Chem* 1997; 43: 1741-8.
33. Carilli I, Kraajenhagen R, Pradella R, D'Okotiro G, Herkner K, Ravan RM et al. In vitro stability of the reticulocyte count. *Clin Lab Haematol*. 1996; 18: 9-11.
34. Cartledge JJ, Thompson D, Verril H, Clarkson P, Eardley I.

- The stability of free and bound prostate-specific antigen. *Brit J Urol* 1999; 84: 810-4.
35. Checkonin VP, Zhirkov YA, Belyaeva A, Ryabakhin IA, Gurina OI, Dimitriyeva TB. Serum time course of two brain specific proteins, alpha (1) brain globulin and neuronspecific enolase, in tickborn encephalitis and Lyme disease. *Clin Chim Acta* 2002; 320: 117-25.
36. Chetty M. The stability of anticonvulsant drugs in whole blood. *Ther Drug Monit* 1994; 16: 491-4.
37. Chiu RWK, Poon LL, Lan TK, Leung TN, Poon PM et al. Quantitative analysis of fetal and total DNA quantitation in material plasma. *Clin Chem* 2001; 47: 1607-13.
38. CLSI, Document H21-A5. Collection, transport and processing of blood specimens for testing plasma-based coagulation assays. Approved standard, 5th ed. Wayne, PA, 2007.
39. CLSI, Document EP7-P. Interference testing in clinical chemistry; proposed guideline. Wayne, PA, 1986.
40. CLSI, Document H3-A6. Procedures for the collection of diagnostic blood specimens by venipuncture. Approved standard, 6th ed. Wayne, PA, 2007.
41. Cobbaert C, Tricarica A. Different effects of IntralipidTM and triacylglycerol rich lipoproteins on Kodak Ektachem serum cholesterol determination. *Eur J Clin Chem Clin Biochem* 1993; 31: 107-9.
42. Colombo JP (ed.). *Klinisch-chemische Urindiagnostik*. Rotkreuz: Labolife, 1994.
43. Conrad K. Autoimmunerkrankungen. In: Guder WG, Nolte J (Hrsg) *Das Laborbuch für Klinik und Praxis*. 2. Aufl München: Elsevier Urban und Fischer 2009.
44. Cooper RG, Sampson EJ, Smith SJ. Preanalytical, including biological variation in lipid and apolipoprotein measurements. *Curr Opin Lipidol* 1992; 3: 365-71.
45. Court JJ, Clark PM, Holder G, Sheppard MC. Collection and storage conditions of serum for the analysis of thyroglobulin. *Proc ACB National Meeting. Ann Clin Biochem* 2002, Suppl: 36-7.
46. Dale JC, Pruef SK. Phlebotomy - a minimalist approach. *Mayo Clin Proc* 1993; 68: 249-55.
47. Dasgupta A, Yared MA, Wells A. Time-dependent absorption of therapeutic drugs by the gel of Greiner Vacuette blood collection tube. *Ther Drug Monit* 2000; 22: 427-31.
48. De Jongh R, Vranken J, Vundelinckx, Bosmans E, Maes M, Heylen R. The effects of anticoagulation and processing on assays of IL-6, sIL-6R, sIL-2R and soluble transferrin receptor. *Cytokine* 1997; 9: 696-701.
49. Delanghe JR, Chapelle JP, Vanderschueren SC. Quantitative nephelometric assay for determining myoglobin evaluated. *Clin Chem* 1990; 36: 1675-8.
50. Develler M, Blanchaert N, Komárek A, Bossuyt X. Can heparin plasma be used instead of serum for nephelometric analysis of serum proteins. *Clin Chem*. 2006; 52: 1609-10.
51. Diver MJ, Hughes JG, Hutton JL, West CR, Hipkin LJ. The long-term stability in whole blood of 14 commonly-requested hormone analytes. *Ann Clin Biochem* 1994; 31: 561-5.
52. Dörner K, Böhler H. *Labor diagnostische Strategien in der Pädiatrie*. Darmstadt: GIT, 1996.
53. Dugan S, Bogema S, Schwartz RW, Lappas NT. Stability of drugs of abuse in urine samples stored at -20 °C. *J Anal Toxicol* 1994; 18: 391-6.
54. Dugué B, Leppänen E, Gräsbeck R. Preanalytical factors and the management of cytokines in human subjects. *Int J Clin Lab Res* 1996; 26: 99-105.
55. Duncanson GO, Worth HJG. Pseudohypophosphataemia as a result of bilirubin interference. *Ann Clin Biochem* 1990; 27: 263-7.
56. Durham BH, Robinson J, Fraser WD. Differences in the stability of intact osteocalcin in serum, lithium heparin plasma and EDTA plasma. *Ann Clin Biochem* 1995; 32: 422-3.
57. Einer G, Zawta B. *Präanalytik-Fibel*. Leipzig, Heidelberg: Barth, 1991, 2nd ed.
58. Engstadt CS, Guttenberg TJ, Osterer B. Modulation of blood cell activation by four commonly used anticoagulants. *Thromb Hemost* 1997; 77: 690-6.
59. Ernst DJ, Szamosi DI. Specimen collection standards complete major revisions. *Med Lab Observ* 2005; 37: 26-9.
60. Erwa W, Bauer FR, Etschwaiger R, Steiner V, Scott CS, Sedlmayr P. Analysis of aged samples with the Abbott CD 400 hematology analyzer. *Eur J Lab Med* 1998; 6: 4-15.
61. Erickson JA, Ashwood EA. Evaluation of a myelin basic protein (MBP) assay and the identification of highly elevated MBP in demyelinating disorders. *Clin Chem*. 2003; 43: A108 (abstract).
62. European IVD Directive for in vitro diagnostics. *Arbeitsblatt der Europäischen Gemeinschaften L331/1 vom 7.12.1998*.
63. Evans K, Mitcheson J, Laker MF. Effect of storage at 4 °C and -20 °C on lipid, lipoprotein, and apoprotein concentrations. *Clin Chem* 1995; 41: 392-6.
64. Evans MJ, Livesey JH, Ellis MJ, Yandle TG. Effect of anticoagulants and storage temperatures on stability of plasma and serum hormones. *Clin Biochem* 2000; 34: 107-12.
65. Evans WE, Oellerich M, Holf DW. *Drug monitoring, Leitfaden für die klinische Praxis*, 2nd ed. Wiesbaden: Abbott Diagnostics Division, 1994.
66. Faynor SM, Rosbinson R. Suitability of plastic collection tubes for cyclosporin measurements. *Clin Chem* 1998; 44: 2220-1.
67. Fiedler GM, Thiery J. Der "fehlerhafte" Laborbefund: 1: Fehlerquellen der prä- und postanalytischen Phase. *Internist* 2004; 45: 315-32.
68. Finney H, Newman DJ, Gruber W, Merle P, Price CP. Initial evaluation of cystatin C measurement by particle-enhanced immunonephelometry of the Behring nephelometer systems (BNA, BNII). *Clin Chem* 1997; 43: 1016-22.
69. Flower L, Ahuja RH, Humphries SE, Mohamed-Ali V. Effects of sample handling on the stability interleukine 6, tumor necrosis factor alpha and leptin. *Cytokine* 2000; 12: 1712-6.
70. Fonseca-Wollheim F da. Ringstudie zur Bilirubininterferenz bei der Bestimmung von Serum-Kreatinin. *Laboratoriumsmedizin* 1988; 12: 317-20.
71. Fonseca-Wollheim F da. Ultrafiltrate analysis confirms the specificity of the selected method for plasma ammonia determination. *Eur J Clin Chem Clin Biochem* 1992; 30: 15 9.
72. Fonseca-Wollheim F da. Deamidation of glutamine by increased plasma α -glutamyltransferase is a source of rapid ammonia formation in blood and plasma specimens. *Clin Chem* 1990; 36: 1479-82.
73. Fonseca-Wollheim F da, Heinze KG, Lomsky K, Schreiner H. Serum ultrafiltration for the elimination of endogenous interfering substances in creatinine determination. *J Clin Chem Clin Biochem* 1988; 26: 523-5.

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amplification by heparin. *J Clin Microbiol* 1991; 29: 676-9.

113. Holodny M, Rainer L, Herman S, Yen-Lieberman B. Stability of plasma immunodeficiency virus load in vacutainer PPT plasma preparation tubes during overnight shipment. *J Clin Microbiol* 2000; 38: 323-6.

114. Huizenga JR, van der Belt K, Gips CH. The effect of storage at different temperatures on cholinesterase activity in human serum. *J Clin Chem Clin Biochem* 1985; 23: 283-5.

115. International council for standardization in haematology. Recommendations of the ICSH for ethylenediamine tetraacetic acid anticoagulation of blood for blood cell counting and syring. Expert panel on cytometry *Am J Clin Pathol* 1993; 100: 371-2.

116. Ishii H, Sarai K, Sanemori H, Kawasaki T. Analysis of thiamine and its phosphate esters by high-performance liquid chromatography. *Anal Biochem* 1979; 97: 191.

117. ISO/EN/DIN 6710. Single-use containers for human venous blood specimen collection. Geneva/Bruelxelles/Berlin, 2007.

118. ISO/EN/DIN 15189. Medical laboratories - Particular requirements for quality and competence Geneva/Bruelxelles/Berlin, 2007.

119. ISO GUIDE 30. Terms and definitions used in connection with reference materials, 2nd ed. Geneva 1992.

120. Johnston A, Cullen G, Holt DW. Quality assurance for cyclosporin assays in body fluids. *Ann Acad Med Singapore* 1991; 20: 3-8.

121. Jung K, Klinggraef P von, Brux B, Sinha P, Schnor L, Loening SA. Preanalytical determinations of total and free prostate-specific antigen and their ratio: blood collection and storage conditions. *Clin Chem* 1998; 44: 685-8.

122. Jung R, Lübcke C, Wogener C, Neumaier M. Reversal of rt-PCR inhibition observed in heparinized clinical specimens. *Biotechniques* 1998; 23: 24-8.

123. Kallner A. Preanalytical procedures in the measurement of ionized calcium in serum and plasma. *Eur J Clin Chem Clin Biochem* 1996; 34: 53-8.

124. Kampa JS, Keffer P. Evaluation of the B-type natriuretic peptide (BNP) and pro-brain natriuretic peptide (NT-pro-BNP) in congestive heart failure (CHF) patients. *Clin Chem* 2003; 49: A 32 (abstract).

125. Kazmierczak SC, Catrou PG, Best AE, Sullivan SW, Briley KP. Multiple regression analysis of interference effects from a haemoglobin-based-oxygen carrier solution. *Clin Chem Lab Med* 1999; 37: 453-64.

126. Keiichi M, Suday A, Kawasaki J. Effect of urine pH, storage time and temperature on stability of catecholamines, cortisol, and creatinine. *Clin Chem* 1998; 44: 1759-62.

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129. Kendall RG, Chapman C, Hartley AE, Norfolk DR. Storage and preparation of samples for erythropoietin radioimmunoassay. *Clin Lab Haematol* 1991; 13: 189-96.

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133. Koch TR, Platoff G. Suitability of collection tubes with separation gels for therapeutic drug monitoring. *Ther Drug Monit* 1990; 12: 277-80.

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136. Koivula T, Groenroos P, Gövert J, Icen A, Irtala K, Penttilä I, et al. Basic urinalysis and urine culture: Finnish recommendations from the Working Group on clean midstream specimens. *Scand J Clin Lab Invest* 1990; 50 (suppl. 200): 26-33.

137. Korte W, Riesen WF. Comparability of serum and plasma concentrations of haemostasis. *Clin Chem Lab Med* 2001; 39: 627-30.

138. Kouri T, Fogazzi G, Gant V, Hallander H, Hofmann W, Guder WG. European urinalysis guidelines. *Scand J Clin Lab Invest* 2000; 60 (suppl 231): 1-96.

139. Kramer KA, O'Brien JF, McConnel JP. Oxidized LDL: Evaluation of an ELISA method for potential use as a marker of cardiovascular risk. *Clin Chem* 2003; 49: A168 (abstract).

140. Kreutzer HJH, Paanakker MPWM. Stability of routine chemistry parameters in lithium heparin gel tubes. *Proc 11th Eur Congr Clin Chem, Tampere* 1995, abstr. 357.

141. Külpmann WR. Determination of electrolytes in serum and serum water. *Wiener Klin Wschr* 1992; suppl: 34-8.

142. Kuhl L, Tamm A. On the stability of vitamin B12, folic acid and ferritin in serum. *Quality of Preanalytical Phase in Europe*. Leuven: Terumo, 2001.

143. Lammers M. Dilution of citrated plasma. *Eur J Clin Chem Clin Biochem* 1996; 34: 369.

144. Lammers M, Gressner AM. Immunonephelometric quantification of free haemoglobin. *J Clin Chem Clin Biochem* 1987; 25: 363-7.

145. Lammers M, Rausch S. Suitability of plasma specimen for immunonephelometric protein assays. *Clin Chem Lab Med* 2002; 40: S114.

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147. Leary NO, Pembroke A, Duggan PF. Measuring albumin and calcium in serum in a dual test with the Hitachi 704. *Clin Chem* 1992; 38: 1342-5.

148. Ledue TB, Collins MF, Craig WY, Ritchie RF. Effect of storage time, temperature and preservative on IgG, transferrin, albumin and alpha-1-microglobulin levels in urine. *Clin Chem* 2000; 46, Suppl.: A 48 (abstract).

149. Lee JS, Carlton E, Ossolinska-Plewnia J, Lu J. Invasive trophoblast antigen (ITA) testing: sample collection and handling. *Clin Chem* 2003; 49:A127 (abstract).

150. Lefevre GF. Analytical performance of the ACB test for ischemia on the Konelab20. *Clin Chem* 2003; 49: A33 (ab-

- stract).
151. Leino A, Koivula MK. Stability of chemical and immunochemical analyses in uncentrifuged plasma samples. *Ann Clin Biochem* 2009; 46: 159-61.
152. Leinonen J, Stenman UH. Reduced stability of prostate-specific antigen after long-term storage of serum at -20°C . *Tumour Biol* 2000; 21: 46-53.
153. Leonard PJ, Persaud J, Motwani R. The estimation of plasma albumin by BCG dye binding on the Technicon SMA 12/60. *Clin Chim Acta* 1971; 35: 409.
154. Leppänen E, Gräsbeck R. The effect of the order of filling tubes after venipuncture on serum potassium, total protein, and aspartate and alanine aminotransferase. *Scand J Clin Lab Invest* 1986; 46: 189-91.
155. Levine B, Smith ML. Stability of drugs of abuse in biological specimens. *Forensic Sci Rev* 1990; 2: 148-56.
156. Lew J, Reichelderfer P, Fowler M, Bremer J, Carrol R, Casols S, et al. Determination of levels of human immunodeficiency virus Type 1 RNA in plasma: reassessment of parameters affecting assay outcome. *J Clin Microbiol* 1998; 38: 1471-9.
157. Lippi G, Brocco G, Salvagno GL, Montagnana M, Guidi GC, Schmidt Gayk H. Influence of sample matrix on the stability of B-CTX at room temperature for 24 and 48 hours. *Clin Lab* 2007; 53:455-9.
158. Lippi G, Giampolo L, Giancresce G. Effect of anticoagulants on lipoprotein(a) measurements with four commercial assays. *Eur J Clin Chem Clin Biochem* 1996; 34: 251-5.
159. Lippi G, Salvagno GL, Solero GP, Franchini M, Guidi GC. Stability of blood cell counts, hematological parameters and reticulocyte indexes on the ADVIA 120 hematological analyzer. *J Lab Clin Med* 2005; 166: 333-40.
160. Lippi U, Schinella M, Nicoli M, Modena N, Lippi G. EDTA induced platelet aggregation can be avoided by a new anticoagulant also suitable for automated blood count. *Haematologia* 1990; 75: 38-41.
161. Lorentz K. Approved Recommendation on IFCC methods for the measurements of catalytic concentration of enzymes. Part 9. IFCC method for α -amylase (1,4- α -D-glucan-4-glucanohydrolase, EC 3.2.1.1). *Clin Chem Lab Med* 1998; 36: 185-203.
162. Luddington R, Peters J, Baker P, Baglin R. The effect of delayed analysis or freeze-thawing on the measurement of natural anticoagulants, resistance to activated protein C and markers of activation of the haemostatic system. *Thromb Res* 1997; 87: 577-81.
163. Lutomski DM, Bower RH. The effect of thrombocytosis on serum potassium and phosphorus concentration. *Am J Med Sci* 1994; 307: 255-8.
164. Maguire GA, Price CP. A continuous monitoring spectrophotometric assay for the measurement of angiotensin-converting enzyme in human serum. *Ann Clin Biochem* 1985; 22: 204.
165. Mair J, Puschendorf B. Aktuelle Aspekte der Labordiagnostik des akuten Myokardinfarktes. *J Lab Med* 1995; 19: 304-18.
166. Matteucci E, Giampetro O. To store urinary enzymes: How and how long? *Kidney Int* 1994; Suppl 47: 58-9.
167. McKenzie D, Henderson A. Electrophoresis of lactate dehydrogenase isoenzymes. *Clin Chem* 1983; 29: 189-95.
168. Meisner M, Tschakowsky K, Schnabel S, Schmidt BJ, Kataliivic A, Schüttler J. Procalcitonin - influence of temperature, storage, anticoagulation and arterial or venous asseration of blood samples on procalcitonin concentrations. *Eur J Clin Chem Clin Biochem* 1997; 35: 597-601.
169. Menssen HD, Melber K, Brandt N, Thiel E. The use of hirudin as universal anticoagulant in haematology, clinical chemistry and blood grouping. *Clin Chem Lab Med* 2001; 39: 1267-77.
170. Merve D, Henley R, Lane G, Field R, Frenneaux M, Dunstan F, McDowell I. Effect of different sample types and stability after blood collection of N-terminal proB type natriuretic peptide as measured with the Roche Elecsys system. *Clin Chem* 2004; 50: 779-80.
171. Meyer T, Monge PK, Sakshaug J. Storage of blood samples containing alcohol. *Acta Pharmacol Toxicol* 1979; 45: 282-6.
172. Morello R, Boos KS, Seidel D. Standardized preanalytics of urinary marker molecules and on line SPE-LC-ECD analysis of catecholamines in urine and plasma. 2nd Pharmaceutical Sciences World Congress (PSWC) 2004, Kyoto, Japan.
173. Morrè SA, Van Valkengoed IGM, De Jong A, Boeke AJP, Van Eijk JTM, Meijer CJLM, et al. Mailed, home-obtained urine specimen: a reliable screening approach for detecting asymptomatic Chlamydia trachomatis infections. *J Clin Microbiol* 1999; 37: 976-80.
174. Mößmer G. Präanalytik von hämosteologischen Laboruntersuchungen. *Der Bay Int* 2004; 24: 124-132.
175. Müller T, Gegenhuber A, Dieplinger B, Poelz W, Haltmayer M. Long term stability of endogeneous B type natriuretic peptide (BNP) and amino terminal pro BNP (NT-proBNP) in frozen plasma samples. *Clin Chem Lab Med* 2004; 42: 942-4.
176. Mussap M, Ruzzante N, Varagnolo M, Plebani M. Quantitative automated particle-enhanced immunonephelometric assay for the routine measurement of human cystatin C. *Clin Chem Lab Med* 1998; 36: 859-65.
177. Nanji AA, Poon R, Hineberg I. Lipaemic interference: Effects of lipaemic serum and Intralipid. *J Clin Pathol* 1998; 41: 1026-7.
178. Narayanan S. Protection of peptidic substrates by protease inhibitors. *Biochim Clin* 1987; 11: 954-6.
179. Narayanan S. Quality control in tumor marker analysis: preanalytical issues. *J Clin Ligand Assay* 1998; 21: 11-7.
180. Narayanan S. The preanalytical phase. An important component of laboratory medicine. *Am J Clin Pathol* 2000; 113: 429-52.
181. Neumaier M, Braun A, Wagener C. Fundamentals of quality assessment of molecular amplification methods in clinical diagnostics. *Clin Chem* 1998; 44: 12-26.
182. Nilson E, Ekholm B, ReesSmith B, Torn C, Hillburan M. Calcium addition to EDTA-Plasma eliminates falsely positive results in the RSR GADA by ELISA. *Clin Chim Acta* 2008; 388: 130-4.
183. Nilsson L, Hedner U, Nilsson M, Robertson B. Shelf-life of bank blood and stored plasma with special reference to coagulation factors. *Transfusion* 1983; 23: 377-81.
184. Numata Y, Doki K, Furukawa Akikuoka S, Asada H, Fukonaga T, et al. Immunoradiometry assay for the N-terminal fragment of proatrial natriuretic peptide in human plasma. *Clin Chem* 1998; 44: 1008-13.
185. Okabe R, Nakatsuka K, Imada M, Miki T, Naka H, Nasaki H, et al. Clinical evaluation of the Elecsys β -crosslabs serum assay, a new assay for degradation products of type I collagen C-telopeptides. *Clin Chem* 2001; 47: 1410-4.
186. Ono T, Kitaguchi K, Takehara M, Shiba M, Hayami K. Serum constituents analysis: effect of duration and temperature of storage of clotted blood. *Clin Chem* 1981; 27: 35-8.

187. Oremek GM, Seiffert UB. Physical activity releases prostate-specific antigen (PSA) from the prostate gland into blood and increases serum PSA concentrations. *Clin Chem* 1996; 42: 691-5.
188. Ossendorf M, Fichtner J, Schroeder S, Thueroff JW, Prellwitz W. In vitro Stabilität des freien prostataspezifischen Antigens (PSA) in Serum- und Vollblutproben. *Anticancer Res* 1997; 17: 4199, No 6C, Abstr. No 58.
189. Pai JK, Curhan GC, Cannuscio CC, Rifai N, Ridker PM, Rimm RB. Stability of novel plasma markers associated with cardiovascular disease: processing within 36 hours of specimen collection. *Clin Chem* 2002; 48: 1781-84.
190. Panteghini M, Pagoni F. Pre-analytical, analytical and biological sources of variation of lipoprotein (a). *Eur J Clin Chem Clin Biochem* 1993; 31: 23-8.
191. Pararo C, Tagini E, Meda R, Scott CS, Novaro O. Improved stability of leucocytes in aged samples: investigation of an alternative anticoagulant strategy in hematology for use with the Abbott Cell-Dyn CD 3500 hematology analyzer. *Eur J Lab Med* 1998; 6: 16-23.
192. Pfund A, Wendland G, Geisen C, Hoepf HW. Reliability of homocysteine measurement. *Herz/Kreislauf* 1999; 33: 381-5.
193. Pfützer A, Pfützer AH, Kann PH, Stute R, Löbig M, Yang JM, Mistry J, Forst T. Clinical and laboratory evaluation of a new specific ELISA for intact proinsulin. *Clin Lab* 2005; 51: 243-9.
194. Probst R, Brandl R, Blümke M, Neumeier D. Stabilization of homocysteine concentration in whole blood. *Clin Chem* 1998; 44: 1567-9.
195. Rainen L, Oelmueller U, Jurgensen S, Wyrich R, Ballas C, Schram J, Herchma C et al. Stabilization of mRNA expression in whole blood samples. *Clin Chem* 2002; 48: 1883-90.
196. Rahmati K, Lernmark A, Becker C, Foltyn-Zadura A, Larsson K, Ivarson S-A, Törn C. A comparison of serum and EDTA-plasma in the measurement of glutamic acid decarboxylase autoantibodies (GADA) and autoantibodies to islet Antigen 2 (IA-2A) using the RSR radioimmunoassay (RIA) and enzyme linked immunosorbent assay (ELISA) kits. *Clin Chem* 2008; 54: 227-35.
197. Ramont L, Thoannes H, Volondat A, Chastang F, Milef M-C, Marquart F-X. Effects of hemolysis and storage conditions on neuron-specific enolase (NSE) in cerebrospinal fluid and serum, in patients in clinical practice. *Clin Chem Lab Med* 2005; 43: 1215-7.
198. Rao LV, Okorodudu AO, Petersen JR, Elghetany MT. Stability of prothrombin time and activated partial prothrombin time tests under different storage conditions. *Clin Chim Acta* 2000; 300: 13-21.
199. Rasbold K, Rendell MS, Goljan E. Simple removal of lipids from serum. *Clin Chem* 1985; 31: 782.
200. Rasmussen K, Møller J, Lyngbak M, Holm Pedersen AM, Dybkaer L. Age- and gender-specific reference intervals for total homocysteine and methylmalonic acid in plasma before and after vitamin supplementation. *Clin Chem* 1996; 42: 630-6.
201. Reisch N, Mauracher B, Bidlingmaier M, Reinke M. ACTH stability dependent on different preanalytical procedures. *Exp. Clin Endocrinol Diabetes* 2006; 114: S 76-7.
202. Reiter W, Göhring P, Stieber P, Pahl H, Kempfer B, Banauch D. The influence of blood taking, storage conditions and interfering factors on a new fibrin monomer test. *J Lab Med* 1996; 20: 112-6.
203. Richardson TD, Wojno KJ, Liang LW, Giachero DA, England BG, Henricks WH, et al. Half-life determination of serum free prostate-specific antigen following radical retropubic prostatectomy. *Urology* 1996; 48: 40-4
204. Richtlinie der Bundesärztekammer zur Qualitätssicherung quantitativer Laboratoriumsmedizinischer Untersuchungen. *Dt Ärztebl* 2008; 105: C301-15.
205. Ricos C, Alvarez V, Cava F, Garcia-Lario JV, Hernandez A, Jimenez CV, et al. Current databases on biological variation: pros, cons and progress. *Scand J Clin Lab Invest* 1999; 59: 491-500.
206. Robinson N, Mangin P, Sangi M. Time and temperature dependent changes in red blood cell analytes used for testing recombinant erythropoietin abuse in sports. *Clin Lab* 2004; 50: 317-23.
207. Rodríguez-Mendizabal M, Lucena MI, Cabello MR, Blanco E, Lopez-Rodríguez B, Sanchez de la Cuesta F. Variations in blood levels of aminoglycosides related to in vitro anticoagulant usage. *Ther Drug Monit* 1998; 20: 88-91.
208. Rosello M, Campos F, Berlanga E, Gavarro A, Acosta A, Torra M. Plasma citrate. An alternative sample for total homocysteine determination. *Clin Chem Lab Med* 2003; 41: S444 (abstract).
209. Roß RS, Eller T, Volbracht L, Paar D. Interferenzen durch Lipämie, Hämolyse und Hyperbilirubinämie am DAX 48-Analysator und ihre klinische Relevanz. *Lab Med* 1994; 18: 233-9.
210. Roß RS, Paar D. Analytisch und klinisch relevante Interferenzen in der Gerinnungsanalytik am Beispiel des MDA 180. *J Lab Med* 1998; 22: 90-4.
211. Rothe G. „Immunphänotypisierung von Leukozyten und Lymphomen“ und „Thrombozytenfunktion (durchflusszytometrisch)“ in Guder WG, Nolte J. *Das Laborbuch für Klinik und Praxis*. 2. Auflage, München: Elsevier, Urban und Fischer, 2009
212. Russell D, Henley R. The stability of parathyroid hormones in blood and serum samples at 4 °C and at room temperature. *Ann Clin Biochem* 1995; 32: 216-7.
213. Sachse C, Baudach A, Tille D, Avenarius HJ, Heller S, Ruby C, et al. Multicentre evaluation of the Cobas Argos blood cell counter. *J Lab Med* 1994; 18: 441-9.
214. Sakurai S, Shiojima I, Tanigawa T, Nakahara K. Aminoglycosides prevent and dissociate the aggregation of platelets in patients with EDTA-dependent pseudothrombocytopenia. *Br J Haem* 1997; 99: 817-23.
215. Sanders GTB. Sample quality. *Klin Chem Mitt* 1992; 23: 11-8.
216. Schiele F, Vincent-Viry M, Herbeth B, Visvikis A, Siest G. Effect of short- or long-term storage on human serum and recombinant apolipoprotein E concentration. *Clin Chem Lab Med* 2000; 38: 525-8.
217. Schmitt UM, Stieber P, Pahl H, Reinmiedl J, Fateh-Moghadam A. Stabilität tumorassoziierter Antigene in Vollblut und Serum. *J Lab Med* 2000; 24: 475.
218. Schmitt Y. Influence of preanalytical factors on the atomic absorption spectrometry determination of trace elements in biological samples. *J Trace Elem Electrolytes Health Dis* 1987; 1: 107-14.
219. Schmitz G, Rothe G, Ruf A, Barlage S, Tschöpe D, Clemenson KJ et al. European working group on clinical cell analysis: Consensus protocol for the flow cytometric characterization of platelet function. *Thrombos Hemostas* 1998; 79: 885-96.
220. Schoenmakers CHH, Kuller L, Lindemans J, Blijenberg BG. Automated enzymatic methods for creatinine measurement

- with special attention to bilirubin interference. *Eur J Clin Chem Clin Biochem* 1993; 31: 861-8.
221. Scholer A, von Rickenbach R, Faffa G, Vetter B. Evaluation of a microparticle enzyme immunoassay for the measurement of sirolimus in whole blood (Abbot IMX Sirolimus). *Clin Lab* 2006; 52: 325-34.
222. Schreiber W, Rausch S, Lammers M. Choice of specimen for immunonephelometric protein assays. *Clin Chim Acta* 2005; 355 suppl. 407 (abstract).
223. Schwinger R, Antoni DH, Guder WG. Simultaneous determination of magnesium and potassium in lymphocytes, erythrocytes and thrombocytes. *J Trace Elem Electrolytes Health Dis* 1987; 1: 88-98.
224. Seitz G, Stickel F, Fiehn W, Werte E, Simanowski UA, Seitz HK. Kohlenhydrat-defizientes Transferrin. Ein neuer, hochspezifischer Marker für chronischen Alkoholkonsum. *Dt Med Wschr* 1995; 120: 391-5.
225. Semjonow A. PSA - unverzichtbar in der Urologie. *Diagnostica Dialog* 1994, 3: 3-4.
226. Serio F, Ruggieri V, Varraso L, DeSario R, Mastrorilli A, Pansini N. Analytical evaluation of the Dade Behring Dimension RXL N-terminal proBNP (NT-proBNP) method and comparison with the Roche Elecsys 2010. *Clin Chem Lab Med* 2005; 43: 1263-73.
227. Sharma A, Anderson K, Baker JW. Flocculation of serum lipoproteins with cyclodextrin: Application to assay of hyperlipidemic serum. *Clin Chem* 1990; 36: 529-32.
228. Shih J, Datwyler SA, Hsu SC, Matias MS, Pacenti DP, Lueders C et al. Effect of collection tube type and preanalytical handling on myeloperoxidase concentrations. *Clin Chem* 2008; 54: 1076-9.
229. Siest G, Bertrand P, Quin B, Herbeth B, Serot JM, Masana L, et al. Apolipoprotein E polymorphism and serum concentration in Alzheimer's disease in nine European centres: the ApoEurope study. *Clin Chem Lab Med* 2000; 38: 721-30.
230. Simó JM, Camps J, Vilella E, Gómez F, Paul A, Joven J. Instability of lipoprotein(a) in plasma stored at -70 °C: Effects of concentration, apolipoprotein(a) genotype, and donor cardiovascular disease. *Clin Chem* 2001; 47: 1673-8.
231. Skopp G, Klingmann A, Pötsch L, Mattern. In vitro stability of cocaine in whole blood and plasma including ecgonine as a target analyte. *Ther Drug Monit* 2001; 23: 174-81.
232. Skopp G, Pötsch L, Klingmann A, Mattern R. Stability of morphine, morphine-3-glucuronide, and morphine-6-glucuronide in fresh blood and plasma and post mortem blood samples. *J Anal Toxicol* 2001; 25: 2-7.
233. Smith AC, Wu AH. Analytical performance evaluation of the Bayer Advia Centaur BNP immunoassay. *Clin Chem* 2003; 49: A71 (abstract).
234. Smith LJ, Long J, Matters DJ, Carr-Smith HD, Bradwell AR. Sample storage and stability for free light chain assays. *Clin Chem* 2003; 49: A05 (abstract).
235. Smoller BA, Kruskall, MS. Phlebotomy for diagnostic laboratory tests in adults: pattern of use and effect on transfusion requirements. *N Engl J Med* 1986; 314: 1233-5.
236. Sokoll LJ, Baum H, Collison PO, Gurr E, Haass M, Luthe H et al. Multicenter analytical performance evaluation of the elecsys proBNP method. *Clin Chem Lab Med* 2004; 42: 965-72.
237. Sölétormos G, Schiöler V, Nielsen D, Skovsgaard T, Dombrowsky P. Interpretation of results for tumor markers on the basis of analytical imprecision and biological variation. *Clin Chem* 1993; 39: 2077-83.
238. Sonntag O. Haemolysis as an interference factor in clinical chemistry. *J Clin Chem Clin Biochem* 1986; 24: 127-39.
239. Sonntag O, Römer M, Haeckel R. Interferences. In Haeckel R (ed.) *Evaluation Methods in Laboratory Medicine*. Weinheim VCH 1993: 101-16.
240. Sonntag O, Glick MR. Serum-Index und Interferogramm – Ein neuer Weg zur Prüfung und Darstellung von Interferenzen durch Serumchromogene. *Laboratoriumsmedizin* 1989; 13: 77-82.
241. Spain MA, Wu AHB. Bilirubin interference with determination of uric acid, cholesterol, and triglycerides in commercial peroxidase-coupled assays. *Clin Chem* 1986; 32: 518-21.
242. Spence N. Differential leucocyte analysis of samples up to 48h old: improved results with citrate pyridoxalphosphate anticoagulant. *Br J Biochem Sci* 1993; 20: 645-63.
243. Stanger O, Herrmann W, Pietrzik K, Fowler B, Geisel J, Di-erkes J, et al. Consensus paper on the rational clinical use of homocysteine, folic acid and B-vitamins in cardiovascular and thromboembolic diseases-guidelines and recommendations. www.dach-liga-homocystein.org
244. Steimer W. Besondere Bedeutung der Präanalytik und Interpretation bei der Bestimmung von Arzneimittelkonzentrationen. *Der Bay Int* 2004; 24: 147-57.
245. Steinbach G, Rau B, Debard A-L, Javarez J-F, Bienvenu J, Ponzio A et al. Multicenter evaluation of a new immunoassay for procalcitonin measurement on the Kryptor®-system. *Clin Chem Lab Med* 2004; 42: 440-9.
246. Stieber P. „Malignes Wachstum“ in Guder WG, Nolte J (Hrsg.) *Das Laborbuch für Klinik und Praxis*. 2. Aufl. München: Elsevier, Urban und Fischer, 2009.
247. Stiegler H, Fischer Y, Vazquez-Jimenez JF, Graf J, Filzmaier K, Fausten B, et al. Lower cardiac troponin T and I results in heparin-plasma than serum. *Clin Chem* 2000; 46: 1338-44.
248. Szasz G, Gerhardt W, Gruber W, Bernt E. Creatinekinase in serum: 2. Interference of adenylate kinase with the assay. *Clin Chem* 1976; 22: 1806-11.
249. Tahara Y, Shima K. Kinetics of HbA1c, glycated albumin, and fructosamine and analysis of their weight functions against preceding plasma glucose level. *Diabetes Care* 1995; 18: 440-7.
250. Tencer J, Thyssel H, Andersson K, Grubb A. Stability of albumin, protein HC, immunoglobulin G, κ -, λ -chain immunoreactivity, orosomuroid and 1-antitrypsin in urine stored at various conditions. *Scand J Clin Lab Invest* 1994; 54: 199-206.
251. Tencer J, Thyssel H, Andersson K, Grubb a. Long-term stability of albumin, protein HC, immunoglobulin G, kappa- and lambda-chain immunoreactivity, orosomuroid and alpha-1-antitrypsin in urine stored at -20 °C. *Scand J Urol Nephrol* 1997; 31: 67-72.
252. Thiery J, Fiedler GM. Der "fehlerhafte" Befund. 2: Häufige Ursachen von Fehlinterpretationen labormedizinischer Befunde. *Internist* 2004; 45: 437-54.
253. Thomas L (ed.). *Clinical Laboratory Diagnosis*. Frankfurt: TH-Books, 1998; Labor und Diagnose, 7th ed. Frankfurt TH-Books, 2008.
254. Tietz NW. *Clinical Guide to Laboratory Tests*, 3rd ed. Philadelphia: Saunders, 1995.
255. Tietz NW. *Textbook of Clinical Chemistry*. Philadelphia: Saunders 1986, 1534-6.

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BD Diagnostics
Preanalytical Systems

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© Niemieckie Towarzystwo Chemii Klinicznej i Medycyny Laboratoryjnej