| INDEX, ÍNDICE, ΕΥΡΕΤΗΡΙΟ | |
|--------------------------|--|
| 1-14 | ENGLISH |
| 15 -29 | ESPAÑOL |
| 30 -44 | ITALIANO |
| 45 -59 | DEUTSCHE |
| 60 - 74 | ΕΛΛΗΝΙΚΟ |
| 75 - 89 | FRANÇAIS |
| 90 -104 | POLSKI |
| 105 -106 | SYMBOLS SIMBOLOS SIMBOLI SYMBOL SYMBOLIS ΣΥΜΒΟΛΟ SIMBOLIS SYMBOLY |

Instructions for Use

Thermo Scientific Sensititre YeastOne Susceptibility Plates

| |
|--------------------------------|
| 029-YEAST – ROW-IVD CID9708 |
| Revision Date: June 12th, 2017 |

THERMO SCIENTIFIC SENSITITRE YEASTONE SUSCEPTIBILITY PLATES

For *in vitro* Diagnostic Use

For more information including plate layout and QC ranges, please refer to www.trekds.com/techinfo. Plate code and batch number required.

Intended Use

The Thermo Scientific™ Sensititre™ YeastOne™ susceptibility plate is an *in vitro* diagnostic microbroth dilution method in a dried 96-well microplate format that provides quantitative minimum inhibitory concentration (MIC) results for non-fastidious yeast including *Candida* species, *Cryptococcus* species, *Aspergillus* species and miscellaneous other rapid-growing yeast species.

Principles of Use

Sensititre YeastOne is a colorimetric microbroth dilution test. Each microplate is dosed with antifungal agents at appropriate dilutions and a colorimetric indicator.

Results are manually read by observing the lowest antifungal concentration that exhibits no growth (as evidenced by no color change).

Precautions

Results should be used as an aid in selecting the drug of choice for treatment.

Only trained personnel should operate the system.

Proper organism handling and disposal methods should be used.

Storage and Shelf Life

Plates should be stored at room temperature (15-25°C) away from direct sunlight and direct heat. Each plate is packaged in foil with a silica gel desiccant. Do not use the plate or broth if (1) the desiccant color is not orange, (2) the expiration date has passed, or (3) the foil pouch is damaged.

Materials

Included:

Sensititre YeastOne plate

Adhesive seal

Not Included [Product Code]:

Thermo Scientific™ Sensititre™ Demineralised Water [T3339]

Sensititre YeastOne Broth [Y3462]

Thermo Scientific™ Sensititre™ Doseheads (for use with Thermo Scientific™ Sensititre™ AIM™ Automated Inoculation Delivery System) [E3010]
Sensititre AIM System [V3020]
Thermo Scientific™ Sensititre™ Vizion™ Digital MIC Viewing System [V2021]
Thermo Scientific™ Sensititre™ Nephelometer [V3011]
Manual Viewer [V4007]
0.5 McFarland Turbidity Standard [E1041]
Bacteriological Loop
20µl Pipette
Sterile Inoculum Reservoir
100ul Pipette and Disposable Tips

Quality Control Strains

Fungal Growth Agar Plates e.g Sabauroud Dextrose Agar (SDA)
Incubator 34 - 36°C, non-CO2
Vortex Mixer

Additional Materials Required for *Aspergillus* spp. Testing

Takashio Agar or Potato Dextrose Agar
Cotton Swab
Sterile Saline with 0.05% Tween™ 20

Specimen Collection and Preparation

Specimens should be collected, transported, stored and plated onto primary isolation medium using standard procedures¹.

Selection of Susceptibility Test Broth

Sensititre System approved broths are performance tested for use with the Sensititre System.

Inoculation Procedure (*Candida* and *Cryptococcus* spp. Testing)

Allow all broths to come to room temperature before use.

Plates should be inoculated within 5 hours of removal from the pouch. A final organism density of approximately $1.5 - 8 \times 10^3$ CFU/ml is recommended.

Steps 1 and 2 should be completed within 15 minutes.

1. Pick several well isolated colonies of >1mm diameter from a pure 24-hour culture (Sabauroud Dextrose Agar) of the yeast isolate and emulsify into demineralised water. Mix well ensuring the suspension is uniform, vortex if required. If clumping occurs, allow

the suspension to settle before adjusting the density. Adjust to a 0.5 McFarland standard visually or with a Sensititre Nephelometer

2. Transfer 20 µl of the suspension into 11 ml of Sensititre YeastOne Broth to give a final inoculum of $1.5 - 8 \times 10^3$ CFU/ml.
3. Transfer 100µl of the final suspension to the Sensititre YeastOne Plate within 15 minutes of completing step 2 by either:
 - a. **Sensititre AIM System** - Replace the tube cap with a Sensititre Single-use Dosehead and inoculate the plate according to the Sensititre AIM user manual. Remove the test tube dosehead combination from the Sensititre AIM System within 30 seconds of dosing a plate and store inverted in a rack or discard
 - b. **Manual pipette** – Pour the broth into a sterile seed trough and inoculate the plate using an appropriate pipette
4. A check of the colony count should be done by removing 10µl from the positive control well and plating onto Sabauroud Dextrose Agar (SDA). A correct inoculum will produce 10-80 colonies
5. Cover all wells with the adhesive seal. Avoid creases prevent well skips

Inoculation Procedure (*Aspergillus* spp. Testing)^{5-7, 22}

Allow all broths to come to room temperature before use.

Plates should be inoculated within 5 hours of removal from the pouch.

A final organism density of approximately $0.5-5 \times 10^4$ CFU/ml is recommended.

1. Subculture from Sabouraud Dextrose Agar onto or Potato Dextrose Agar.
2. Incubate for 7 days at 35°C to obtain adequate sporulation.
3. Collect conidia with a cotton swab and suspend in sterile saline with Tween™. Allow heavy particles to settle for 3 to 5 minutes.
4. Collect supernatant and mix with a vortex mixer
5. Adjust turbidity of supernatant to 80 to 82% transmittance at 530nm measured with a spectrophotometer equivalent to an inoculum of $0.6 - 5 \times 10^6$ cfu/ml. Alternatively, adjust to a 0.5 McFarland Standard.
6. Add 100µl of the suspension to 11ml of Sensititre YeastOne broth to give a final inoculum of $0.5-5 \times 10^4$ cfu/ml.
7. Transfer 100µl of the final suspension to the Sensititre YeastOne plate within 15 minutes of completing step 2 by either:

a. Sensititre AIM System.

Replace the tube cap with a Sensititre Single-use Dosehead and inoculate the plate according to the Sensititre AIM user manual.

Remove the test tube dosehead combination from the Sensititre AIM System within 30 seconds of dosing a plate and store inverted in a rack or discard.

b. **Manual pipette.** Pour the broth into a sterile seed trough and inoculate the plate using an appropriate pipette.

8. Check the colony count by removing 10µl of the inoculum from the positive control well and plating the inoculum onto Sabauroud Dextrose Agar (SDA). A correct inoculum will produce 50-500 colonies.
9. Cover all wells with an adhesive seal. Avoid creases to prevent well skips.

Incubation

Incubate plates at 35°C in a non-CO₂ incubator.

Incubation temperatures over 35°C may affect performance.

- *Candida* species should be incubated for 24 to 25 hours (refer to step 1 under READING TEST RESULTS section)
- *Cryptococcus* species should be incubated for 72 hours.
- *Aspergillus* species should be incubated for 48 to 72 hours.

Reading Test Results

Plates may be read visually under normal laboratory lighting using a manual mirror viewer or by using the Sensititre Vizion System. Refer to the Sensititre Vizion System user manual for additional information. Yeast growth in the antifungal solutions will be evident as a change in the colorimetric growth indicator from blue (negative) to red (positive). Some yeast species may not change the indicator completely to red but display more of a purpling. Some organisms may show a slight purpling in posaconazole, voriconazole, fluconazole, itraconazole and ketoconazole.

1. Examine the positive growth well after incubation
 - a. *Candida* species should be incubated for 24 to 25 hours (refer to step 1 under READING TEST RESULTS section)
 - b. *Cryptococcus* species should be incubated for 72 hours
 - c. *Aspergillus* species should be incubated for 48 to 72 hours
2. If the growth well is red, the endpoints for the antifungals can be determined. For *Candida* species, if the well is blue or faint purple, re-incubate for an additional 24 hours and re-examine.

DO NOT READ TURBIDITY IN THE SENSITITRE YEASTONE PLATE. Read Only Color Change.

3. The MIC is the lowest concentration of an antifungal agent that substantially inhibits growth of the organism as detected by a color change.

Interpretation of Results

TABLE 1. Illustration and the interpretation of test results that may occur

| Well Concentration $\mu\text{g/ml}$ | | | | | | | R = RED: Positive growth indication B = BLUE: Negative growth indication |
|-------------------------------------|---|---|---|----|----|---|--|
| 1 | 2 | 4 | 8 | 16 | 32 | | |
| A. | R | R | R | B | B | B | Typical growth pattern; MIC endpoint is 8 $\mu\text{g/ml}$. |
| B. | R | R | R | R | R | R | Growth in all wells; MIC endpoint is >32 $\mu\text{g/ml}$. |
| C. | B | B | B | B | B | B | No growth in any well; MIC endpoint is ≤ 1 $\mu\text{g/ml}$. |
| D. | R | R | R | B | R | R | "Skipped Well". MIC endpoint is >32 $\mu\text{g/ml}$. Disregard "skip" when wells on either side have growth. If more than one "skip" should occur in a column, the test results are invalidated ¹ |
| E. | R | R | B | B | R | R | Double "Skipped Well". The test should be repeated ¹ |

¹With careful techniques, these occurrences are not common.

Reading Notes

Amphotericin B.

For amphotericin B at 24 hours, the endpoints are typically easily defined and the MIC is read as the lowest drug concentration that prevents any discernible color change. Trailing endpoints with Amphotericin B are not usually encountered.

The first well showing a distinct color change as compared to the positive growth well is the MIC.

Flucytosine and Azole Antifungals

Candida albicans, *C. glabrata* and *C. tropicalis* with flucytosine and azoles, such as fluconazole, itraconazole, ketoconazole, voriconazole and posaconazole may give endpoints that are

typically less sharp because of trailing growth, and may be a significant source of variability. Trailing occurs when a slight colour change persists and it is often identical for all drug concentrations above the MIC. The MIC should be read as the first well showing a less intense colour change compared to the positive growth control well. Reference strains of defined susceptibility may also help to train personnel. Isolates of *Candida krusei* are assumed to be intrinsically resistant to fluconazole and their MICs should not be interpreted, (1) a comment should accompany the test result reported.

Trailing endpoint: This occurs when a slight color change persists and is often identical in several concentrations. The MIC should be read as the first well showing a less intense color change compared to the positive growth control well.

Echinocandins

MIC end points should be determined after 24 hours of incubation at 35°C. The MIC should be read as the first well showing a less intense color change as compared to the positive control well.


Itraconazole

Itraconazole can occasionally come out of solution at concentrations of ≥ 4 $\mu\text{g/ml}$. This can result in the affected well exhibiting growth and turning red.


From time to time we are encountering paradoxical growth in the higher concentrations of Itraconazole on the Sensititre YeastOne Panels which results in pinking of these wells.

The paradoxical effect also known as the Eagle phenomenon refers to an observation where an increase in the antimicrobial concentration beyond a certain point paradoxically results in an increase to the number of bacteria that survive. An explanation could be that as the concentration is too high, the agent might be self-antagonising the receptor with which it binds (penicillin binding proteins, for example, in the case of a penicillin).

Resolution

The growth in the high concentration should be ignored unless you have growth in all of the other concentrations of Itraconazole. In the example below the 0.5µg/ml well highlighted with the black square is where the MIC result should be recorded.

If you have any other questions or concerns please contact the technical support department on Phone: +44 (0) 1256 694287 | Fax: +44 (0) 1256 463388 or see page 12 for list of emails.

Contamination/ Skips

Alternatively, a pink (growth) well between blue (no growth) wells could be indicative of contamination. Sub-culture well contents to ascertain the cause.

A blue well in a series of red growth wells indicates a “skip” and should be ignored. The MIC should be read above any skip wells. If there is more than one skipped well, the antifungal should not be reported.

Quality Control

Frequency of quality control testing should be conducted according to local guidelines¹

Inoculum should be cultured onto a suitable medium to check for purity. Test results are invalid if a mixed culture is detected.

All Sensititre Plates include positive control wells. Tests are invalid unless there is distinct growth in all positive control wells

The following cultures from the American Type Culture Collection (ATCC™) are recommended for user quality control:

Issatchenkia orientalis (*Candida krusei*) ATCC 6258

Candida parapsilosis ATCC 22019

Results should **not** be reported if QC results are not in range.

The inoculation, reading and interpretation of Sensititre YeastOne susceptibility plates when testing for user quality control should be performed as described on the preceding section.

Table 2. Recommended 24 and 48 hour MIC limits for two quality control strains as per Broth Microdilution CLSI M27 (Ref.2). Ranges that are different or additional to published quality control ranges are underlined.

| Antifungal Agent | <i>Issatchenkia orientalis</i> ATCC6258 | | <i>Candida parapsilosis</i> ATCC 22019 | |
|------------------|--|------------|---|------------|
| | 24 hour | 48 hour | 24 hour | 48 hour |
| 5 – Flucytosine | 4-16 | 8-32 | <u>0.12-0.5</u> | 0.12-0.5 |
| Amphotericin B | 0.5-2 | 1-4 | 0.25-2 | 0.5-4 |
| Anidulafungin | 0.03-0.12 | - | 0.25-2 | - |
| Caspofungin | 0.12-1 | 0.25-1 | 0.25-1 | 0.5-4 |
| Fluconazole | 8-64 | 16-128 | 0.5-4 | <u>2-8</u> |
| Itraconazole | 0.12-1 | 0.25-1 | 0.06-0.5 | 0.12-0.5 |
| Micafungin | 0.06-0.5* | 0.12 – 0.5 | 0.5-2 | 0.5-4 |
| Posaconazole | 0.06-0.5 | 0.12-1 | 0.03-0.25 | 0.06-0.25 |
| Voriconazole | 0.06-0.5 | 0.12-1 | 0.015-0.12 | 0.03-0.25 |

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Expected QC values are provided in CLSI M27; for *Aspergillus*, refer to guidance document M38.

Isolates of *Issatchenkia orientalis* (*Candida krusei*) are assumed to be intrinsically resistant to fluconazole and their MICs should not be interpreted using this scale

NOTE 1: If minimal inhibitory concentrations (MICs) for *Candida* spp. are measured using a scale that yields results falling between categories, the next higher category is implied. Thus an isolate with a fluconazole MIC of 12.5 ug/mL would be placed in the S-DD category.

Please refer to CLSI¹ for more information concerning interpretation of results.

Limitations

1. Sensititre YeastOne Plates are for use with non-fastidious yeast including *Candida* species, *Cryptococcus* species, *Aspergillus* species and miscellaneous rapid growing yeast species. They are not intended for fastidious or slow growing yeast such as *Histoplasma* or *Blastomyces*, and filamentous fungi.
2. Comparison between the Sensititre YeastOne System at 24 hours and the CLSI reference method at 48 hour was evaluated. However due to the difficulty in correlating end points of trailing organisms (*C. albicans*) at 48 hours incubation, high error rates are observed.
3. Testing of fungi and antifungal agents is inherently less precise than testing bacteria.
4. Some investigators believe the 24-hour reading is more appropriate than the 48-hour reading because of the problem with trailing with certain isolates. The CLSI official

standard indicates that readings should be accomplished at 48 hours. Until sufficient data is collected and analyzed, the question of most clinically relevant time of reading remains unanswered. Reporting of results should indicate clearly the times of reading.

5. For additional guidance, refer to CLSI Antifungal Susceptibilities for yeast Standard M27 and for *Aspergillus* refer to M38.
6. Colour change is the indicator of the end point, not turbidity. (This fact alleviates some major concerns with the interpretation of certain *Candida* species because of 'trailing'. Trailing is more commonly seen with isolates other than those of blood and other sterile body fluids.)
7. Do not read at 24 hours if the control well has not completely turned positive.
8. The performance of Voriconazole with *Cryptococcus* species and rapid growing yeast species has not been determined. Voriconazole MIC's should therefore only be reported for *Candida* species.
9. Use only with Sensititre System approved yeast susceptibility inoculum broth. The use of other broths could result in error.
10. As with any in-vitro susceptibility testing method, the results of testing should be correlated with the patient's clinical response to prescribed therapy.
11. Performance has only been established for Amphotericin B, Itraconazole Posaconazole and Voriconazole with *Aspergillus* species.
12. Reference 5 showed high level (>99% agreement) with CLSI method for amphotericin B and *Aspergillus* species. Lower agreement was seen with Itraconazole, *A. fumigatus*, *A. flavus* and *A. terreus* had >90% agreement whilst *A. nidulans* had 85% and *A. ustus* 33% on panels incubated for 48 hours with a 10^3 cfu/ml inoculum. Agreement over 90% was observed with an inoculum of 10^4 cfu/ml and 72 hours incubation with itraconazole and amphotericin B (M38).
13. Correlation of the MIC for Caspofungin to the treatment outcome following caspofungin use has not been fully established. (9)
14. Performance of the Sensititre YeastOne system with Anidulafungin, Caspofungin, and Micafungin with *Cryptococcus*, *Aspergillus* species and rapid growing yeast species **other than *Candida* species** has not been established. Anidulafungin, Caspofungin and Micafungin MIC's should therefore only be reported for *Candida* species.
15. Performance of the Sensititre YeastOne Plate with Posaconazole with *Cryptococcus* has not been established.
16. Only instruments supported by the Sensititre system i.e. a simple mirror viewer, , Sensititre Vizion, , must be used to report results with CE IVD and FDA cleared Sensititre products, any other system used will not be supported.

Performance

Panels read manually are designed to give comparable performance to CLSI reference microbroth procedure. Comparable performance is defined as > 90% essential and categorical agreement to within a doubling dilution of the reference MIC (1).

Bibliography

- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087 (M27).
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holliday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock, D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadiis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS, Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., *etal.* (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology*. **35**: 139-143.
- (7) Pemán, J., *etal.* (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares, M., J., G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., *etal.* (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623
- (10) Pfaller, M., A., *etal.* (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117

- (11) Pfaller, M, A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posiconazole and ravuconazole. *Journal of Clinical Microbiology* **42**: 4577-4580
- (12) Espinel-Ingroff, A., *etal.* (2004) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS, M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721
- (13) Linares, M, J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595
- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459
- (16) Canton, E., *etal* (2005) .Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS, M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607
- (17) Holliday, N.M., *etal* (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract* **C-191**
- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS, M38-A , Sensititre YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 -1129
- (19) Espinel-Ingroff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622
- (20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne® plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001
- (21) Alexander, B. D., T.C. Byrne, K.L. Smith, K.E. Hansen, K.J. Anstrom, J.R. Perfect and L. B. Reller. (2007) Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706
-

(22). Patel, R., *et al.* (2007). Clinical evaluation of the Sensititre YeastOne plate for testing susceptibility of filamentous fungi to posiconazole. *Journal of Clinical Microbiology*. **45**: 2000-2001

DISCLAIMER

The information provided in this technical insert is current at the time of printing and may change without notice.

The latest information can be downloaded from www.trekds.com\techinfo or by
Contacting Thermo Fisher Scientific Microbiology Technical Services.

Manufactured by TREK Diagnostic Systems (Part of Thermo Fisher Scientific Inc)
Units 17-19 Birches Industrial Estate, East Grinstead,
West Sussex. RH19 1XZ, UK.
Tel: +44- 1342-318777

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalog marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. ATCC is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Technical Support email addresses

| | |
|-----------------------------|---|
| Austria | microbiology.techsupport.de@thermofisher.com |
| Belgium | microbiology.techsupport.benelux@thermofisher.com |
| Czech Republic | mikrobiologie.tech.podpora.cz@thermofisher.com |
| Denmark | microbiology.techsupport.nordic@thermofisher.com |
| France | microbiologie.techsupport.fr@thermofisher.com |
| Finland | mikrobiologia.tekninentuki.fi@thermofisher.com |
| Germany | microbiology.techsupport.de@thermofisher.com |
| Italy | microbiologia.supportotecnico.it@thermofisher.com |
| Netherlands | microbiology.techsupport.benelux@thermofisher.com |
| Norway | microbiology.techsupport.nordic@thermofisher.com |
| Spain | microbiologia.soporte.es@thermofisher.com |
| Sweden | microbiology.techsupport.nordic@thermofisher.com |
| UK Technical Support | microbiology.techsupport.uk@thermofisher.com |

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- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087 (M27).
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holliday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock, D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS, Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., *et al.* (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology*. **35**: 139-143.
- (7) Pemán, J., *et al.* (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares, M., J., G. Charriel, F. Solís, F. Rodríguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253

- (9) Kartsonis, N., et al. (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623
- (10) Pfaller, M, A., et al. (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117
- (11) Pfaller, M, A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posiconazole and ravuconazole. *Journal of Clinical Microbiology* **42**: 4577-4580
- (12) Espinel-Ingroff, A., et al. (2004) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS, M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721
- (13) Linares, M, J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595
- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459
- (16) Canton, E., et al (2005) .Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS, M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607
- (17) Holliday, N.M., et al (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract* **C-191**
- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS, M38-A , Sensititre YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 -1129
- (19) Espinel-Ingroff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622

(20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne® plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001

(21) Alexander, B. D., T.C. Byrne, K.L. Smith, K.E. Hansen, K.J. Anstrom, J.R. Perfect and L. B. Reller. (2007) Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706

(22). Patel, R., *etal.* (2007). Clinical evaluation of the Sensititre YeastOne plate for testing susceptibility of filamentous fungi to posiconazole. *Journal of Clinical Microbiology*. **45**: 2000-2001

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Bibliografia

- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087 (M27).
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holliday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock, D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadiis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS, Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., *etal.* (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology*. **35**: 139-143.
- (7) Pemán, J., *etal.* (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares, M., J., G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., *etal.* (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623
- (10) Pfaller, M., A., *etal.* (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117
-

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- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087 (M27).
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holliday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock, D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadiis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS, Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., et al. (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology*. **35**: 139-143.
- (7) Pemán, J., et al. (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares, M., J., G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253.
- (9) Kartsonis, N., et al. (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623.
- (10) Pfaller, M., A., et al. (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50**: 113-117.

- (11) Pfaller, M, A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne[®] antifungal susceptibility testing of the new triazoles voriconazole, posiconazole and ravuconazole. *Journal of Clinical Microbiology* **42**: 4577-4580.
- (12) Espinel-Ingroff, A., *et al.* (2004) Multicenter comparison of the Sensititre YeastOne[®] colorimetric antifungal panel with the NCCLS, M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718-721.
- (13) Linares, M, J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902.
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne[®] colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595.
- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459.
- (16) Canton, E., *et. al* (2005) Sensititre YeastOne[®] Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS, M27-A2 multicenter study. *Antimicrobial Agents and Chemotherapy*. **49**: 1604-1607.
- (17) Holliday, N.M., *et al.* (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne[®] susceptibility plate. *American Microbiology Abstract* **C-191**.
- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS, M38-A , Sensititre YeastOne[®] and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 -1129.
- (19) Espinel-Ingroff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622.
- (20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne[®] plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001.
- (21) Alexander, B. D., T.C. Byrne, K.L. Smith, K.E. Hansen, K.J. Anstrom, J.R. Perfect and L. B. Reller. (2007) Comparative evaluation of Etest and Sensititre YeastOne[®] panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706.
-

(22). Patel, R., *et al.* (2007). Clinical evaluation of the Sensititre YeastOne plate for testing susceptibility of filamentous fungi to posiconazole. *Journal of Clinical Microbiology*. **45**: 2000-2001.

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Das ATCC Licensed Derivative®-Emblem, die ATCC Licensed Derivative®-Wortmarke und die ATCC-

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- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087 (M27).
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holliday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock, D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS, Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., *etal.* (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology*. **35**: 139-143.
- (7) Pemán, J., *etal.* (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares, M., J., G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., *etal.* (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623

- (10) Pfaller, M, A., *etal.* (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117
- (11) Pfaller, M, A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posiconazole and ravuconazole. *Journal of Clinical Microbiology* **42**: 4577-4580
- (12) Espinel-Ingroff, A., *etal.* (2004) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS, M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721
- (13) Linares, M, J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595
- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459
- (16) Canton, E., *etal* (2005) .Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS, M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607
- (17) Holliday, N.M., *etal* (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract* **C-191**
- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS, M38-A , Sensititre YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 -1129
- (19) Espinel-Ingroff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622
- (20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne® plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001

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Bibliographie

- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087 (M27).
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holliday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.

- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock, D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*. *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadiis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS, Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., *etal.* (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology*. **35**: 139-143.
- (7) Pemán, J., *etal.* (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares, M., J., G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., *etal.* (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623
- (10) Pfaller, M., A., *etal.* (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117
- (11) Pfaller, M., A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posiconazole and ravuconazole. *Journal of Clinical Microbiology* **42**: 4577-4580
- (12) Espinel-Ingroff, A., *etal.* (2004) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS, M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721
- (13) Linares, M., J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595

- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459
- (16) Canton, E., *etal* (2005) .Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS, M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607
- (17) Holliday, N.M., *etal* (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract* **C-191**
- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS, M38-A , Sensititre YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 -1129
- (19) Espinel-Ingroff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622
- (20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne® plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001
- (21) Alexander, B. D., T.C. Byrne, K.L. Smith, K.E. Hansen, K.J. Anstrom, J.R. Perfect and L. B. Reller. (2007) Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706
- (22). Patel, R., *etal*. (2007). Clinical evaluation of the Sensititre YeastOne plate for testing susceptibility of filamentous fungi to posiconazole. *Journal of Clinical Microbiology*. **45**: 2000-2001

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- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087 (M27).
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holliday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock, D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS, Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., *etal.* (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology*. **35**: 139-143.
- (7) Pemán, J., *etal.* (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares, M., J., G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., *etal.* (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623
- (10) Pfaller, M., A., *etal.* (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117

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Producent: TREK Diagnostic Systems (część firmy Thermo Fisher Scientific Inc)