diagnostics BX, BV, BXV - faktor

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Dátum poslednej revízie: 1.2.2018

PURPOSE OF USE AND DESCRIPTION

Diagnostic discs BX-factor, BV-factor a BXV-factor are designed to detection and isolation of *Haemophilus sp.* bacteria.

PRINCIPLE

Paper discs are impregnated with bacitracin B (5 IU) and growth factors in form of hemine (BX-factor), nicotinamide adenine dinucleotide - NAD (BV-factor) or with both (BXV-factor). Growth factors are releasing from disc to agar. Different species of *Haemophilus* create zone around particular disc according to nutrition requirements.

DIAGNOSTIC KIT CONTAINS

- tube with diagnostic discs
- BV-factor (Ref. 6037-500 - min. 500pc, Ref. 6037-50 - min. 50pc)
- BX-factor (Ref. 6038-500 - min. 500pc, Ref. 6038-50 - min. 50pc)
- BXV-factor (Ref. 6039-500 - min. 500pc, Ref. 6039-50 - min. 50pc)

WARNINGS AND SECURITY PRECAUTIONS

- For in vitro diagnostics use.**
- For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products has to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Do not use product out of expiration.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of overexpired or devaluated product liquidate according to regulations of your country. Packaging materials are recyclable.

STORAGE

Store diagnostic discs in dark, discs BX-factor at temperature from +2 to +8 °C, discs BV-factor and BXV-factor at temperature from -18 to -10 °C. Opened packages of these discs can be stored in refrigerator.

REQUIRED MATERIALS

- Petri dish with agar
- tampons, loops, burner, tubes and other basic equipment of microbiological laboratory

RECOMMENDED PROCEDURE

Spread 24 hour culture of tested strain on agar without blood (for example Mueller Hinton agar). Place diagnostic disc with sterile loop on agar, so that whole area of disc is in contact with medium. Distance between discs must be at least 2,5 cm. Incubate 18-24 h in 5-7% CO₂ atmosphere at temperature 35-37 °C.

EVALUATION AND INTERPRETATION

Transparent colonies grown around the disk BXV indicate the presence of bacteria of the *Haemophilus* genus. To identify the type use BV and BX discs. Factor requirements

Bacterial strain	BX	BV	BXV
<i>Haemophilus aegyptius</i>	-	-	+
<i>Haemophilus aphrophilus</i>	-	-	+
<i>Haemophilus ducreyi</i>	+	-	+
<i>Haemophilus haemolyticus</i>	-	-	+
<i>Haemophilus influenzae</i>	-	-	+
<i>Haemophilus parahaemolyticus</i>	-	+	+
<i>Haemophilus parainfluenzae</i>	-	+	+

+ = zone creation, - = without zone

Functionality of diagnostic discs can be verified with this control strain:

Bacterial strain	CCM / ATCC	BX + BV	BXV
<i>Haemophilus influenzae</i>	4457 / 49777	+	+
<i>Streptococcus pneumoniae</i>	4501 / 49619	> 17 mm	> 17 mm

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

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