NADAL® FOB Test (test cassette)

REF 272001

DE	Gebrauchsanweisung	2	CZ	Návod k použití	23
EN	Instruction for use	5	FI	Käyttöohje	26
FR	Instructions d'utilisation	8	SE	Användarinstruktioner	29
ES	Instrucciones de uso	11	NO	Bruksanvisning	32
IT	Istruzioni per l'uso	14		Symbols	35
PL	Sposób użycia	17		Our Teams	36
PT	Instruções de Utilização	20			

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14. Referenzen

- Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
- Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
- Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
- Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
- Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
- Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5): 107-110
- St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
- Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

Rev. 3, 2015-10-15 OM/UJ

1. Intended Use

The NADAL® FOB Test is a rapid, visual immunoassay for the qualitative, presumptive detection of human haemoglobin in faecal specimens. The test is intended for use as an aid in the diagnosis of lower gastrointestinal (g.i.) pathologies. The test is designed for professional use only.

2. Introduction and Clinical Significance

Colorectal cancer is one of the most commonly diagnosed types of cancer and a leading cause of cancer-related deaths. Screening for occult blood in faeces is likely to improve the odds of the detection of colorectal cancer at an early stage, thus reducing mortality.

Previously, commercially available FOB tests utilised the Guaiac Method, which requires a special diet prior to testing in order to avoid false positive and false negative results. The NADAL® FOB Test is designed to detect human haemoglobin in faecal samples. The test is based on an immunochemical method that improves specificity of detection of lower gastrointestinal disorders including colorectal cancers and adenomas without any dietary restrictions.

3. Test Principle

The NADAL® FOB Test is used to detect human haemoglobin through the visual interpretation of colour development on the internal strip. Antibodies to human haemoglobin are immobilised in the test line region of the membrane. During testing, the specimen reacts with antibodies to human haemoglobin conjugated with coloured particles and pre-coated on the sample pad of the test. The mixture then migrates along the membrane by capillary action and interacts with components on the membrane. If there is sufficient human haemoglobin in the specimen, a coloured line will form in the test line region of the membrane. The presence of this coloured line indicates a positive test result, while its absence indicates a negative test result. The appearance of a coloured line in the control line region serves as a procedural control indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

Each test kit is suitable for 20 test evaluations.

- 20 individually wrapped test cassettes
- 20 specimen collection tubes with specimens diluent buffer
- 1 package insert

5. Additional Materials Required

- Timer
- Stool collector (Ref. 272004)
- Patient information block with short instructions regarding stool specimen collection (Ref. 272040BL)
- Specimen collection container (optional, if an undiluted stool sample is required and the dilution step is carried out only immediately before testing)

6. Storage & Stability

The test should be stored at 2-30°C until the expiry date printed on the sealed foil pouch. The test should remain in the sealed foil pouch until use. Do not freeze the test. Care should be taken to protect the components of the test kit from contamination. Do not use the test if there is evidence of

microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- For single use only.
- Do not use the test beyond the expiration date indicated on the foil pouch.
- Do not use the test if the foil pouch is damaged.
- The test contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is, therefore, recommended that tests be treated as potentially infectious and handled by observing usual safety precautions (e.g., do not ingest or inhale).
- Avoid cross-contamination of samples by using a new specimen collection tube for each sample.
- Read the entire package insert carefully prior to testing.
- Do not eat, drink or smoke in the area where specimens and test kits are being handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for the proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- The specimens diluent buffer contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of the specimens diluent buffer or extracted samples, always flush with copious quantities of water to prevent azide buildup.
- Bring all reagents to room temperature (15-30°C) before use.
- Do not add any solution to the reaction area (result area).
- Patients should not collect samples during their menstrual period or 3 days before and after it, or if they have bleeding haemorrhoids, blood in urine, or if they have straining during bowel movements.
- Do not substitute or mix components from different test kits.
- Humidity and high temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation

Store the specimen collection tube at room temperature. It is recommended to collect stool specimen using a stool collector. Avoid dilution of stool specimens with urine or water from the toilet bowl.

1. Hold the specimen collection tube upright and remove the light-blue cap. Remove the applicator stick (spiral rod). Be careful not to spill or splash buffer solution from the collection tube.
2. Collect specimens by inserting the applicator stick in three different sites of the stool sample.

Note: Please make sure to insert the applicator stick three times in a row into the stool sample. Do NOT put the applicator

cator stick back into the tube in between collecting the 3 samples. Be careful not to spill fluid from the tube. The nature of the sample and compliance with these instructions will affect your test result.

- Place the applicator stick, along with the sample, back into the specimen collection tube and close it properly.
- Shake the specimen collection tube vigorously to mix the specimen and specimens diluent buffer.

Note:

The NADAL® FOB Test is intended for use with human faecal specimens only. Patients should not collect samples during their menstrual period or 3 days before and after it, or if they have bleeding hemorrhoids, blood in urine or if they experience strain during their bowel movement.


Alcohol, aspirin and other medications taken in excess may cause gastrointestinal irritation, resulting in occult bleeding. Such substances should be discontinued at least 48 hours prior to testing.

No dietary restrictions are necessary before testing.

Perform the test as soon as possible after sample collection. Do not leave specimens at room temperature for prolonged periods of time. On the whole, the storage time of samples at 2-8°C should not exceed 7 days. The transportation of samples to the clinic can be carried out at room temperature without any problems; however it should not exceed 2 hours. Bring specimens to room temperature prior to testing.

9. Test Procedure

- Bring test cassettes and patient's samples (extracted samples) to room temperature (15°C to 30°C) prior to testing.
- Remove a test cassette from the sealed foil pouch and place it on a clean, level surface. Label the test cassette with a patient or control identification. For best results, the assay should be performed within one hour.
- Shake the specimen collection tube thoroughly to ensure proper mixing of the faecal sample with the specimens diluent buffer.

- Unscrew the white protection cap from the specimen collection tube. Using a piece of tissue paper, break the tip of the collection tube using a twisting motion.

- Hold the specimen collection tube vertically and dispense 3 drops of solution into the sample well (S) of the test cassette.

Avoid trapping air bubbles in the sample well (S) and do not add any solution to the result area.

- Read the result after 5 minutes. Do not interpret results after more than 10 minutes.

10. Result Interpretation

Positive:

Two coloured lines appear on the membrane. One line appears in the control line region (C) and the other one in the test line region (T).

Negative:

Only one coloured line appears in the control region (C). No coloured line appears in the test line region (T).

Invalid:

The control line (C) fails to appear.

Results from any test which has not produced a control line at the specified reading time should be discarded. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Note:

The intensity of colour in the test line region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of colour in the test line region should be considered positive. Note that this is a qualitative test only and it cannot be used to determine the concentration of analytes in the specimen.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

An internal procedural control is included in the test cassette.

A coloured line appearing in the control line region (C) is considered an internal positive procedural control, confirming sufficient specimen volume and correct procedural technique.

However, when the faecal sample is tested, the background may appear slightly yellowish due to the original colour of the sample. This is acceptable as long as it does not interfere with the interpretation of test result. The test is invalid if the background fails to clear and obscures the reading of the result.

12. Limitations

- The NADAL® FOB Test is for professional *in-vitro* diagnostic use only. It should only be used for the qualitative detection of human haemoglobin.
- The presence of blood in stool specimens may be due to causes other than colorectal bleeding, such as hemorrhoids, blood in urine or stomach irritation.
- Negative results do not exclude bleeding, since some polyps and colorectal region cancers can bleed intermittently or not at all. Additionally, blood may not be uniformly distributed in faecal samples and colorectal polyps at an early stage may not bleed.
- Urine and excessive dilution of specimens with toilet water may cause erroneous test results.
- The NADAL® FOB Test may show decreased sensitivity for upper gastrointestinal bleeding, as blood degrades as it passes through the gastrointestinal tract.
- Not all colorectal bleeding is due to precancerous or cancerous polyps. As with all diagnostic tests, a confirmed diagnosis should only be made by the physician after all clinical and laboratory findings have been evaluated.

13. Performance Characteristics

Analytical Sensitivity

Specimens containing human haemoglobin with concentration equal to or higher than 40 ng/ml produce positive results. In some cases specimens containing human haemoglobin with concentrations of less than 40 ng/ml can also produce positive results.

Hook or Prozone effect:

The working range of the NADAL® FOB Test is 2 µg/g to 25 mg/g faeces (= 40 ng/ml up to 500,000 ng/ml). At higher concentrations the test shows a "high dose Hook or Prozone effect". In case of presumable Hook-Effect, please dilute the sample and repeat the measurement.

Analytical Specificity:

The NADAL® FOB Test is specific for human haemoglobin and shows no cross-reaction with haemoglobin from bovine, turkey, pig, rabbit, chicken, horse and goat blood at concentrations up to 1 mg/mL.

Interfering Substances

The NADAL® FOB Test shows no cross-reaction with the listed substances:

Analytes	Concentration	Analytes	Concentration
Ascorbic acid	20 mg/dL	Urea	2000 mg/mL
Oxalic acid	60 mg/dL	Glucose	2000 mg/dL
Bilirubin	100 mg/dL	Caffeine	40 mg/dL
Uric acid	60 mg/dL	Albumin	2000 mg/dL
Acetylsalicylic acid	20 mg/dL		

		NADAL® FOB Test		
Another rapid test		+	-	Total
	+	325	9	334
	-	16	1024	1040
	Total	341	1033	1374

Relative Sensitivity:

97.3% (95.56%-99.04%)*

Relative Specificity:

98.4% (97.64%-99.16%)*

Overall Agreement:

98.2% (97.47%-98.89%)*

*95% Confidence Interval

Note:

A recently published study by the Shinshu-University School of Medicine in Japan examined the cost-value ratio of multiple measurements. It shows that the relative sensitivity increases with the amount of tests and its relative specificity slightly decreases.

Results:

Amount of tests	sensitivity	specificity
1	58%	96%
2	89%	95%
3	100%	94%

14. References

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5): 107-110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

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14. Bibliographie

- Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
- Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
- Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7): 2023-2027
- Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1): 41-48
- Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3): 189-190
- Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5): 107-110
- St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
- Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

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14. Reference

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5): 107-110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(3):396-9.

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14. Referenser

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5): 107-110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

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nal von minden GmbH