NADAL[®] FOB Test (test cassette)

REF 272001

DE	Gebrauchsanweisung	2	PL	Sposób użycia	17
EN	Instruction for use	5		Symbols	20
FR	Instructions d'utilisation	8		Our Teams	20
ES	Instrucciones de uso	11			
IT	Istruzioni per l'uso	14			

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Friedenstrasse 32
93053 Regensburg
Germany

Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

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15. Referenzen

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5): 107-110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

Rev. 1.1, 2013-10-21 OM

1. Intended Use

The NADAL® Faecal Occult Blood (FOB) Test Device is a rapid, visual immunochromatographic test for the qualitative detection of human haemoglobin in faecal samples. The test is for professional and *in vitro* use only. This test was created to detect indications of gastrointestinal (g.i.) disorders like bowel cancer and major haemorrhoids.

2. Introduction and/or Clinical Significances

Intestinal cancer is one of the most commonly diagnosed types of cancer and one of the most common causes of cancer related deaths. By testing for occult blood in faeces, indications of intestinal cancer can be detected at an early stage.

Previously, tests for occult blood in faeces were based on the Guaiac Method, which required a special diet prior to testing in order to avoid false positive and false negative results. The high-specificity NADAL® FOB test was designed to detect human haemoglobin in faecal samples. The test is based on an immunochemical method that improves the detection of diseases in the colonic area as opposed to the classical Guaiac tests due to the higher sensitivity and specificity of the NADAL® test.

3. Principle of the Test

The immunological NADAL® FOB rapid test is based upon an immunochromatographic method with a specific antibody reaction, detecting human haemoglobin.

The stool sample containing haemoglobin reacts with the specific monoclonal antibodies, which are bound to the gold particles. This complex spreads over the membrane and reaches the test line (T) which is pre-coated with anti-haemoglobin antibodies.


For a positive result, the haemoglobin molecules, loaded with the gold marked antibodies from the stool sample, are bound onto the test line and become visible through a pink colouring. For a negative result, due to the absence of the haemoglobin molecules, there will be no colouring in the test region. A coloured band will always appear at the control region (C) to serve as a procedural indicator for the proper performance of the test and the device.

4. Reagents and Materials Supplied

Every test package is suitable for 20 test evaluations.

- 20 individually wrapped test cassettes
- 20 sample collection tubes: each contains 2 ml of 0.1 M Phosphate buffered saline (PBS), with Bovine Serum Albumin (BSA) and 0.02% sodium azide.
- 20 patient cards (with patient instructions manual on the back)
- 1 package insert

The plastic case of the test cassette encloses one test strip, with at the right end the sample well and at the left end the window of the reaction zone. In the reaction zone there is the test (T) region and the control (C) region where lines indicate the presence or absence of the analyte after performing the test procedure.

5. Additional Required Materials

- Stopwatch
- Stool container

6. Storage & Stability

The test kit is to be stored refrigerated (2-8°C) or at room temperature (up to 30°C) in the sealed pouch for the duration of the shelf-life.

7. Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- For single use only.
- Do not use the test kit after expiration date.
- Do not use the sample collection tubes from different lots.
- Do not open the test cassette from foil pouch until ready to perform the test.
- All patient samples should be treated as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- Buffered saline contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of buffered saline or extracted samples, always flush with copious quantities of water to prevent azide build up.
- Patients should closely follow the specimen collection procedures.
- Patients should not collect samples during their menstrual period, if they have bleeding haemorrhoids, blood in the urine, or if they have strain during bowel movement.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Do not use test if pouch is damaged.
- Do not spill solution into the reaction zone.
- Do not touch the reaction zone of the device to avoid contamination.
- Avoid cross-contamination of samples by using a new specimen collection container and specimen collection tube for each sample (optional additional materials: Art. No. 272004).
- Bring all reagents to room temperature (15-30°C) before use.
- Humidity and high temperature can adversely affect results.

8. Specimen Collection and Preparation

The sample collection as well as the pre-treatment should be performed by the patient at home. Give the patient the tubes and the patient card (possibly with sample collector). Label

the patient card with the patient's name, date of birth and practice stamp.

The following points are for the patient's attention:

1. Store the test tube at room temperature.
2. Collect the stool, possibly with the use of a stool collector. If this is not possible, avoid dilution of the stool sample with urine or water from the toilet bowl.

3. Keep the collection tube up straight and open the bright blue cap. Remove the bright blue applicator stick.

4. Insert the applicator stick (spiral rod) in three different places of the stool sample.

Note: Please make sure to insert the spiral rod three times in a row into the stool sample. Do NOT put the sample collector back into the tube in between collecting the 3 samples. Continue to watch that no fluid escapes from the tube. The nature of the sample and the compliance with these instructions will affect your test result.

5. Replace the applicator stick containing the sample in the collection tube and close properly.
6. Indicate on the patient's card the date and time of sampling.
7. If the test is not performed immediately, i. e. the test tube is not brought to the clinic on the same day, the closed test tube should be stored refrigerated (2-8°C) and shielded from light.
8. Bring the test tube and patient card to the clinic within 1 or 2 days.

Storage of the sample


If the test is not performed immediately or the sample is taken at the patient's home, the test tube should be refrigerated (2-8°C). After storage (2-8°C), the sample should be brought to room temperature prior to testing. Transportation of the sample can take place at room temperature without problems.

9. Procedure of the Test

1. Test device and patient's samples (extracted sample) should be brought to room temperature (15°C to 30°C) prior to testing.
2. Remove the test device from its pouch when ready to perform the test. Bring the device to room temperature to avoid condensation on the membrane. Label the device with patient or control identification.
3. Shake the collection tube thoroughly to ensure proper mixing of the faecal sample with the extraction solution.
4. Unscrew the white protection cap from the stool sample tubes. Using a piece of tissue paper, break the tip of the collection tube using a twisting motion.

5. Hold the collection tube vertically and dispense 3 drops of solution (app. 120 µL) into the sample well of the test device.
6. Observe the result in 5 minutes. Strong positive results may be observed sooner. Do not interpret results after 8 minutes.

10. Interpretation of the Results

Positive:

Two pink-red coloured bands appear, one in the control region (C) and one in the test region (T). When testing with strong positive samples, the intensity of the control band may be lighter than expected. Comparison of the line intensities is not recommended.

Negative:

Only one pink-red coloured band appears in the control region (C). No apparent faint pink to red coloured band in the test region (T).

Invalid:

A total absence of pink coloured bands in both regions is an indication of procedural error or that test reagents may have deteriorated. Repeat the test with a new test device and if condition persists, contact the manufacturer for technical assistance.

The test lines may get darker after some time. This does not have any affect on the result.

11. Quality Control

An internal procedural control is included in the test. A coloured band appearing on the Control region (C) of the membrane indicates proper performance of the test and the device. However, when the faecal sample is tested, the background may appear slightly yellowish due to the original colour of the sample. This is acceptable as long as it does not interfere with the interpretation of test result. The test is invalid if the background fails to clear and obscures the reading of the result.

12. Expected Values

A positive result appears if haemoglobin is present within the sample. Gastrointestinal blood in the stool can also be due to other factors including hemorrhoids or the mixing of the stool sample with urine containing blood.

Negative results do not exclude the existence of polyps and bowel cancer, some polyps or types of bowel cancer are intermittent or never bleed. In addition, blood is not uniformly distributed within stool samples. Bowel polyps in early stages do not always necessarily bleed.

13. Limitations

- This test kit is to be used for the qualitative detection of human haemoglobin in faecal samples. A positive result indicates the presence of human haemoglobin in faecal samples.
- The presence of blood in stools may be due to several causes, besides colorectal bleeding, such as haemorrhoids, blood in urine or stomach irritations.
- All colorectal bleedings may not be due to precancerous or cancerous polyps. The data obtained by this test should be used in conjunction with other clinical findings and testing methods, such as barium enema, sigmoidoscopy or colonoscopy performed by the physician.
- Negative results do not exclude bleeding since some polyps and colorectal region cancers can bleed intermittently or not at all. Additionally, blood may not be uniformly distributed in faecal samples. Colorectal polyps at an early stage may not bleed.
- Urine and excessive dilution of sample with water from toilet bowl may cause erroneous test results.
- This test may be less sensitive for detecting upper g.i. bleeding because blood degrades as it passes through the g.i. track.

14. Performance Characteristics

Analytical Sensitivity

A sample containing human haemoglobin at concentration equal to or higher than 40 ng/ml produces a positive result. In some cases sample containing human haemoglobin at concentrations less than 40 ng/ml can also be tested as positive.

Hook or Prozone effect:

The working range of the NADAL® FOB test is 2 µg/g to 125 µg/g faeces (= 40 ng/ml up to 2500 ng/ml). At higher concentrations the test shows a "high dose Hook or Prozone effect". In the case of a presumable Hook-Effect, please dilute the sample and repeat the measurement.

Test Specificity:

The NADAL® Faecal Occult Blood Test Device is specific for human haemoglobin and does not show any cross-reaction with the haemoglobin from bovine, pig, rabbit, horse and sheep concentrations up to 0.5 mg/ml.

The NADAL® Faecal Occult Blood Test Device also does not show any cross reaction with bilirubin, vitamin C and horse radish peroxidase.

Note:

In a recently published study of the Shinshu-University School of Medicine in Japan the cost-value ratio of multiple measurements was examined. It shows that the relative sensitivity increases with amount of test and its relative specificity decreases slightly.

Results:

Amount of tests	sensitivity	specificity
1	58%	96%
2	89%	95%
3	100%	94%

15. References

1. Dam, J.V., et. al.; Faecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Faecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5): 107-110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

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15. Referencias

1. Dam, J.V., et. al.; Faecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402.
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094.
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7):2023-2027.
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1): 41-48.
5. Ransohoff, D.F. and Lang, C.A.; Improving the Faecal Occult-Blood Test; The New England Journal of Medicine; (1996) 334 (3):189-190.
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5): 107-110.
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):369-9

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6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5): 107-110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

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15. Bibliografia

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190

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nal von minden GmbH