# NADAL<sup>®</sup> H. pylori Ag Test (test cassette)

REF 262002

<b>DE</b>	Gebrauchsanweisung	2	<b>PL</b>	Sposób użycia	17
<b>EN</b>	Instruction for use	5	<b>PT</b>	Instruções de Utilização	20
<b>FR</b>	Instructions d'utilisation	8		Symbols	23
<b>ES</b>	Instrucciones de uso	11		Our Teams	24
<b>IT</b>	Istruzioni per l'uso	14			

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## 14. Referenzen

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

Rev. 1, 2014-08-05 OM/UJ

## 1. Intended Use

The NADAL® H. pylori Ag Test (Faeces) is a rapid visual immunoassay for the qualitative presumptive detection of *Helicobacter pylori* antigens in human faecal specimens. The NADAL® H. pylori Ag Test is intended to be used as an aid in the diagnosis of *H. pylori* infection.

## 2. Introduction and Clinical Significances

*Helicobacter pylori*, also known as *Campylobacter pylori*, is a spiral-shaped, gram-negative bacterium with typical flagella. The bacterium infects the gastric mucosa by secreting the vacuolating cytotoxin A (VacA) and causes several gastroenteric diseases such as non-ulcerous dyspepsia, active gastritis or gastric and duodenal ulcer. *H. pylori* is classified as carcinogen agent type I and can increase the risk of stomach adenocarcinoma.

Many *H. pylori* strains have been isolated. The one expressing the cytotoxin associated gene A (CagA) antigen is found to be highly immunogenic and is of utmost clinical importance because of its association with cytotoxic factors. It is widely reported in the specialist literature that infected patients possessing antibodies against CagA are at five times higher risk of gastric cancer than the reference group infected with a CagA negative bacterial strain. The presence of the CagA antibodies determines the persistence of *H. pylori* infection and ulceration.

CagA antigen associated to others, such as CagII and CagC seems to act as a starting agent of a sudden inflammatory response, which may provoke peptic ulceration, allergic episodes and decrease the efficacy of therapy.

At present several invasive and non-invasive methods are available to detect *H. pylori* infection. Invasive methods which require endoscopy of the gastric mucosa with a histologic, cultural and urease investigation are cost- and time-consuming.

Non-invasive methods include classical ELISA, immunoblotting assays and breath tests. But the latter are extremely complicated and highly unselective.

## 3. Principle of the Test

The NADAL® H. pylori Ag Test (Faeces) is designed to detect *Helicobacter pylori* through visual interpretation of colour development on the internal strip. *H. pylori* antibodies are immobilised in the test line region of the membrane. During the test, the specimen reacts with coloured *H. pylori* antibody colloidal gold conjugates which are precoated on the sample pad of the test. The mixture then moves along the membrane by capillary action and interacts with reagents on the membrane.

If there are enough *H. pylori* antigens in the specimen, a coloured line appears in the test line region (T) of the membrane. This coloured line indicates a positive result, while its absence indicates a negative result. A coloured line should always appear in the control line region (C) and serves as a procedural control indicating that sufficient volume of sample has been added and that membrane wicking has occurred.

## 4. Reagents and Materials Supplied

- 20 individually packed test cassettes
- 20 specimen dilution tubes with buffer

- 1 package insert

## 5. Additional Required Materials

- Timer
- Patient cards
- Centrifuge (if necessary)
- Micropipette (if necessary)

## 6. Storage & Stability

The test cassette can be stored at 2-30°C until the expiry date printed on the sealed pouch. The test cassette must remain in the sealed pouch until use. Do not freeze. Care should be taken to protect the components of the test kit from contamination. Do not use the test if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

## 7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Do not use the test after the expiration date indicated on the package. Do not use the test if the foil pouch is damaged. Do not reuse tests.
- This test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore, recommended that these products be treated as potentially infectious and handled by observing usual safety precautions (e.g., do not ingest or inhale).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Read the entire procedure carefully prior to testing.
- Do not eat, drink or smoke in the area where the specimens and kits are handled. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for the proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Humidity and temperature can adversely affect results.
- Used testing materials should be discarded according to local regulations.

## 8. Specimen Collection and Preparation




- The NADAL® H. pylori Ag Test (Faeces) is intended only for use with human faecal specimens.
- Perform testing immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Specimens can be stored at 2-8°C for up to 72 hours.
- Bring specimens to room temperature prior to testing.
- Pack specimens in compliance with applicable regulations for transportation of etiological agents, in case they need to be shipped.

## 9. Procedure of the Test

**Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.**

**Specimen collection and pre-treatment:**

1. Use the specimen dilution tubes provided in the test kit for specimen collection. Best results will be obtained if the assay is performed within 6 hours after collection.
2. Unscrew and remove the dilution tube applicator. Be careful not to spill or spatter solution from the tube. Collect specimens by inserting the applicator stick into at least 3 different sites of the faeces to collect approximately 50 mg of faeces (equivalent to 1/4 of a pea).
3. Place the applicator back into the tube and screw the cap tightly. Be careful not to break the tip of the dilution tube.
4. Shake the specimen dilution tube vigorously to mix the specimen and the extraction buffer.

**Testing**

1. Remove the NADAL® H. pylori Ag test cassette from the sealed pouch and place it on a clean, level surface. Label the test cassette with patient or control identification. To obtain the best result, the assay should be performed within one hour.
2. Using a piece of tissue paper, break the tip of the dilution tube. Hold the tube vertically and dispense 2 drops of the solution into the specimen well (S) of the test cassette. Avoid trapping air bubbles in the specimen well (S). Do not drop any solution into the result window. As the test begins to work, you will see colour moving along the membrane.
3. Wait for the coloured line(s) to appear. The result should be read at 10 minutes. Do not interpret the result after 20 minutes.

**Note:** If the specimen does not migrate (presence of particles), centrifuge the extracted specimens in the tubes. Collect 80 µL of supernatant, dispense it into the specimen well (S) of a new test cassette and start afresh following the instructions mentioned above.

**10. Interpretation of the Results****Positive:**

Two coloured lines appear on the membrane. One line appears in the control line region (C) and the other line appears in the test line region (T).

**Negative:**

Only one coloured line appears in the control line region (C). No coloured line appears in the test line region (T).

**Invalid:**

The control line fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a

new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

**Note:**

The intensity of colour in the test line region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of colour in the test line region should be considered as positive result. Note that this is a qualitative test only and cannot be used to determine the concentration of analytes in the specimen.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

**11. Quality Control**

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

External controls are not supplied with this test kit. It is recommended that positive and negative controls be tested as a *Good Laboratory Practice* to confirm the test procedure and to verify proper test performance.

**12. Limitations**

- The NADAL® H. pylori Ag Test (Faeces) is for professional *in-vitro* diagnostic use and should be used for the qualitative detection of *Helicobacter pylori* only.
- During certain antibiotic treatments, the concentration of *H. pylori* antigens can decrease below the minimum detection level of the test. Therefore, diagnosis should be made with caution during antibiotic treatment.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test but should only be made by the physician after all clinical and laboratory findings have been evaluated.

**13. Performance Characteristics**

The NADAL® H. pylori Ag Test vs. Biopsy/Histology/RUT

Biopsy/ Histology/ RUT	NADAL® H. pylori Ag Test		
	+	-	Total
+	132	0	132
-	0	154	154
Total	132	154	286

Relative Sensitivity: >99.9% (97.3% - 100.0%)\*

Relative Specificity: >99.9% (97.6% - 100.0%)\*

Overall Agreement: >99.9% (98.7% - 98.8%)\*

\*95% Confidence Interval

**Cross-reactivity:**

Cross-reactivity of the NADAL® H. pylori Ag Test has been studied with the following organisms at 10<sup>9</sup> organisms/ml. No cross-reactivity was observed with these organisms when tested with the NADAL® H. pylori Ag Test (Faeces).

<i>Staphylococcus aureus</i>	<i>Pseudomonas aeruginosa</i>	<i>Enterococcus faecalis</i>
Group <i>C Streptococcus</i>	<i>Klebsiella pneumoniae</i>	<i>Branhamella catarrhalis</i>
<i>Candida albicans</i>	<i>Proteus mirabilis</i>	<i>Acinetobacter spp</i>
<i>Salmonella choleraesuis</i>	<i>Gardnerella vaginalis</i>	<i>Acinetobacter calcoaceticus</i>
<i>E. coli</i>	<i>Chlamydia trachomatis</i>	<i>Neisseria gonorrhea</i>
Group <i>B Streptococcus</i>	<i>Proteus vulgaris</i>	<i>Enterococcus faecium</i>
<i>Hemophilus influenzae</i>	<i>Neisseria meningitidis</i>	<i>Rotavirus</i>

#### 14. References

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.

3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.

4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996; 100:355-415.

5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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#### 14. Bibliografia

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Gancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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nal von minden GmbH