nal von minden

UrineCheck 7

DOA Adulteration Tests

REF 630925

DE	Gebrauchsanweisung	2	FI	Käyttöohje	18
EN	Instructions for use	4	SE	Användarinstruktioner	20
FR	Instructions d'utilisation	6	DK	Brugervejledning	22
ES	Instrucciones de uso	8	NL	Instructies voor gebruik	24
IT	Istruzioni per l'uso	10		Literature	26
PL	Sposób użycia	12		Symbols	27
PT	Instruções de utilização	14		Our Teams	28
CZ	Návod k použití	16			

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Rev.03, 2018-10-16 UITI

1. Intended Use

The nal von minden UrineCheck 7 is intended for use in the pre-screening procedure of urine samples for drugs of abuse, in order to prevent any chance of adulteration.

2. Introduction

The nal von minden UrineCheck 7 is a ready to use fixed plastic strips, offering up to seven different test zones which can be then disposed after use. Only fresh and non centrifugal urine samples should be used with the test. With the nal von minden UrineCheck 7, urine samples can be tested for commercially available adulteration substances such as Creatinine, Nitrite, Glutaraldehyde, Bleaching agents, Pyridinium chlorochromate and other oxidants as well adulteration via acids or bases through a pH Value. The results from these tests can be used to indicate dilution of the urine or the presence of chemicals prior to the actual drug test being carried out. Dilution for example can be detected through the presence of Creatinine or abnormalities the specific weight of a sample.

3. Principle of test

All seven tests are based on a colour reaction of the indication reagent on the test zone with the substances in the urine sample. The results are read by comparison of the test zone with an added colour card.

Creatinine: Test for dilution

The Creatinine reacts with the creatinine indicator under alkaline conditions to form a violet-brown complex. The concentration level is directly proportional to the colour intensity of the test zone.

Nitrite: Test for added Nitrite

The Nitrite reacts with the aromatic amine in the acid medium. The diazonium reacts with the indicator to form a pink/purple colour.

Glutaraldehyde: Test for Aldehyde

The Aldehyde part of Glutaraldehyde reacts with the indicator to form a pink/purple colour.

pH: Detection of acids or bases

The test is based on the known method of two different pH-indicators leading to easily distinguishable colour values within a wide pH range. The colour ranges from orange (low pH) to yellow and green to blue (high pH).

Specific Gravity: Test for dilution

This test is based on the apparent pH change of certain polyelectrolytes in relation to the ion concentration. In presence of the indicator, where there is low ionic concentration there will be a blue or blue-green colour whilst in urine with a higher ionic concentration there will be a green\ yellow colour present.

Bleach: Testing for the presence of bleach in urine

The presence of bleaching agents is shown by a blue-green, brown, or orange colour complex.

Pyridinium: test for the presence of Pyridinium Chlorochromate in urine

The presence of chromate results in a blue-green colour complex.

4. Material supplied

- 25 test strips in a round box with desiccant
- Package insert

5. Additional required material (not supplied)

- Urine cup
- Paper towel
- Timer

6. Precautions and Warnings

- Do not use test beyond the expiry date.
- For single use only.
- Store tests between 15-30 °C.
- Keep out of direct sunlight.
- Keep test strips in closed box and sealed pouch. Remove only as many strips as needed and close the box immediately after use.
- Do not dispose the desiccant out of the box.
- Test should be disposed according to legal regulations
- Avoid cross-contamination by using a new urine cup for every new sample
- Urine samples can be infectious. Therefore, clean handling and secure storage of test strips is recommended.
- Do not touch the reaction fields/pads.

7. Storage and Stability


The nal von minden UrineCheck 7 can be stored at room temperature (15-30 °C). The test must be stored in original pouch/box. Do not dispose the desiccant out of the box. Keep out of direct sunlight. Remove only as many strips as needed and close the box immediately after use.

8. Sample collection and preparation

The urine sample must be collected in a clean, dry container. The urine sample may be collected at any time. After sampling, the urine should be tested as soon as possible. Samples should not be centrifuged or spiked with preservatives. Fill a small part of the urine sample into another vessel in order to avoid possible contamination of the total urine sample. Execute the drug test as specified in the respective package insert of your test system. Freeze the urine sample, if you cannot test it within an hour. Urine samples can be stored for up to 48 hours at 2-8 °C. For longer storage keep the specimens frozen at -20 °C. Frozen specimens must be thawed prior to testing, brought to room temperature and mix thoroughly.

9. User information

- Remove only as many strips as you can use immediately and close the vial tightly.
- Immerse all test fields completely in fresh, well-mixed urine. Remove the strip immediately to avoid tripping the reagents from the test fields.
- When removing the strip wipe off excess urine from the vessel edge. Gently tap the long edge on a paper towel to further remove excess urine between the wells and thus avoid overflow (contamination from adjacent test areas).
- Compare each test pad with the corresponding colour blocks at the specified time. The correct reading time is critical for optimal results.
- **The results are obtained by direct colour comparison of the colour of the test pads with the respective colour block on the colour card on the round box.**

**Note:**

All adulteration parameters should be read after **1-2 minutes**, to distinguish positive from negative urine. **Colour changes that occur later after 2 minutes, have no diagnostic value.**

10. Interpretation of results

Semi-quantitative results are obtained by visually comparing the colour of each test zone with the appropriate reference colour block, which is shown on the container label or the drug test. No equipment/accessories required.

11. Quality control

For best results, the performance of the test strips should be controlled with positive and negative control samples, tested whenever a new test is made out or a new bottle is opened. Each laboratory should establish its own performance standards and review all operating parameters, if these standards are not met.

12. Limitations

Any comparison made with the colour card will depend on the interpretation of the individual user. It is recommended that all laboratory personnel that interpret the test results are tested for colour blindness. Like all laboratory tests, definitive diagnostic or therapeutic decisions should be based not only on a single test method. Some compounds or physical conditions that can affect the test results are listed below. Drugs that discolour the urine can also lead to anomalous results, because they impair the reactions of the reagents on the test strip.

13. Expected values**Creatine**

The daily creatine excretion of the human body is usually constant depending on muscle mass. The DOT policy states that a Creatinine level less than 20 mg/dl in samples is an indication of manipulation. Although variance can occur due differences in age, sex, diet and muscle mass, samples with Creatinine values of less than 20 mg/dl are considered to be adulterated.

Specific gravity

The specific gravity of urine varies from 1.003 to 1.03. Urine of adults with normal diets and normal fluid intake has an average specific gravity of 1.016 - 1.022. An elevated value can be obtained in the presence of low amounts of protein. The DOT guideline states that a specific gravity <1.003 of the urine sample is proof of manipulation. Specific gravity and creatinine should be tested together in order to obtain a better decision about adulteration.

pH

The normal urine pH-value ranges from 4.5 to 8. Values below 4 or above 9 indicate adulteration

Nitrite

Although nitrite is not a natural component of urine, nitrite levels up to 3.6 mg/dl can be found in some urine samples. This value can be caused by urinary tract infections, bacterial contamination or improper storage. The nal von minden UrineCheck 7 consider nitrite levels above 7.5 mg/dl as abnormal.

Glutaraldehyde

Glutaraldehyde is not a natural component of human urine and therefore should not occur in normal urine. Its presence in urine indicates a possible manipulation. However, a false positive result can be detected when there are ketones in the urine. Ketones can occur in the urine when a person is suffering from ketoacidosis, malnourished or has other metabolic abnormalities.

Bleach

The presence of bleaching agents in the urine indicates manipulation, as they are not natural components of urine.

Pyridinium

A blue-green or grey discolouration indicates adulteration with Pyridinium because is not a natural component of urine.

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

Rev.03, 2018-10-16 HuBu

Literature

1. U.S. Department of Transportation, Drug Testing Procedures Handbook.
2. Young, D.S. et al, Clinical Chemistry, 21 (9), 1975.
3. Friedman, R.B. et.al. Clinical Chemistry: 26 (4), 1980.
4. Tietz, N.W. Clinical Guide to Laboratory Tests, Third Edition, W.B. Saunders Company, p566, 1995.
5. Chu, S.Y. and Sparks. Clinical Chemistry, 17, 1984.
6. Tingsrud, K.M. and Linne, J.J., Urinalysis and body Fluids A color Test and Atlas, Mosby-Year Book, Inc., 1995.
7. Henry, J.B. et al.: Clinical Diagnosis and Management by Laboratory Methods, 16th Ed. Philadelphia: Saunders; (1979)

Rev.03, 2018-10-16 ReWi

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consultense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	Temperaturbegrenzung	Temperature limitation	Limites de température	Limitación de temperatura	Limiti di temperatura	Temperatura przechowywania
LOT	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
REF	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für	Sufficient for	Suffisant pour	Válido para	Sufficiente per	Wystarczający na
	Achtung,	Caution	Attention	Precaución	Attenzione	Uwaga
	Vor Sonnenlicht schützen	Keep away from sunlight	Conserver à l'abri de la lumière du soleil	Manténgase fuera de la luz del sol	Conservare al riparo dalla luce solare	Chronić przed promieniowaniem słonecznym

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegrænsning
LOT	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges
	Prazo de validade	Spotřebujte do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato
REF	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Bestil listingsnummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant
	Suficiente para	Dostačuje pro	Lukumäärä	Räcker till	Voldoende voor	Tilstrækkeligt til
	Atenção	Upozornění	Varoitus	Viktigt	Voorzichtig	Vigtigt
	Manter afastado da luz solar	Chraňte před slunečním světlem	Säilytettävä auringonvalolta suojattuna	Utsätt inte för direkt solljus	Niet blootstellen aan zonlicht	Opbevares mørkt

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel –UK: 0808 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel. Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax: 0800 024 9519

Denmark:

Tel: +31 703075 603
Free Tel: +45 80 88 87 53
Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH