# NADAL® Calprotectin Test (test cassette)

REF 1212001

|           |                      |    |           |                          |    |
|-----------|----------------------|----|-----------|--------------------------|----|
| <b>DE</b> | Gebrauchsanweisung   | 2  | <b>PT</b> | Instruções de Utilização | 16 |
| <b>EN</b> | Instructions for use | 5  |           | Symbols                  | 19 |
| <b>ES</b> | Instrucciones de uso | 7  |           | Our Teams                | 20 |
| <b>IT</b> | Istruzioni per l'uso | 10 |           |                          |    |
| <b>PL</b> | Sposób użycia        | 13 |           |                          |    |

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#### 15. Referenzen

1. Vieira, A. et al., "Inflammatory bowel disease activity assessed by fecal calprotectin and lactoferrin: correlation with laboratory parameters, clinical, endoscopic and histological indexes", BMC Research Notes 2009, 2:221.
2. Hanai, H. et al. «"Relationship Between Fecal Calprotectin, Intestinal Inflammation, and Peripheral Blood Neutrophils in Patients with Active Ulcerative Colitis" Digestive Diseases and Sciences, Sept. 2004, Vol 49, No 9, pp 1438-1443.
3. Bonnin Tomas, A, et al. "Calprotectina fecal como marcador diferencia entre patología gastrointestinal orgánica y funcional". Rev. Esp. de Enf. Dig. 2007, Vol 99, No 12, pp. 689-693.

Rev. 0, 2016-01-05 OM/UJa

### 1. Intended Use

The NADAL® Calprotectin Test is a rapid chromatographic immunoassay (non-invasive assay) for the qualitative detection of calprotectin (hCP) in human faecal specimens. The test is intended for use as an aid in the diagnosis of inflammatory gastrointestinal disorders and is designed for professional use only.

### 2. Introduction and Clinical Significance

Calprotectin is a calcium-binding protein that makes up 5% of the total protein and 60% of the cytosolic protein of neutrophils. It has bacteriostatic and fungistatic properties and is found in faeces at levels six times higher than that in plasma. That faecal biomarker is useful for assessing the activity of inflammatory bowel disease (IBD). IBD includes Crohn's disease (CD) and ulcerative colitis (UC) and is associated with elevated neutrophils.

The NADAL® Calprotectin Test is useful in differentiating organic IBD from functional gastrointestinal disease, such as irritable bowel syndrome (IBS). It is a simple, non-invasive assay that is especially useful for children, who may require general anesthesia for colonoscopy. Moreover, this faecal calprotectin detection can be used to predict relapse.

### 3. Test Principle

The NADAL® Calprotectin Test is a qualitative lateral flow immunoassay for the detection of calprotectin in human faecal samples. Monoclonal antibodies against calprotectin are pre-coated in the test line region of the membrane. During testing, the sample reacts with other anti-human calprotectin antibodies coated onto coloured particles which are pre-dried on the internal test strip. The mixture migrates along the membrane by capillary action. In the case of a positive result, specific antibodies present on the membrane react with the mixture conjugate and generate a coloured line.

A green line should always appear in the control line region (C) and serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

### 4. Reagents and Materials Supplied

- 10 NADAL® Calprotectin test cassettes
- 10 specimen collection tubes with buffer
- 1 package insert

### 5. Additional Materials Required

- Specimen collection container
- Pipette (for liquid stool samples)
- Disposable gloves
- Timer

### 6. Storage & Stability

Store the test as packaged in the sealed foil pouch at 2-30°C. The test is stable until the expiration date printed on the foil pouch. The test cassette should remain in the sealed foil pouch until use. Do not freeze tests.

### 7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.

- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

### 8. Specimen Collection and Preparation

#### Specimen Collection and Storage:

A sufficient quantity of feces should be collected (1-2 g or 1-2 mL for liquid sample). Stool samples should be collected without preservatives or transport media in clean and dry containers. Samples can be stored in the refrigerator (2-8°C) for up to 7 days prior to testing. For longer storage (maximum 1 year) the specimen must be kept frozen at -20°C. In this case, the sample should be completely thawed and brought to room temperature before testing. Repeated freezing and thawing cycles are not recommended.

#### Specimen preparation:

Use a separate specimen collection tube with buffer for each sample. Unscrew the cap of the tube and insert the applicator stick into the faecal specimen four times to collect the sample. Close the tube with the buffer and stool sample. The tube with the sample can be stored at 2-8°C for up to 7 days.

Shake the tube to ensure good sample dispersion. For liquid stool samples, collect the faecal specimen with a pipette and add 15 µL to the specimen collection tube with buffer.


### 9. Test Procedure

**Bring the tests and samples with buffer to room temperature (15-30°C) prior to testing. Do not open foil pouches until you are ready to perform the assay. Use a separate test cassette for each sample.**

1. Remove the NADAL® Calprotectin test cassette

from the sealed foil pouch and use it as soon as possible.

- Shake the specimen collection tube to ensure good sample dispersion. Break off the tip of the tube.
- Dispense exactly 4 drops into the specimen well (S) of the test cassette. Start the timer.
- Read the result at exactly 10 minutes after dispensing the sample.

## 10. Result Interpretation

### Positive:

Two lines appear in the result window. A green line appears in the control line region (C) and a red line appears in the test line region (T).

### Negative:

Only one green line appears in the control line region (C).

### Invalid:

If no green control line (C) develops, the assay is invalid, regardless of the presence or absence of the red test line.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure. In this case, review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit and contact your distributor.

**Note:** The colour intensity of the red line in the test line region (T) may vary depending on the concentration of antigens in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined with this qualitative test.

## 11. Quality Control

An internal procedural control is included in the test cassette:

A green line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

*Good laboratory practice (GLP)* recommends the use of control materials to ensure proper test kit performance.

## 12. Limitations

- The NADAL® Calprotectin Test only indicates the presence of calprotectin in the specimen (qualitative detection) and should be used for the detection of calprotectin in faecal specimens only. Neither the quantitative value nor the rate of increase in calprotectin can be determined with this test.

- An excess of sample may cause an inaccurate result (brown lines appear). If so, dilute the sample with buffer and repeat the test.
- Some stool samples can decrease the intensity of the control green line.
- In the case of patients with active neutrophilic inflammatory bowel diseases such as Crohn's disease and ulcerative colitis, the result may be positive for faecal calprotectin. The NADAL® Calprotectin Test may be used for patients with chronic diarrhea.
- Positive results confirm the presence of calprotectin in faecal samples (nevertheless, it can be due to several causes, inflammatory bowel disease, colorectal cancer and some enteropathies). Positive results should be followed up with additional diagnostic procedures by a physician to determine the exact cause of inflammation.
- Stool samples from patients with non-steroidal anti-inflammatory drug (NSAID) treatment could show a positive result.
- Neonatal faecal calprotectin levels have been reported to be higher than those of older children with a median of 167 µg/g.

## 13. Expected Values

Higher levels of calprotectin in the stool are associated with an increased risk of relapse in patients with inflammatory bowel disease (IBD). Some studies established the value equal or higher than 50 µg hCP/g faeces as a cut-off value to detect adult patients with GI inflammatory disorders.

## 14. Performance Characteristics

### Cut-off

The cut-off of the NADAL® Calprotectin Test is 500 ng/mL (50 µg hCP/g faeces) for human calprotectin.

### Sensitivity and specificity

An evaluation was performed to compare the results obtained with the NADAL® Calprotectin Test and those obtained with a commercially available immunoassay (Calprest®, Eurospital).

The following results were obtained using the NADAL® Calprotectin Test:

Sensitivity >94% and specificity >93%

### Cross-reactivity

An evaluation was performed to determine the cross-reactivity of the NADAL® Calprotectin Test. There was no cross-reactivity found with faecal markers occasionally present in faeces:

|                              |                   |
|------------------------------|-------------------|
| Bovine and swine hemoglobin  | Human hemoglobin  |
| Bovine and swine transferrin | Human lactoferrin |
| Bovine lactoferrin           | Human transferrin |

## 15. References

- Vieira, A. et al., "Inflammatory bowel disease activity assessed by fecal calprotectin and lactoferrin: correlation with laboratory parameters, clinical, endoscopic and histological indexes", BMC Research Notes 2009, 2:221.
- Hanai, H. et al. «Relationship Between Fecal Calprotectin, Intestinal Inflammation, and Peripheral Blood Neutrophils in Patients with Active Ulcerative Colitis" Digestive Diseases and Sciences, Sept. 2004, Vol 49, No 9, pp 1438-1443.
- Bonnin Tomas, A. et al. "Calprotectina fecal como marcador diferencia entre patologia gastrointestinal orgánica y funcional". Rev. Esp. de Enf. Dig. 2007, Vol 99, No 12, pp. 689-693.

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## 15. Referencias

1. Vieira, A. et al., "Inflammatory bowel disease activity assessed by fecal calprotectin and lactoferrin: correlation with laboratory parameters, clinical, endoscopic and histological indexes", BMC Research Notes 2009, 2:221.
2. Hanai, H. et al. «"Relationship Between Fecal Calprotectin, Intestinal Inflammation, and Peripheral Blood Neutrophils in Patients with Active Ulcerative Colitis" Digestive Diseases and Sciences, Sept. 2004, Vol 49, No 9, pp 1438-1443.
3. Bonnin Tomas, A. et al. "Calprotectina fecal como marcador diferencia entre patologia gastrointestinal orgánica y funcional". Rev. Esp. de Enf. Dig. 2007, Vol 99, No 12, pp. 689-693.

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#### 15. Bibliografia

1. Vieira, A. et al., "Inflammatory bowel disease activity assessed by fecal calprotectin and lactoferrin: correlation with laboratory parameters, clinical, endoscopic and histological indexes", BMC Research Notes 2009, 2:221.
2. Hanai, H. et al. « "Relationship Between Fecal Calprotectin, Intestinal Inflammation, and Peripheral Blood Neutrophils in Patients with Active Ulcerative Colitis" Digestive Diseases and Sciences, Sept. 2004, Vol 49, No 9, pp 1438-1443.
3. Bonnin Tomas, A. et al. "Calprotectina fecal como marcador diferencia entre patologia gastrointestinal orgánica y funcional". Rev. Esp. de Enf. Dig. 2007, Vol 99, No 12, pp. 689-693.

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


## 15. Referências

1. Vieira, A. et al., "Inflammatory bowel disease activity assessed by fecal calprotectin and lactoferrin: correlation with laboratory parameters, clinical, endoscopic and histological indexes", BMC Research Notes 2009, 2:221.
2. Hanaj, H. et al. «"Relationship Between Fecal Calprotectin, Intestinal Inflammation, and Peripheral Blood Neutrophils in Patients with Active Ulcerative Colitis" Digestive Diseases and Sciences, Sept. 2004, Vol 49, No 9, pp 1438-1443.
3. Bonnin Tomas, A. et al. "Calprotectina fecal como marcador diferencia entre patologia gastrointestinal orgánica y funcional". Rev. Esp. de Enf. Dig. 2007, Vol 99, No 12, pp. 689-693.

Rev. 0, 2016-01-05 MNL

| Symbol   | Deutsch                        | English                               | Français                                     | Español  | Italiano                                    | Polski                             |
|--|--------------------------------|---------------------------------------|--|--|---|------------------------------------|
|  | CE<br>Konformitätszeichen      | CE marking of<br>conformity           | Conformité aux<br>normes<br>européennes      | Conformidad europea                                | Conformità europea                          | Znak zgodności CE                  |
|  | Gebrauchsanweisung<br>beachten | Consult instructions<br>for use       | Consulter la notice<br>d'utilisation         | Consúltense las<br>instrucciones de uso            | Consultare le<br>istruzioni per l'uso       | Przestrzegać instrukcji<br>obsługi |
|  | In-vitro-Diagnostika           | In-vitro diagnostic<br>medical device | Dispositif médical de<br>diagnostic in vitro | Producto sanitario<br>para diagnóstico<br>in vitro | Dispositivo medico-<br>diagnostico in vitro | Tylko do diagnostyki in<br>vitro   |
|  | Temperaturbegrenzung           | Temperature<br>limitation             | Limites de<br>température                    | Límite de temperatura                              | Limiti di temperatura                       | Temperatura<br>przechowywania      |
|  | Chargenbezeichnung             | Batch code                            | Code du lot                                  | Código de lote                                     | Codice lotto                                | Numer serii                        |
|  | Nicht zur<br>Wiederverwendung  | Do not reuse                          | Ne pas réutiliser                            | No reutilizar                                      | Non riutilizzare                            | Tylko do jednorazowego<br>użytku   |
|  | Verwendbar bis                 | Use by                                | Utiliser jusqu'au                            | Fecha de caducidad                                 | Utilizzare entro                            | Data ważności                      |
|  | Bestellnummer                  | Catalogue Number                      | Référence du<br>catalogue                    | Número de catálogo                                 | Riferimento di<br>Catalogo                  | Numer katalogowy                   |
|  | Hersteller                     | Manufacturer                          | Fabricant                                    | Fabricante   | Fabbricante                                 | Producent                          |
|  | Ausreichend für <n><br>Ansätze | Sufficient for<br><n> tests           | Suffisant pour "n"<br>tests                  | Suficiente para <n><br>utilizaciones               | Sufficiente per<br>"n" saggi                | Wystarczający na <n><br>Powtórzeń  |

| Symbol   | Português  | Český   | Suomi  | Svenskt   | Nederlands  | Dansk                                       | Norsk                                  |
|--|--|---|--|---|---|---|--|
|    | Conformidade<br>com as normas<br>européias         | CE certifikát                                       | CE-merkitty  | CE-märkning   | CE-markering                                      | CE-mærkning                                 | CE standardisert                       |
|   | Consultar as<br>instruções de<br>utilização        | Viz návod k použití                                 | Katso<br>käyttöohjetta   | Läs bruksanvisningen  | Raadpleeg de<br>gebruiksaanwijzing                | Se bruksanvisningen                         | Les<br>bruksanvisning<br>nøye          |
|  | Dispositivo médico<br>para diagnóstico in<br>vitro | Diagnostický<br>zdravotnický<br>prostředek in vitro | In vitro -<br>diagnostiikkaan<br>tarkoitettu<br>lääkinnällinen laite | Medicinteknisk<br>produkt avsedd för<br>in vitro-diagnostik | Medisch<br>hulpmiddel voor<br>in-vitrodiagnostiek | Medicinsk udstyr til<br>in vitro-diagnostik | In-vitro diagnostic<br>medisinsk enhet |
|  | Limites de<br>temperatura                          | Teplotní omezení                                    | Lämpötilarajat   | Temperatur-<br>begränsning                                  | Temperatuurlimiet                                 | Temperatur-<br>begrænsning                  | Temperatur<br>begrænsning              |
|  | Código do lote                                     | Kód řáže  | Eräkoodi   | Satsnummer  | Code van de partij                                | Batchcode                                   | Merking                                |
|  | Não reutilizar                                     | Pro jednorázové<br>použití                          | Kertakäyttöinen  | Får inte<br>återanvändas                                    | Niet opnieuw<br>gebruiken                         | Må ikke genbruges                           | Må ikke brukes om<br>igjen             |
|  | Prazo de validade                                  | Spotřebuje do                                       | Käytettävä<br>viimeistään  | Används före  | Houdbaar tot                                      | Udløbsdato                                  | Tidtaking                              |
|  | Número de<br>catálogo                              | Katalogov číslo                                     | Luettelonumero   | Listnummer  | Catalogus nummer                                  | Best il l ingsnummer                        | Katalog nummer                         |
|  | Fabricante   | Výrobce   | Valmistaja   | Tillverkare   | Fabrikant   | Fabrikant                                   | Produsent                              |
|  | Suficiente para<br><n> test                        | Dostačuje pro <n><br>testů                          | Lukumäärä <n><br>test  | Räcker till <n> test  | Voldoende voor<br><n> test                        | Tilstrækkeligt til <n><br>test              | Tilstrækkelig<br>for<n> tester         |

## Our Teams

### Germany:

#### Regensburg

Tel: +49 941 290 10-0  
Fax: +49 941 290 10-50

#### Moers

Tel: +49 2841 99820-0  
Fax: +49 2841 99820-1

#### Austria:

Tel: +49 941 290 10-29  
Free Tel: 0800 291 565  
Fax: +49 290 10-50  
Free Fax: 0800 298 197

#### UK & Ireland:

Tel: +49 941 290 10-18  
Free Tel – UK: 0808 234 1237  
Free Tel – IRE: 1800 555 080  
Fax: +49 290 10-50

#### France:

France Tel: 0800 915 240  
France Fax: 0800 909 493

### Switzerland

Swiss Tel: 0800 564 720  
Swiss Fax: 0800 837 476

### Belgium

Belgium Tel: 0800 718 82  
Belgium Fax: 0800 747 07

### Luxembourg

Lux. Tel: 800 211 16  
Lux. Fax: 800 261 79

### Spain:

Tel: +49 941 290 10-759  
Free Tel: 900 938 315  
Fax: +49 941 290 10-50  
Free Fax: 900 984 992

### Italy:

Tel: +49 941 290 10-34  
Fax: +49 941 290 10-50

### Poland:

Tel: +49 941 290 10-44  
Free Tel: 00 800 491 15 95  
Fax: +49 941 290 10-50  
Free Fax: 00 800 491 15 94

### Portugal:

Tel: +49 941 290 10-735  
Tel. Verde: 800 849 230  
Fax: +49 941 290 10-50  
Fax Verde: 800 849 229

### Netherlands:

Tel: +31 30 75 600  
Free Tel: 0800 0222 890  
Fax: +31 70 30 30 775  
Free Fax: 0800 024 9519

### Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607  
Free Tel: +45 80 88 87 53  
Tax: +31 703030 775

### Laboratory Diagnostics Team:

Tel: +49 941 290 10-40  
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH