NADAL® hCG pregnancy test (test cassettes)

REF 142000/ 152000
142002/ 152002

DE	Gebrauchsanweisung	2	PL	Sposób użycia	24
EN	Instruction for use	7		Symbols	28
FR	Instructions d'utilisation	11		Our Teams	28
ES	Instrucciones de uso	16			
IT	Istruzioni per l'uso	20			

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NADAL® hCG Pregnancy Test (test cassette)

(Ref. 142000/ 152000/ 142002/ 152002)

This Package Insert contains information on the following nal von minden NADAL® hCG Pregnancy Tests:

Product Nr.:	Cut-off-Level	Packaging
142000	10 mIU/ml in Urine (Early Test)	Test Cassette, individually sealed
152000	10 mIU/ml in Urine or Serum (Early Test)	Test Cassette, individually sealed
142002	25 mIU/ml in Urine	Test Cassette, individually sealed
152002	25 mIU/ml in Urine or Serum	Test Cassette, individually sealed

1. Intended Use

The NADAL® hCG Pregnancy Test (142002/152002) or the Early Pregnancy Test (142000/152000) is a one-step test kit for the rapid detection of the pregnancy hormone hCG in urine or in serum. This test provides a visual, qualitative result and has been developed for professional use.

2. Introduction and Clinical Significances

Human chorionic gonadotropin (hCG) is a glycoprotein hormone produced by the developing placenta shortly after fertilization and can already be detected in the early stage of a pregnancy in urine and serum.

In normal pregnancy, hCG can be detected in both serum and urine as early as 7 days after conception. The hCG level continues to rise very rapidly, frequently exceeding 100 mIU/ml by the first missed menstrual period, and culminates in the 10th to 12 week of pregnancy in a range 11500-300000 mIU/ml.

The early appearance of hCG in urine and blood soon after conception and its rapid rise in concentration make it an ideal marker for the early detection of pregnancy.

Elevated hCG levels in urine, comparable to the early stages of pregnancy may be associated with a hydatidiform mole or choriocarcinoma. Therefore, these possibilities should be ruled out before a pregnancy is diagnosed.

The hCG Cassette Test is a rapid chromatographic immunoassay for the detection of hCG at a concentration of 10 mIU/ml or 25 mIU/ml in urine (or serum with 152000 and 152002).

The test uses monoclonal antibodies to detect hCG in the urine (or serum). The specificity of this immune reaction screens out a cross-reactivity interference from structurally related hormones (hFSH, hLH and hTSH) at physiological levels.

3. Principle of the Test

The NADAL® hCG Cassette Test is a chromatographical immunoassay test for the rapid qualitative detection of hCG in urine and serum samples from a concentration of 10 mIU/ml or 25 mIU/ml in urine or urine/serum.

In the test highly specific antibodies are used against hCG. The specificity of this immunological linkage shows no cross-reactivity interference from structurally related hormones (hFSH, hLH and hTSH) at physiological levels.

Antibodies against the α -chain of hCG are bound to the membrane in the test region and goat anti-mouse antibodies are bound to the control region. During the test, the sample is

first mixed with colloidal gold particles, which have monoclonal antibodies against β -chain of hCG on its surface. All particles are then transported through the capillary and over the membrane through the liquid. If a positive result occurs, i.e. showing the presence of hCG, a coloured line will appear in the test region. If there is no coloured line in the test region, the test result is negative. A coloured line should always appear in the control region. It shows the successful execution of the test.

4. Reagents and Materials Supplied

Within the plastic case of the test cassette a testing strip is located. On the right end of the cassette, the sample opening and on the left end, the reactionfield window can be found. After addition of the sample in the round sample opening (contact point for urine), coloured lines in the test and the control region will become visible. To facilitate the allocation, the control region is labeled with C. The test region – where a line will appear only for positive results – is marked with T.

Materials supplied:

- 1 User Manual
- 30 Test cassettes, individually wrapped
- 1 single use pipette per Test cassette

5. Additional Required Materials

- Collection container
- Stop-watch

6. Storage and Stability

The test should be stored up until the expiration date in the original packaging at 2-30°C.

7. Warnings and Precautions

- Only for professional *in-vitro* use.
- Only for single use.
- Sample must not touch the reaction field (results field).
- Avoid cross-contamination of samples by using a new specimen collection container and specimen pipette for each sample.
- Do not touch the reaction field to avoid contamination.
- Only urine samples (or Serum samples with 152000/152002) should be used – no other liquids.
- Bring the Urine- (or Serum-) sample to room temperature (18-30°C).
- Do not use after the expiration date.
- Do not use test if pouch is damaged.
- Use the test immediately after removing it from the foil packaging or the box as it is very sensitive to humidity.
- Be aware of necessary sample volume.
- Be aware of the specified evaluation times.
- Be aware of cross-reactivity.
- Test should be stored and transported at the correct temperature.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological

NADAL® hCG Pregnancy Test (test cassette)

(Ref. 142000/ 152000/ 142002/ 152002)

hazards throughout testing and follow standard procedures for proper disposal of specimens.

For serum samples:

- Using serum samples that contain traces of HAMA (human anti-mouse antibodies) can lead to false-positive results. That is why, in case of a positive result, you should test the serum sample for the presence of HAMA. (e.g. LUCIO® Medical HAMA ELISA, Art.Nr. ELI-6060). Likewise other heterophilic antibodies can affect the test result.
- Different kinds of hCG can be predominant in a serum sample depending on the stage of pregnancy. That is why, with group results, it is recommended to verify the result by urine sample.
- When using serum, results should only be interpreted in conjunction with another test method or clinical finding.

8. Specimen Collection and Preparation**Urine samples:**

- The urine specimen must be collected in a clean and dry plastic or glass container without any preservatives. The first morning urine is preferred since it generally contains the highest concentration of hCG. However, urine collected at any time of day may be used.
- The specimens and test must be at room temperature for the correct execution of the procedure. Urine samples exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain clear supernatant for testing. Only the supernatant may be used for testing.
- If testing is not immediate, the urine specimens may be stored at 2-8°C for up to 48 hours. For long-term storage, specimens should be kept below -20°C.

Serum sample (only for 152000 and 152002):

- Collect blood aseptically by venipuncture, avoiding hemolysis, into plain, heparinised or EDTA tubes.
- Separate serum or plasma from cruro as soon as possible to avoid hemolysis. Use only clear, non-hemolysed specimens. Specimens containing particulate matter should be clarified by centrifugation prior to assay.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- Ideally, testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept below -20°C.

9. Procedure of the Test

1. Please follow the instructions for sampling. Let the sealed test cassettes and the urine-/serum samples or the controls reach room temperature (18-30°C) before use.
2. Remove the test cassette from the packaging and label it with a characteristic for the sample which is to be used.
3. Put the test cassette on a clean and even surface. Pipette around 0.12 mL (ca. 120 µl) of the sample liquid (marked volume or 3 drops) into the round sample opening (S).

Avoid the formation of bubbles. While pipetting, the pipette should be used vertically and use a new tip for every new sample. At this point set the timer.

4. Wait for the coloured lines to appear. The result must be interpreted after the following reaction times:

- 10 mIU/ml Early test (142000/ 152000):
→ 5-10 minutes
- 25 mIU/ml (142002/ 152002):

3 min for urine

3 min

5 min for serum

5 min

5-10 min

5. Do not interpret results after the appropriate read time. It is important that the background is clear before the result is read.

Note: Results which are interpreted later might detect significantly less amounts of hCG. This increases the risk for false-positive results. Low hCG values can appear if natural abortions, certain tumours, a high physiological hCG concentration of not pregnant women and the use of medicaments containing hCG occurs.

10. Interpretation of the Results


To interpret the test, the reaction field is observed. Within this field one or two red lines appear.

POSITIVE (pregnant):

Two coloured lines appear. In addition to the line in control region (C), a second line becomes visible in the test region (T). The hCG concentration corresponds to one appearing during a pregnancy.

**NEGATIVE (not pregnant):**

Only one red line appears in the result window – the control line (C). In the test region (T) no line will appear.

**INVALID:**

If no line appears the test has not been performed correctly or the reagents reached their expiry date. If no red line appears even after repeated dropping of the sample liquid, discard the test; it is invalid.

**Notes for the Interpretation:**

- In the case of a strong T-line, the (upper) control line C can be very pale. The test result will still be valid.
- In case of a negative test result (under circumstances with strongly diluted samples) even though a pregnancy is suspected, the test should be repeated 48 to 72 hours and in addition, a quantitative test should be performed since the

hCG concentration might be not high enough. The first morning urine contains the highest hCG concentration.

- Already shortly after nidation in the womb hCG-concentrations (less than 50 mIU/ml) in urine or serum are very low. However since a significant number of pregnancies (ca. 31%) end naturally after the first 3 months, it is recommended to repeat the test at a later point in time even with a very weakly positive test result. A repeatedly negative result can be caused by a spontaneous or induced abortion. In that case the hCG concentration will decrease. 22% of the clinically not detected pregnancies are ended by an unnoticed abortion and become apparent as a late period. If however during a repetition an intensely coloured line is visible in the test region, an increase of the hCG-concentration can be detected.
- The colour intensity in the test region is proportional to the hCG concentration in the sample. However, it should not be used for the determination of the hCG concentration or the growth rate of the hCG concentration since the test can only be performed qualitatively.
- A sample hCG concentration below the cut-off level of this test might result in a weak line appearing in the test region (T) after an extended period of time. A line in the test region (T) seen after the read time could be indicative of a low hCG level in the sample. If such results are seen, it is recommended that the test be repeated with a new sample in 48-72 hours or that an alternate confirmation method is used.

11. Quality Control

A procedural control is included in the test. The appearance of a coloured line in the control region (C) confirms the correct performance of the test and the reactivity of the reagents.

12. Limitations

Factors taking influence on the test result:

Various factors can – even if the patient is not pregnant – cause a hCG concentration of over 10 mIU/ml. Among these factors are tumour diseases (e.g. chorion carcinomas), in individual cases changing years and hCG containing medication.

Genito-urinary infections can lead to a false-positive hCG test result. An increased number of leucocytes, a basic pH, an increased amount of nitrite, erythrocytes and proteins in urine among others are indicative of genito-urinary infections. In case of suspicion, please verify this by means of urinalysis strips or specific swab tests. Very suspicious urine samples should be verified by laboratory testing.

Therefore, in case of the presence of hCG in urine or serum – determined by the hCG cassette test – no pregnancy should be diagnosed until the possibilities mentioned before are excluded.

The hCG concentration of a strongly diluted urine sample (i.e. low specific weight) can be below the cut-off. If a pregnancy is still suspected, the test should be repeated after 48-72 hours with the first morning urine of the patient. Even a very weak line in the test region should not be diagnosed as “pregnant” immediately but should be reason for a repetition of the test.

Furthermore, the final diagnosis should not be based on the result of one test alone. The diagnosis should be made after evaluation of all relevant findings from the clinical and laboratory examination by a doctor.


Immunologically effective substances, as used within an antibody therapy can lead to an invalid test result.

Expected hCG concentrations:

With healthy men and healthy, not pregnant women, no hCG can be detected with the hCG-cassette test.

The hCG concentration in urine or serum of pregnant women increases rapidly after nidation in the womb.

A concentration of 25 mIU/ml or higher is present 7-10 days after conception (equal to 4-5 days before absences of the menstruation). In that sense the hCG cassette test already can show hCG on the first day after the absence of the period – the hCG concentration can reach values of up to 100 mIU/ml. The hCG concentration reaches maximum values in a range 11500-300000 mIU/ml in the 10th to 12 week of pregnancy (2-4 months after the first missed menstrual period).

Finally, the hCG concentration will decrease to low values during the remaining time of the pregnancy. After giving birth, the hCG concentration decreases rapidly and will reach its normal value only a few weeks after birth.

13. Performance Characteristics

The lower cut-off of the hCG cassette test was set at 10 mIU/ml (for 142000/152000) or 25 mIU/ml (for 142002/152002) (Calibration with the 4th IRP of hCG).

Sensitivity

The NADAL® pregnancy test does not show any prozone effect up to a maximal concentration (600000 mIU/ml).

To evaluate the sensitivity of the hCG cassette test, a study with > 100 test persons was performed. It was proved that after 5 minutes the early test (142002/152002) showed a sensitivity of 10 mIU/ml or the standard test (142002/152002) of 25 mIU/ml. With samples with low hCG concentrations a weak red line may appear in the test region.

To evaluate the sensitivity of the hCG cassette test, the following trials were performed. Urine samples of five not-pregnant persons were set by addition of hCG at concentrations of 0, 10, 25, 40 and 100 mIU/ml. 100 of these samples were tested blindly with the hCG cassette test.

NADAL® hCG Pregnancy Test (test cassette)

(Ref. 142000/ 152000/ 142002/ 152002)

hCG (mIU/ml)	0	10	25	40	100
Samples	20	20	20	20	20
negative	20	5	0	0	0
positive	0	15	20	20	20

Specificity**Cross Reactivity**

The specificity of the hCG cassette test was determined by the use of cross reactivity studies with known amounts of LH (300 mIU/ml), FSH (1000 mIU/ml) and TSH (1 mIU/ml). This is equal to a multiple amount of the physiologically maximum concentration to be expected. A cross reactivity could not be detected.

Interferences

hCG free urine samples and urine samples with 25 mIU/ml hCG were spiked with numerous connections. These connections did not interfere with the test results.

Substance	Concentration
Paracetamol	0.2 mg/ml
Acetylsalicyl acid	0.2 mg/ml
Ascorbin acid	0.2 mg/ml
Atropin	0.2 mg/ml
Coffeine	0.2 mg/ml
Gentisin acid	0.2 mg/ml
Glucose	20 mg/ml
Hemoglobine	10 µg/ml

A full list of tested substances is available on request.

Analytical Specificity

In a clinical multicenter study the results of the NADAL® hCG Combo pregnancy test were compared with another commercially available hCG membranetest for urine and serum testing. For the urine study, 159 samples were examined and both test methods showed 88 negative and 71 positive results. For the serum study, 73 samples were examined and both test methods showed 51 negative, 21 positive and 1 invalid result. The results showed a trueness of >99% of the NADAL® hCG pregnancy test Combo (test cassette for urine/serum) and the second hCG membrane test (urine/serum).

hCG-Reference method (Urine)

Method		other hCG rapid test		Result
hCG Test Cassettes	Result	Positive	Negative	
	Positive	71	0	71
	Negative	0	88	88
Result		71	88	159

Sensitivity: 100% (95-100%)*

Specificity: 100% (96-100%)*

Trueness: 100% (98-100%)*

*95%-confidence interval

hCG-Reference method (Serum)

Method		other hCG rapid test		Result
hCG Test Cassettes	Result	Positive	Negative	
	Positive	21	0	21
	Negative	0	51	51
Result		21	51	72

Sensitivity: 100% (84-100%)*

Specificity: 100% (93-100%)*

Trueness: 100% (95-100%)*

*95%-confidence interval

Standardisation

The hCG strip test was standardised with the 4th international WHO reference standard preparation (IRP 75/589).

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nal von minden GmbH