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
Medtronic

Resterilization Guideline of AF Solutions Catheter Cables

CryoAblation			Phased RF Ablation	EP Systems
Cables	Coaxial Umbilical (203CX) Electrical Umbilical (2035U)	Achieve™ Electrical Cable (990066)	Catheter Interface Cable (990004)	RF Catheter Cable (05116) RF Generator Adaptor (04821) Catheter Connecting Cables (05504SP, 05514SP, 05518SP, 05508SP, 4819)
Single-Use	Supplied Sterile Do not reuse or reprocess a sterile device for purpose of reuse. Reuse or reprocessing may compromise the structural integrity of the device or create a risk of contamination from the device that could result in patient injury, illness, or death.	Supplied Sterile Do not reuse or reprocess a sterile device for purpose of reuse. Reuse or reprocessing may compromise the structural integrity of the device or create a risk of contamination from the device that could result in patient injury, illness, or death.	Supplied Sterile	Supplied sterile OR non-sterile 05116 is supplied sterile 04821 is supplied non-sterile 05504SP, 05514SP, 05518SP, 05508SP, and 4819 are supplied sterile
Resterilization	Do not resterilize a single-use device for purpose of reuse. Resterilization may compromise the structural integrity of the device or create a risk of contamination from the device that could result in patient injury, illness, or death.	Do not resterilize a single-use device for purpose of reuse. Resterilization may compromise the structural integrity of the device or create a risk of contamination from the device that could result in patient injury, illness, or death.	Can be resterilized If repeat sterilization is required, it is recommended that the Catheter Interface Cable be steam sterilized using the following parameters: Method of sterilization: Steam • Steam cycle temperature: 127 °C (260 °F) • Steam cycle time at temperature: 35 ± 5 minutes • Cool-down time: 30 minutes minimum If steam sterilization parameters cannot be followed, it is recommended that a sterile probe cover be used to fully enclose the catheter interface cable before use.	Can be resterilized Do not steam or dry heat sterilize any component. The cables may be sterilized with ethylene oxide (EtO). No other components may be sterilized; otherwise, the equipment may be damaged. The cable must be completely dry before being sterilized. 1. Remove the cable from its original package prior to sterilization. The original package is not a suitable sterilization packaging material. 2. Wrap the cable in packaging permeable to ethylene oxide. 3. Sterilize the cable using procedures validated for effectiveness. Typical conditions for sterilization are: • 500 mg/L ethylene-oxide • 50%-70% relative humidity • 50-55 °C with exposure and aeration time dependent on the particular vessel 4. Use an acceptable method for determining sterilization effectiveness, such as biological indicators. 5. Visually inspect the cable and connectors. 6. Do not use the cable if it appears damaged. Damage includes, but is not limited to, deterioration of the cable insulation (brittleness, cracking, thinning, or bare spots).
# of Times Cable Can Be Reused	NA	NA	Up to 10 times if the recommended sterilization parameters are followed. OR if a sterile probe cover is used in place of steam sterilization, there is no restriction for the reuse of the Catheter Interface Cable.	Not specified in labeling. Do not use the cable if it appears damaged.

www.medtronic.eu

Europe

Medtronic International Trading Sàrl
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00

United Kingdom/Ireland

Medtronic Limited
Building 9
Croxley Green Business Park
Hatters Lane
Watford
Herts WD18 8WW
www.medtronic.co.uk
Tel: +44 (0)1923 212213
Fax: +44 (0)1923 241004

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UC201200648 EE © Medtronic 2011. Wszelkie prawa zastrzeżone. Wydrukowano w Europie