CleanCart® cartridges

Efficient and safe care of the fluid path

With CleanCart A and CleanCart C cartridges, regular dialysis machine maintenance and disinfection assure proper dialysis machine function and the best treatment delivery requiring less machine operation from clinical staff.

- Less space for storage
- Less heavy lifting for staff
- No risk from liquid chemicals
- Easy to safely streamline disinfection combining heat disinfection with decalcification and cleaning
- Removes organic deposits, fats and proteins with weekly use of CleanCart A cartridge
- Removes precipitated calcium and magnesium salts with daily use of CleanCart C cartridge
- Recyclable polypropylene container for less environmental impact

Applying dialysis disinfectants becomes a convenient and simple operation. So you can make the most of your time with patients.

CleanCart® cartridge

CleanCart® A

The CleanCart A cartridge contains 13 grams of anhydrous sodium carbonate powder (in compliance with Ph.Eur., USP). It is used in combination with the heat disinfection program to clean the dialysis machine fluid path. When circulated throughout the fluid path, the sodium carbonate solution removes organic deposits, fats and proteins.

CleanCart® C

The CleanCart C cartridge contains 32 grams of citric acid anhydrate powder (in compliance with Ph.Eur., USP). The CleanCart C cartridge is used in combination with the heat disinfection program to both decalcify and disinfect the dialysis machine fluid path simultaneously. When circulated throughout the fluid path, the citric acid solution removes precipitated salts of calcium and magnesium.

Technical data		
	CleanCart A	CleanCart C
Action (in combination with the heat disinfection program)	Disinfection & cleaning of the fluid path from organic deposits, fats, and proteins	Disinfection & decalcification of the fluid path
Cartridge content (PH. Eur., USP)	13 grams	32 grams
Total cartridge weight	Approx. 32 grams	Approx. 50 grams
Cartridge and label material	Recyclable polypropylene (PP)	Recyclable polypropylene (PP)
Recommended usage	Refer to operator's manual for respective dialysis machine	Refer to operator's manual for respective dialysis machine
PH (prepared solution)	Approx. 11	Approx. 2
Residual test needed	No	No
Biodegradable content	Yes	Yes
Storage conditions	Store below +30°C	Store below +30°C
Shelf life	24 months from date of manufacture	24 months from date of manufacture
Product code	114010	114011
Packaging	Box of 33 cartridges	Box of 33 cartridges
Box dimensions	279 x 247 x 81 mm	279 x 247 x 81 mm
Total box weight	1200 grams	1800 grams

Gambro® is a registered trademark of Gambro Lundia AB
CleanCart® is a trademark of Gambro Hospal (Switzerland) Ltd registered in Europe and in other countries
Not for use in the US.

Gambro Lundia AB
PO Box 10101
SE-22010 Lund
Sweden
Phone + 46 46 16 90 00
partner@gambro.com
www.gambro.com

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Gambro and CleanCart are trademarks belonging to the Gambro Group.