# Alexis® O

## Wound Protector/Retractor

### Wound Protection / Incision Retraction Device

#### Instructions For Use

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#### Mode d'emploi

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#### Gebruiksaanweisung

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#### Brugsanvisning

#### Brugsanvisning

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
#### Brugsanvisning

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#### Kullanma Talimatı

ENGLISH
<p><b>INDICATIONS:</b> The Alexis O wound protector/retractor is a device that allows the surgeon to access the abdominal cavity through an atraumatically retracted wound, providing maximum exposure with minimum incision size. In addition to incision retraction, it is intended to protect against wound contamination during both laparoscopic and open surgery. The small and medium Alexis O wound protectors/retractors are capable of temporarily closing an incision to maintain pneumoperitoneum during laparoscopic surgery, or to serve as an additional trocar port site. In addition, indicated sizes may be used to access the thoracic cavity or other soft tissue retraction during cardiac and general surgical procedures.</p> <p><b>DESCRIPTION:</b> The Alexis O wound protector/retractor is a one piece design, consisting of a wound sheath with a ring at both ends.</p> <p><b>CONTRAINDICATIONS:</b> Use of the Alexis O wound protector/retractor is contraindicated when, in the judgment of the physician, use of this device would be contrary to the best interest of the patient.</p> <p><b>PRECAUTIONS:</b> Sweep area to ensure tissue is not trapped between ring and tissue layer.</p> <p><b>WARNINGS:</b></p> <ul style="list-style-type: none"> <li>Do not use if the product or sterile package is damaged.</li> <li>This device was designed and tested for single patient use only. Do not reuse, reprocess or resterilize this device. Reuse, reprocessing or resterilization may alter the structural and/or functional integrity of this device which may result in patient injury, infection, illness or death. Risk of residual contamination and resterilization failure may lead to patient injury, infection, illness or death.</li> <li>Avoid sharp instruments contacting the device, as it may damage the device or cause patient injury.</li> </ul> <p><b>INSTRUCTIONS FOR USE:</b></p> <ul style="list-style-type: none"> <li><b>PREPARE SURGICAL SITE</b> <ul style="list-style-type: none"> <li>Prepare surgical site according to the hospital's standard procedures, making sure the skin is clean and dry.</li> </ul> </li> <li><b>MARK INCISION</b> <ul style="list-style-type: none"> <li>Mark an appropriate length incision line using a sterile skin marker (Illustration <b>1</b>).</li> </ul> </li> <li><b>MAKE AN INCISION</b> <ul style="list-style-type: none"> <li>Make an incision along marked incision line (Illustration <b>2</b>).</li> </ul> </li> <li><b>INSERT WOUND RETRACTOR</b> <ul style="list-style-type: none"> <li>Where applicable, dissect the appropriate space to prepare for proper placement of the green ring. Collapse the green ring by squeezing the two sides together and inserting it through the incision until desired position is reached (Illustration <b>3</b>).</li> <li>For applicable models, keep the tether hanging outside of the sheath and patient (Illustration <b>3</b>).</li> <li>After retraction, sweep area to ensure tissue is not trapped between ring and tissue layer.</li> </ul> </li> <li><b>DEPLOY WOUND RETRACTOR</b> <ul style="list-style-type: none"> <li>Grasp the white ring at opposite ends, perpendicular to the incision. Pull up until the green ring sits tightly against the tissue layer.</li> <li>Grasp the white ring at the 10 and 2 o'clock position. Roll inward or outward until desired retraction is achieved, and the outer ring sits against the skin (Illustration <b>4</b>).</li> <li>Gently pull up on the white ring to ensure it is securely anchored in place (Illustration <b>5</b>).</li> </ul> </li> <li><b>USE SMALL OR MEDIUM WOUND RETRACTOR TO CLOSE A WOUND AND ESTABLISH INSUFFLATION</b> <ul style="list-style-type: none"> <li>To close the wound and establish insufflation, unroll and rotate the white ring so that it twists the sheath until it forms an airtight seal (Illustration <b>6</b>).</li> <li>Grasp the white ring at the 10 and 2 o'clock position. Roll inward or outward until desired retraction is achieved, and the outer ring sits against the skin (Illustration <b>7</b>).</li> <li>Gently pull up on the white ring to ensure it is securely anchored in place (Illustration <b>8</b>).</li> </ul> </li> <li><b>USE SMALL OR MEDIUM WOUND RETRACTOR AS A TROCAR PORT SITE</b> <ul style="list-style-type: none"> <li>Place a 5mm trocar cannula in the opening of the retractor (Illustration <b>9</b>).</li> <li>While holding the trocar cannula in place, unroll and rotate the white ring so that it twists the sheath around the cannula until it forms an airtight seal (Illustration <b>9</b>).</li> <li>Grasp the white ring at the 10 and 2 o'clock position. Roll inward or outward until desired retraction is achieved, and the outer ring sits against the skin (Illustration <b>9</b>).</li> <li>Gently pull up on the white ring to ensure it is securely anchored in place (Illustration <b>9</b>).</li> </ul> </li> <li><b>REMOVE WOUND RETRACTOR</b> <ul style="list-style-type: none"> <li>If necessary, loosen and unwind the Alexis O wound protector/retractor.</li> <li>Grasp the green ring and pull ring up and out of the incision (Illustration <b>10</b>).</li> <li>For applicable models, the Alexis O wound protector/retractor can be removed by grasping the tether and pulling upward, allowing the green ring to be removed from the incision (Illustration <b>10</b>).</li> </ul> </li> </ul>

#### HOW SUPPLIED:

This device is supplied sterile and for single patient use.

**STORAGE AND HANDLING:**  
Handle with care. Store at room temperature. Avoid prolonged temperature excursions.

#### SPECIFICATIONS:

Model Number	Size	Description	Incision Size	Wound Retraction	Thoracic Retraction	Maintain Pneumoperitoneum	Trocar Site
C8401	S	Small	2.5-6 cm	x	x	x	x
C8402	M	Medium	5-9 cm	x	x	x	x
C8403	L	Large	9-14 cm	x			
C8404	XL	Extra Large	11-17 cm	x			
C8405	XXL	Extra Extra Large	17-25 cm	x			

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ENGLISH
<p><b>INDICATIONS:</b> The Alexis O wound protector/retractor is a device that allows the surgeon to access the abdominal cavity through an atraumatically retracted wound, providing maximum exposure with minimum incision size. In addition to incision retraction, it is intended to protect against wound contamination during both laparoscopic and open surgery. The small and medium Alexis O wound protectors/retractors are capable of temporarily closing an incision to maintain pneumoperitoneum during laparoscopic surgery, or to serve as an additional trocar port site. In addition, indicated sizes may be used to access the thoracic cavity or other soft tissue retraction during cardiac and general surgical procedures.</p> <p><b>DESCRIPTION:</b> The Alexis O wound protector/retractor is a one piece design, consisting of a wound sheath with a ring at both ends.</p> <p><b>CONTRAINDICATIONS:</b> Use of the Alexis O wound protector/retractor is contraindicated when, in the judgment of the physician, use of this device would be contrary to the best interest of the patient.</p> <p><b>PRECAUTIONS:</b> Sweep area to ensure tissue is not trapped between ring and tissue layer.</p> <p><b>WARNINGS:</b></p> <ul style="list-style-type: none"> <li>Do not use if the product or sterile package is damaged.</li> <li>This device was designed and tested for single patient use only. Do not reuse, reprocess or resterilize this device. Reuse, reprocessing or resterilization may alter the structural and/or functional integrity of this device which may result in patient injury, infection, illness or death</li></ul>