
Alexis® Laparoscopic System

Instructions For Use Instrucciones de uso Mode d'emploi Istruzioni per l'uso Gebrauchsanweisung Gebruiksaanwijzing Bruksanvisning Brugsanvisning Käyttöohjeet Bruksanvisning Instruções de utilização Instrukcja używania Kullanma Talimatı

INDICATIONS:
The Alexis O wound protector/retractor is a device that allows the surgeon to access the abdominal cavity through an atraumatically retracted wound, providing maximum exposure with minimum incision size. In addition to incision retraction, it is intended to protect against wound contamination during both open and laparoscopic surgery. The Alexis O wound protector/retractor may be used to access the thoracic cavity or other soft tissue retraction during cardiac and general surgical procedures. The laparoscopic cap is intended to be used during laparoscopic surgery to seal the retracted wound or to convert the retracted wound site into an additional trocar port for 12mm trocars.

DESCRIPTION:
The Alexis Laparoscopic System includes an Alexis O wound protector/retractor and laparoscopic cap. The Alexis O wound protector/retractor consists of a wound sheath with a ring at both ends.

KEY - ILLUSTRATIONS AND TERMINOLOGY (Illustration 1):
A. Alexis laparoscopic cap
B. Alexis O wound protector/retractor

CONTRAINDICATIONS:
Use of the Alexis Laparoscopic System is contraindicated when, in the judgement of the physician, use of this device would be contrary to the best interest of the patient.

PRECAUTIONS:
• Sweep area to ensure tissue is not trapped between ring and tissue layer.

WARNINGS:
• Do not use if the product or sterile package is damaged.
• This device was designed and tested for single patient use only. Do not reuse, reprocess or resterilize this device. Reuse, reprocessing or resterilization may alter the structural and/or functional integrity of this device which may result in patient injury, infection, illness or death. Risk of residual contamination and resterilization failure may lead to patient injury, infection, illness or death.
• Avoid sharp instruments contacting the device, as it may damage the device or cause patient injury.

INSTRUCTIONS FOR USE:
• **PREPARE SURGICAL SITE**
◦ Prepare the surgical site according to the hospital's standard procedures, making sure the skin is clean and dry.
• **MARK INCISION**
◦ Mark an appropriate length incision line using a sterile skin marker (Illustration 2).
• **MAKE AN INCISION**
◦ Make an incision along marked incision line (Illustration 3).
• **INSERT WOUND RETRACTOR**
◦ Where applicable, dissect the appropriate space to prepare for proper placement of the green ring. Collapse the green ring by squeezing the two sides together and inserting it through the incision until desired position is reached (Illustration 4).
◦ For applicable models, keep the tether hanging outside the sheath and patient.
◦ Before retraction, sweep area to ensure tissue is not trapped between green ring and tissue layer.
• **DEPLOY WOUND RETRACTOR**
◦ Grasp the white ring at opposite ends, perpendicular to the incision. Pull up until the green ring sits tightly against the tissue layer.
◦ Grasp the white ring at the 10 and 2 o'clock position. Roll inward or outward until desired retraction is achieved, and the outer ring sits against the skin (Illustration 5).
◦ Gently pull up on the white ring to ensure it is securely anchored in place (Illustration 6).
• **ATTACH LAPAROSCOPIC CAP**
◦ Place the laparoscopic cap on top of the Alexis O wound protector/retractor's white ring (Illustration 7).
◦ Push the outer edge of the laparoscopic cap and white ring of the Alexis O wound protector/retractor together until the white ring is fully enclosed by the laparoscopic cap (Illustration 8).
◦ Use as a trocar port site or to maintain insufflation (Illustration 9).
• **REMOVE LAPAROSCOPIC SYSTEM**
◦ Remove the laparoscopic cap from the Alexis O wound protector/retractor (Illustration 10). If necessary, loosen and unwind the Alexis O wound protector/retractor.
◦ Grasp the green ring and pull up and out of the incision (Illustration 11).
◦ For applicable models, the Alexis O wound protector/retractor can be removed by grasping the tether and pulling upward, allowing the green ring to be removed from the incision.

HOW SUPPLIED:
This device is supplied sterile and for single patient use.
STORAGE AND HANDLING:
Handle with care. Store at room temperature. Avoid prolonged temperature excursions.
SPECIFICATIONS:

Model Number	Size	Description	Incision Size	Tether
C8501	S	Small	2.5-6cm	X
C8502	M	Medium	5-9cm	

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